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

Figure 1
Kinetics of light emission for uninfected and infected cells using the MycoAlert™ PLUS Kit. The B/A ratio indicates the presence or absence of mycoplasma.


Figure 2
Performance of MycoAlert™ PLUS Assay compared to MycoAlert™ Assay. A dilution series of MycoAlert™ Assay Control demonstrates the increased sensitivity of various luminometer models when using MycoAlert™ PLUS Assay, compared to the first generation MycoAlert™ Assay.

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Table

Mycoplasma species validated down to 1 CFU/ml at a 95% confidence level
[≥23 out of 24 positive]

Species	Claimed Sensitivity	Probit Analysis of Validation Data
	1 CFU/ml	CFU/ml
<i>A. laidlawii</i> *	+	<1.0
<i>M. arginini</i>	+	0.15
<i>M. fermentans</i>	+	0.29
<i>M. hominis</i> *	+	<1.0
<i>M. hyorhinis</i>	+	0.51
<i>M. orale</i>	+	0.11
<i>M. pneumoniae</i>	+	0.02
<i>M. salivarium</i>	+	1.41
<i>S. citri</i>	+	0.37
<i>M. gallisepticum</i>	+	not yet validated
<i>M. synoviae</i>	+	not yet validated

*For 1 CFU/ml 24 out of 24 positive; 0.5 CFU/ml and 0.1 CFU/ml not tested.

For use in quality control/manufacturing process only.
MYCOTOOL is a trademark of Roche.

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
Figure 3
Effect of MycoZap™ Elimination Reagent on cell viability and mycoplasma removal. The MycoZap™ Reagent treatment eliminates mycoplasma in as few as 4 days with minimal impact on cell viability. Mycoplasma detection was performed with MycoAlert™ Assay and viability was determined by ViaLight™ Plus Assay.

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	Mycoplasma Only Solution	All-in-one Solutions	
	Prophylactic	Plus-CL	Plus-PR
Prevention against mycoplasma	■	■	■
Prevention against - Gram (+) bacteria - Gram (-) bacteria - Fungi - Yeast	No; but can be used in combination with other antibiotic formula of choice detection and prevention	■	■
Suited for primary cells	■	■	■
Suited for cell lines	■	■	

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*MycoZap™ Spray is only available in Austria, Belgium, Denmark, France, Germany, Italy, Luxembourg, Netherlands, Poland, Spain, Sweden, Switzerland, UK and distributed countries on the European continent. For availability in other countries, please contact our Customer Service (for contact details please refer to back page).

 For more information, please visit www.lonza.com/mycoplasma

References

- Drexler HG, Uphof CC (2002) Mycoplasma contamination of cell cultures: Incidence, sources, effects, detection, elimination, prevention. *Cytotechnology* 39: 75–90
- Koshimizu K, Kotani H (1981). In: Procedures for the Isolation and Identification of Human, Animal and Plant Mycoplasmas (Nakamura, M., ed.), *Saikou, Tokyo*, 87–102.
- Gong H, Zölzer F, von Recklinghausen G, Rössler J, Breit S, Havers W, Fotsis T, Schweigerer L (1999) Arginine deiminase inhibits cell proliferation by arresting cell cycle and inducing apoptosis. *Biochem Biophys Res Comm* 261: 10–14.
- Ben-Menachem G, Mousa A, Brenner T, Pinto F, Zähringer U, Rottem S (2001) Choline deficiency induced by Mycoplasma fermentans enhances apoptosis of rat astrocytes. *FEMS Microbiol Letters* 201: 157–162
- McGarrity MF, Vanaman V, Sarama J (1984) Cyto genetic effects of mycoplasma infection of cell cultures: a review. *In Vitro* 20: 1–18
- Sokolova IA, Vaughan ATM, Khodarev NN (1998) Mycoplasma infection can sensitize host cells to apoptosis through contribution of apoptotic-like endonuclease[s]. *Immunol Cell Biol* 76: 526–534
- Doersen CJ, Stanbridge EJ (1981) Effects of mycoplasma contamination on phenotypic expression of mitochondrial mutants in human cells. *Mol Cell Biol* 1: 321–329
- Stanbridge EJ (1971) Mycoplasmas and cell cultures. *Bacteriological Reviews* 35: 206–227
- Darin N, Kadhom N, Brière JJ, Chretien D, Bébéar CM, Rötig A, Munnich A, Rustin P (2003) Mitochondrial activities in human cultured skin fibroblasts contaminated by Mycoplasma hyorhinis. *BMC Biochem* 4: 15
- Rottem S (2003) Interaction of mycoplasmas with host cells. *Physiol Rev* 83: 417–432
- Miller CJ, Kassem HS, Pepper SD, Hey Y, Ward TH, Margison GP. (2003) Mycoplasma infection significantly alters microarray gene expression profiles. *Biotechniques* 35: 812–814

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Lonza Cologne GmbH – 50829 Cologne, Germany

For research use only. Not for use in diagnostic procedures.
MycoTOOL is a trademark of Roche.

Unless otherwise noted, all trademarks herein are marks of the Lonza Group or its affiliates. The information contained herein is believed to be correct and corresponds to the latest state of scientific and technical knowledge. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information and no warranty is expressed or implied concerning the use of these products. The buyer assumes all risks of use and/or handling. Any user must make his own determination and satisfy himself that the products supplied by Lonza Group Ltd or its affiliates and the information and recommendations given by Lonza Group Ltd or its affiliates are (i) suitable for intended process or purpose, (ii) in compliance with environmental, health and safety regulations, and (iii) will not infringe any third party's intellectual property rights.

© Copyright 2013, Lonza Cologne GmbH. All rights reserved.
BR-MycoProd Chinese 06/13 CD-BR022_CN