Nucleofector™ Technology

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Nucleofector™ 技术

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Nucleofector™ Technology – the Superior Non-viral Method


Figure 1. Transfection of the human natural killer cell line NKL using traditional electroporation and Nucleofection. 5×10^6 NKL cells were transfected with $2.5 \mu\text{g}$ of pmaxGFP™ Vector. Nucleofection: Nucleofector™ Solution V; Program 0-017. Competitor B electroporation: 25 mV, $96 \mu\text{F}$. Transfection efficiency was monitored by flow cytometry after 24 hours. Cells transfected by Nucleofection show a significantly better transfection efficiency compared to cells transfected by traditional electroporation. Cell viability, as measured 18 hours after transfection, was also superior using Nucleofection. (Data courtesy of Dr. John Coligan, Laboratory of Immunogenetics, NIH/NIAID, Rockville, MD, USA. *J Immunol Methods* [2004] 284: 133-140.)

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DNA Delivery Straight Into the Nucleus


Figure 2. Normal human dermal fibroblasts (neonatal) were transfected with $2.5 \mu\text{g}$ TMR-labeled plasmid DNA encoding eGFP. After 2 hours, cells were fixed with 3.5% PFA and analyzed by confocal microscopy. TMR label is shown in (A), GFP fluorescence in (B), DAPI nuclear staining in (C) and a merge of all three fluorescent labels in (D).

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
Next Generation 4D-Nucleofector™ Platform

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
Figure 3. Comparison of 4D-Nucleofector® Units of first (Gen 1) and next (Gen 2) generation for all functional units. Bars represent transfection efficiencies for maxGFP™ protein and viabilities for different cell types: Jurkat cells (X Unit and 96-well Unit), NHDF-neo (Y Unit and LV Unit with 1 mL cartridge) and suspension HEK293 cells (LV Unit with LV cartridge).

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Transferability Between Nucleofection Conditions Between Different Formats


Figure 4. Various primary cells were transfected in the two Nucleocuvette™ vessel formats (20μL and 100μL) using the indicated programs. Twenty-four hours post Nucleofection cells were analyzed for transfection efficiency (flow cytometry) and viability (cell number normalized to no program control).

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
Figure 5. Efficient adherent Nucleofection of neurons in 24-well culture plates. Mouse cortical neurons were seeded into poly-D-lysine coated 24-well plates (1×10^5 cells/well) After 6 DIV, cells were transfected with pmaxGFP™ Vector using the AD1 4D-Nucleofector™ Y Kit. One day post Nucleofection, cells were stained by MAP2 antibody (red) and analyzed by fluorescence microscopy for maxGFP™ protein expression.



Figure 6. Human umbilical vein endothelial cells (HUVEC) were isolated and plated in passage 1 into collagen-coated 24-well plates at a density of 50,000 cells/well. After 1DIV cells were transfected with 16µg pmaxGFP™ Vector using AD1 4D-Nucleofector™ Y Solution and program CA-215. Cells were analyzed for maxGFP™ Protein expression after 24h. (Data kindly provided by M. Sauvage, Pharmaceutical Industry, FR)

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
Optimization of Nucleofection Conditions Within Just One Experiment

Examples of cell lines that have been optimized by customers using the Cell Line Optimization 96-well Nucleofector™ Kit.

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

Same conditions used for the 96-well Shuttle™ and the HT Nucleofector™ Systems. The HT Nucleofector™ System works with 96-well Shuttle™ Parameters, thus the full spectrum of already optimized protocols is available for the HT Nucleofector™ System.

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
Two formats available. (A) 1mL Nucleocuvette™ Cartridge: 1mL filling volume for up to 1×10^8 cells (manual filling via sterile injection port) (B) LV Nucleocuvette™ Cartridge: Up to 20mL processing volume (in 1 mL steps) for up to 1×10^9 cells (automatic filling via reservoirs or bags)

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货号	产品名	描述	规格
Nucleofector™ Devices and Systems			
384-well Nucleofector™ System			
AAU-1001	384-well Nucleofector™ System	Includes power supply, plate handler, laptop, and software	
4D-Nucleofector™ System			
AAF-1003B	4D-Nucleofector™ Core Unit		
AAF-1002L	4D-Nucleofector™ LV Unit	Including 2 LV reservoir racks. Requires the core unit to build complete system	
AAF-1003X	4D-Nucleofector™ X Unit	Requires the Core Unit to build complete system	
AAF-1003Y	4D-Nucleofector™ Y Unit	Requires the Core Unit to build complete system	
AAF-1003S	4D-Nucleofector™ 96-well Unit	Requires the Core Unit to build complete system	
Nucleofector™ Kits			
V4XC-1024	SE Cell Line 4D-Nucleofector™ X Kit L	100 µL Nucleocuvette™ Vessel	24 reactions
V4XC-2024	SF Cell Line 4D-Nucleofector™ X Kit L	100 µL Nucleocuvette™ Vessel	24 reactions
V4XC-2032	SF Cell Line 4D-Nucleofector™ X Kit S	20 µL Nucleocuvette™ Strip	32 reactions (16-well)
V4XC-3024	SG Cell Line 4D-Nucleofector™ X Kit L	100 µL Nucleocuvette™ Vessel	24 reactions
V4XC-9064	Cell Line Optimization 4D-Nucleofector™ X Kit	20 µL Nucleocuvette™ Strip	64 reactions (16-well)
V4XP-2024	P2 Primary Cell 4D-Nucleofector® X Kit L	100 µL Nucleocuvette™ Vessel	24 reactions
V4XP-2032	P2 Primary Cell 4D-Nucleofector® X Kit S	20 µL Nucleocuvette™ Strip	32 reactions (16-well)
V4XP-3024	P3 Primary Cell 4D-Nucleofector® X Kit L	100 µL Nucleocuvette™ Vessel	24 reactions
V4XP-3032	P3 Primary Cell 4D-Nucleofector® X Kit S	20 µL Nucleocuvette™ Strip	32 reactions (16-well)
V4XP-4024	P4 Primary Cell 4D-Nucleofector® X Kit L	100 µL Nucleocuvette™ Vessel	24 reactions
V4XP-9096	Primary Cell Optimization 4D-Nucleofector™ X Kit	20 µL Nucleocuvette™ Strip	96 reactions (16-well)
V4YP-1A24	AD1 4D-Nucleofector™ Y Kit	24-well Dipping Electrode	24 reactions
V4YP-2A24	AD2 4D-Nucleofector™ Y Kit	24-well Dipping Electrode	24 reactions
V4SC-2096	SF Cell Line 96-well Nucleofector™ Kit	20 µL Nucleocuvette™ Plate	96 reactions (96-well)
V4SP-3096	P3 Primary Cell 96-well Nucleofector™ Kit	20 µL Nucleocuvette™ Plate	96 reactions (96-well)
V4LP-3002	P3 Primary Cell 4D-Nucleofector™ LV Kit L	1 mL Nucleocuvette™ Cartridge	2 reactions
V4LP-3020	P3 Primary Cell 4D-Nucleofector™ LV Kit XL	LV Nucleocuvette™ Cartridge	1 reaction
V5SP-3002	P3 Primary Cell 384-well Nucleofector™ Kit	20 µL Nucleocuvette™ Plate	768 reactions (384-well)
VCA-1003	Cell Line Nucleofector™ Kit V	100 µL aluminum cuvette	25 reactions
VPA-1002	Human T Cell Nucleofector™ Kit	100 µL aluminum cuvette	25 reactions
GMP-grade consumables			
G4LP3-22500	TheraPEAK™ P3 Nucleofector® Solution Set		22.5 mL + 5 mL
G4LP3-126000	TheraPEAK™ P3 Nucleofector® Solution Set		126 mL + 28 mL
G4LN-7501	TheraPEAK™ 1 mL Nucleocuvette® Cartridge - coming soon		5 reactions
G4LN-7520W	TheraPEAK™ LV Nucleocuvette® Cartridge, weldable tubing		5 reactions
G4LR-1501W	TheraPEAK™ 4D-Nucleofector® Reservoir, weldable tubing		5 reactions

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