


## 小动物活体光学成像技术在心血管疾病研究中的应用

Revvity 小动物活体光学成像技术已在生命科学基础研究、临床前医学研究及药物研发等领域得到广泛应用。在众多应用领域中，心血管疾病是目前研究的热点领域之一。常用于心血管研究的光学标记方法包括：1. 利用荧光虫荧光素酶（Firefly Luciferase）或荧光蛋白作为报告基因，通过标记特定基因而研究心血管疾病相关基因在心血管疾病中的作用；2. 通过外源注射功能性荧光探针，观测心血管疾病发展过程中的分子变化。应用小动物活体光学成像技术进行心血管疾病研究主要集中在三个方面：1. 应用生物发光技术，研究细胞治疗心血管疾病效果；2. 应用功能性荧光探针，了解疾病发展的分子机理和药物治疗心血管疾病的效果；3. 心血管疾病相关基因的作用及其治疗。下面将对一些具体实例进行阐述：

### 一. 应用生物发光技术，研究细胞治疗心血管疾病效果

通过细胞示踪，可以了解细胞在心脏的分布并且了解心脏的调控机制。新型的报告基因（triple-fusion, TF）包括萤火虫荧光素酶，单体红色荧光蛋白和缩短的胸苷激酶。TF 报告基因体外转染鼠胚胎干细胞，然后注射胚胎干细胞到裸鼠心肌层。Revvity 的 IVIS 成像系统显示的生物发光成像和正电子发射层断成像（PET）有很强相关性，并且能够监测移植后胚胎干细胞存活，增殖和迁移（下图）。




(Lee & Wu, Stem Cell Migration, 2011)

上图：移植含有 TF 报告基因的胚胎干细胞到无胸腺的裸鼠心脏。使用生物发光（a，上图）和 PET（a，下图）纵向监测胚胎干细胞在心脏内移植后的存活和增殖。（b & c），定量分析生物发光和 PET 成像的强度，之间有很强的关联性。


## 二. 应用功能性荧光探针，了解疾病发展的分子机理和药物治疗心血管疾病效果

血管疾病发展过程中伴随着分子表达的改变，比如基质金属蛋白酶的表达特异性升高，因此通过观测体内分子的改变能够判断心血管疾病的发展和治疗效果。研究发现基质金属蛋白酶 (Metalloproteinases, MMPs) 在动脉粥样硬化中高度聚集，应用 Revvity 的 MMPsense 探针并通过 IVIS Spectrum 仪器成像显示人动脉粥样硬化中金属蛋白酶的活性。下图说明，应用 MMPsense 探针能够检测人颈动脉内膜切除手术的动脉粥样硬化斑块。




(Wallis et al, Circulation. 2009)

在 apoE<sup>-/-</sup> 心血管疾病小鼠模型中，大动脉根的组织蛋白酶活性 (Protease Sensor, PS) 升高 (左下图)。使用 Revvity 的 Prosense, IntegrinSense 和 CatBFAST 探针分别对心血管疾病中大动脉根的组织蛋白酶和整联蛋白  $\alpha\beta3$  活性进行观测。如右下图所示，活体和体外组织的成像结果相一致 (使用 Revvity 的 FMT 小动物活体荧光断层成像系统)。




(Nahrendorf et al., Arterioscler. Thromb. Vasc. Biol., 2009) 2012)



(Li et al., International Journal of Molecular Imaging,

Atorvastatin (阿托伐他汀)可以降低炎症性单核细胞的募集和血管细胞粘附分子的表达,而且斑块的组织蛋白酶主要来源于炎症性单核细胞和血管细胞。如下图所示,Revvity的FMT小动物活体荧光断层成像系统显示出Atorvastatin能有效治疗apoE<sup>-/-</sup>心血管疾病小鼠模型。通过Atorvastatin治疗后,主根动脉的组织蛋白酶活性下降。




(Nahrendorf et al., Arterioscler. Thromb. Vasc. Biol., 2009)

上图: Atorvastatin 降低 apoE<sup>-/-</sup>心血管疾病小鼠模型的组织蛋白酶活性,基因表达和胆固醇水平。


### 三. 心血管疾病相关基因的作用及其治疗

利用小动物活体光学成像,能够研究在活体动物水平上心血管疾病相关基因在其发展过程中的作用。利用生物发光技术进行心血管疾病相关基因的研究,主要用基因-荧光素酶的共表达载体,通过荧光素酶产生的生物发光信号反映该基因的表达情况,研究该基因的相关作用。

NF-κB 是细胞核内的转录因子,能被炎症性的细胞因子(TNF, IL-1β, IL-6, IL-8)激活。下图所示,科研人员使用荧光素酶 NF-κB 转基因小鼠研究心肌梗塞,冠状动脉结扎小鼠(Myocardial infarction)中的 NF-κB 基因表达水平相对于假性手术处理的对照小鼠(Sham)显著升高。研究表明 NF-κB 基因在心肌梗塞的发展中起了重要的作用。




使用 NF-κB 转基因小鼠进行同种异体血管化心脏移植的研究。Revvity 的 IVIS 成像系统显示野生型 Balb/c 心脏同种异种移植到萤火虫荧光素酶标记 NF-κB 转基因 C57BL/6 小鼠，NF-κB 的基因表达水平明显提高，说明 NF-κB 在同种异体移植中被激活（下图 A）。CD154 抗体（MR1）治疗可以有效地提高同种异体血管化心脏移植后受体小鼠的耐受性和存活率，并且抑制 NF-κB 基因表达水平。然而服用 CpG DNA 结合 Toll 样受体 9，提高 NF-κB 表达水平并且抑制 CD154 抗体治疗效果（下图 B）。



(Ma et al, Transplantation, 2008)

荧光素酶标记的低氧报告基因载体系统（Double oxygen-sensing vector system, DOSVS）心肌内注射到小鼠心脏。通过间断性缺氧实验，研究人员使用 Revvity 的 IVIS 成像系统观测低氧报告基因表达的变换。下图所示在心脏缺氧（On）情况下，低氧报告基因表达升高；然而在常氧条件（Off）下，低氧报告基因表达下降。实验证明低氧报告基因可以用来监测心脏缺氧或者局部缺血的产生和发展。



Fomicheva et al, Molecular Therapy, 2008

## Cardiovascular Disease PAPERS IMAGED ON THE IVIS

1. Alter J, Sennoga CA, Lopes DM, Eckersley RJ, Wells DJ. Microbubble stability is a major determinant of the efficiency of ultrasound and microbubble mediated in vivo gene transfer. *Ultrasound in medicine & biology*. Jun 2009;35(6):976-984.
2. Bai X, Yan Y, Song YH, Seidensticker M, Rabinovich B, Metzle R, Bankson JA, Vykoukal D, Alt E. Both cultured and freshly isolated adipose tissue-derived stem cells enhance cardiac function after acute myocardial infarction. *European heart journal*. Feb 2010;31(4):489-501.
3. Barth AS, Kizana E, Smith RR, Terrovitis J, Dong P, Leppo MK, Zhang Y, Miake J, Olson EN, Schneider JW, Abraham MR, Marban E. Lentiviral vectors bearing the cardiac promoter of the Na<sup>+</sup>-Ca<sup>2+</sup> exchanger report cardiogenic differentiation in stem cells. *Mol Ther*. May 2008;16(5):957-964.
4. Boddington SE, Henning TD, Jha P, Schlieve CR, Mandrussow L, DeNardo D, Bernstein HS, Ritner C, Golovko D, Lu Y, Zhao S, Daldrop-Link HE. Labeling human embryonic stem cell-derived cardiomyocytes with indocyanine green for noninvasive tracking with optical imaging: an FDA-compatible alternative to firefly luciferase. *Cell transplantation*. 2010;19(1):55-65.
5. Bogdanov AA, Jr., Lin CP, Kang HW. Optical imaging of the adoptive transfer of human endothelial cells in mice using anti-human CD31 monoclonal antibody. *Pharmaceutical research*. Jun 2007;24(6):1186-1192.
6. Cao F, Lin S, Xie X, Ray P, Patel M, Zhang X, Drukker M, Dylla SJ, Connolly AJ, Chen X, Weissman IL, Gambhir SS, Wu JC. In Vivo Visualization of Embryonic Stem Cell Survival, Proliferation, and Migration After Cardiac Delivery. *Circulation*. February 21, 2006 2006;113(7):1005-1014.
7. Cao YA, Bachmann MH, Beilhack A, Yang Y, Tanaka M, Swijnenburg RJ, Reeves R, Taylor-Edwards C, Schulz S, Doyle TC, Fathman CG, Robbins RC, Herzenberg LA, Negrin RS, Contag CH. Molecular imaging using labeled donor tissues reveals patterns of engraftment, rejection, and survival in transplantation. *Transplantation*. Jul 15 2005;80(1):134-139.
8. Chimenti I, Smith RR, Li TS, Gerstenblith G, Messina E, Giacomello A, Marban E. Relative Roles of Direct Regeneration Versus Paracrine Effects of Human Cardiosphere-Derived Cells Transplanted Into Infarcted Mice. *Circulation research*. Jan 28 2010.
9. Duysen EG, Lockridge O. Whole body and tissue imaging of the butyrylcholinesterase knockout mouse injected with near infrared dye labeled butyrylcholinesterase. *Chemico-biological interactions*. Mar 7 2008.
10. Fomicheva EV, Turner, II, Edwards TG, Hoff J, Arden E, D'Alecy LG, Metzger JM. Double Oxygensensing Vector System for Robust Hypoxia/Ischemia-regulated Gene Induction in Cardiac Muscle In Vitro and In Vivo. *Mol Ther*. Jun 24 2008.
11. Fujishiro J, Kawana H, Inoue S, Shimizu H, Yoshino H, Hakamata Y, Kaneko T, Murakami T, Hashizume K, Kobayashi E. Efficiency of adenovirus-mediated gene transduction in heart grafts in rats. *Transplantation proceedings*. Jan-Feb 2005;37(1):67-69.
12. Gruber PJ, Li Z, Li H, Worrada D, Huang B, Abdullah I, Wang W, El-Deiry W, Ferrari VA, Zhou R. In vivo imaging of MLC2v-luciferase, a cardiac-specific reporter gene expression in mice. *Acad Radiol*. September 2004;11(9):1022-1028.
13. Hendry SL, 2nd, van der Bogt KE, Sheikh AY, Arai T, Dylla SJ, Drukker M, McConnell MV, Kutschka I, Hoyt G, Cao F, Weissman IL, Connolly AJ, Pelletier MP, Wu JC, Robbins RC, Yang PC. Multimodal evaluation of in vivo magnetic resonance imaging of myocardial restoration by mouse embryonic stem cells. *The Journal of thoracic and cardiovascular surgery*. Oct 2008;136(4):1028-1037 e1021.
14. Hingtgen SD, Li Z, Kutschke W, Tian X, Sharma RV, Davisson RL. Superoxide Scavenging and AKT Inhibition in the Myocardium Ameliorate Pressure Overload-induced NF- $\kappa$ B Activation and Cardiac Hypertrophy. *Physiol Genomics*. Jan 26 2010.
15. Hsiang CY, Chen YS, Ho TY. Nuclear factor- $\kappa$ B bioluminescence imaging-guided transcriptomic analysis for the assessment of host-biomaterial interaction in vivo. *Biomaterials*. Jun 2009;30(17): 3042-3049.
16. Hsu M, Muchova L, Morioka I, Wong RJ, Schroder H, Stevenson DK. Tissue-specific effects of statins on the expression of heme oxygenase-1 in vivo. *Biochemical and biophysical research communications*. May 12 2006;343(3):738-744.
17. Huang M, Chan DA, Jia F, Xie X, Li Z, Hoyt G, Robbins RC, Chen X, Giaccia AJ, Wu JC. Short hairpin RNA interference therapy for ischemic heart disease. *Circulation*. Sep 30 2008;118(14 Suppl):S226-233.
18. Jearawiriyapaisarn N, Moulton HM, Buckley B, Roberts J, Sazani P, Fucharoen S, Iversen PL, Kole R. Sustained dystrophin expression induced by peptide-conjugated morpholino oligomers in the muscles of mdx mice. *Mol Ther*. Sep 2008;16(9):1624-1629.

19. Karlsen A, Retterstol L, Laake P, Paur I, Kjolsrud-Bohn S, Sandvik L, Blomhoff R. Anthocyanins inhibit nuclear factor-kappaB activation in monocytes and reduce plasma concentrations of pro-inflammatory mediators in healthy adults. *The Journal of nutrition*. Aug 2007;137(8):1951-1954.
20. Kutschka I, Chen IY, Kofidis T, Arai T, von Degenfeld G, Sheikh AY, Hendry SL, Pearl J, Hoyt G, Sista R, Yang PC, Blau HM, Gambhir SS, Robbins RC. Collagen Matrices Enhance Survival of Transplanted Cardiomyoblasts and Contribute to Functional Improvement of Ischemic Rat Hearts. *Circulation*. July 4, 2006;114(1\_suppl):I-167-173.
21. Kutschka I, Chen IY, Kofidis T, von Degenfeld G, Sheikh AY, Hendry SL, Hoyt G, Pearl J, Blau HM, Gambhir SS, Robbins RC. In vivo optical bioluminescence imaging of collagen-supported cardiac cell grafts. *J Heart Lung Transplant*. Mar 2007;26(3):273-280.
22. Kutschka I, Kofidis T, Chen IY, von Degenfeld G, Zwierzchoniewska M, Hoyt G, Arai T, Lebl DR, Hendry SL, Sheikh AY, Cooke DT, Connolly A, Blau HM, Gambhir SS, Robbins RC. Adenoviral Human BCL-2 Transgene Expression Attenuates Early Donor Cell Death After Cardiomyoblast Transplantation Into Ischemic Rat Hearts. *Circulation*. July 4, 2006;114(1\_suppl):I-174-180.
23. Li Z, Lee A, Huang M, Chun H, Chung J, Chu P, Hoyt G, Yang P, Rosenberg J, Robbins RC, Wu JC. Imaging survival and function of transplanted cardiac resident stem cells. *Journal of the American College of Cardiology*. Apr 7 2009;53(14):1229-1240.
24. Li Z, Wu JC, Sheikh AY, Kraft D, Cao F, Xie X, Patel M, Gambhir SS, Robbins RC, Cooke JP, Wu JC. Differentiation, survival, and function of embryonic stem cell derived endothelial cells for ischemic heart disease. *Circulation*. Sep 11 2007;116(11 Suppl):I46-54.
25. Ma L, Xiang Z, Sherrill TP, Wang L, Blackwell TS, Williams P, Chong A, Chari R, Yin DP. Bioluminescence Imaging Visualizes Activation of Nuclear Factor-kappaB in Mouse Cardiac Transplantation. *Transplantation*. Mar 27 2008;85(6):903-910.
26. Martinez-Fernandez A, Nelson TJ, Yamada S, Reyes S, Alekseev AE, Perez-Terzic C, Ikeda Y, Terzic A. iPS programmed without c-MYC yield proficient cardiogenesis for functional heart chimerism. *Circulation research*. Sep 25 2009;105(7):648-656.
27. Min JJ AY, Moon S, Kim YS, Park JE, Kim SM, Le UN, Wu JC, Joo SY, Hong MH, Yang DH, Jeong MH, Song CH, Jeong YH, Yoo KY, Kang KS, Bom HS. In vivo bioluminescence imaging of cord blood derived mesenchymal stem cell transplantation into rat myocardium. *Annals of Nuclear Medicine*. 2006;20(3): 165-170.
28. Nelson TJ, Martinez-Fernandez A, Yamada S, Perez-Terzic C, Ikeda Y, Terzic A. Repair of acute myocardial infarction by human stemness factors induced pluripotent stem cells. *Circulation*. Aug 4 2009;120(5):408-416.
29. Parker AL, Waddington SN, Buckley SM, Custers J, Havenga MJ, van Rooijen N, Goudsmit J, McVey JH, Nicklin SA, Baker AH. Effect of neutralizing sera on factor x-mediated adenovirus serotype 5 gene transfer. *Journal of virology*. Jan 2009;83(1):479-483.
30. Partridge J, Carlsen H, Enesa K, Chaudhury H, Zakkar M, Luong L, Kinderlerer A, Johns M, Blomhoff R, Mason JC, Haskard DO, Evans PC. Laminar shear stress acts as a switch to regulate divergent functions of NF-kappaB in endothelial cells. *Faseb J*. Nov 2007;21(13):3553-3561.
31. Peterson JR, Infanger DW, Braga VA, Zhang Y, Sharma RV, Engelhardt JF, Davisson RL. Longitudinal Non-invasive Monitoring of Transcription Factor Activation in Cardiovascular Regulatory Nuclei Using Bioluminescence Imaging. *Physiol Genomics*. Jan 29 2008.
32. Prasad KM, Xu Y, Yang Z, Toufektsian MC, Berr SS, French BA. Topoisomerase Inhibition Accelerates Gene Expression after Adeno-associated Virus-mediated Gene Transfer to the Mammalian Heart. *Mol Ther*. Apr 2007;15(4):764-771.
33. Redding KM, Chen BL, Singh A, Re RN, Navar LG, Seth DM, Sigmund CD, Tang WW, Cook JL. Transgenic mice expressing an intracellular fluorescent fusion of angiotensin II demonstrate renal thrombotic microangiopathy and elevated blood pressure. *American journal of physiology*. Jun 2010;298 (6):H1807-1818.
34. Rodriguez-Porcel M, Gheysens O, Chen IY, Wu JC, Gambhir SS. Image-Guided Cardiac Cell Delivery Using High-Resolution Small-Animal Ultrasound. *Molecular Therapy*. 2005;12(6):1142-1147.
35. Sheikh AY, Lin SA, Cao F, Cao YA, van der Bogt KE, Chu P, Chang CP, Contag CH, Robbins RC, Wu JC. Molecular Imaging of Bone Marrow Mononuclear Cell Homing and Engraftment in Ischemic Myocardium. *Stem Cells*. Jul 12 2007.
36. Swijnenburg RJ, Govaert JA, van der Bogt KE, Pearl JI, Huang M, Stein W, Hoyt G, Vogel H, Contag CH, Robbins RC, Wu JC. Timing of bone marrow cell delivery has minimal effects on cell viability and cardiac recovery after myocardial infarction. *Circ Cardiovasc Imaging*. Jan 2010;3(1):77-85.

37. Tanaka M, Swijnenburg R-J, Gunawan F, Cao Y-A, Yang Y, Caffarelli AD, de Bruin JL, Contag CH, Robbins RC. In Vivo Visualization of Cardiac Allograft Rejection and Trafficking Passenger Leukocytes Using Bioluminescence Imaging. *Circulation*. August 30, 2005 2005;112(9\_suppl):I-105-110.
38. Tillmanns J, Carlsen H, Blomhoff R, Valen G, Calvillo L, Ertl G, Bauersachs J, Frantz S. Caught in the act: in vivo molecular imaging of the transcription factor NF-kappaB after myocardial infarction. *Biochemical and biophysical research communications*. Apr 14 2006;342(3):773-774.
39. Tolar J, Wang X, Braunlin E, McElmurry RT, Nakamura Y, Bell S, Xia L, Zhang J, Hu Q, Panoskaltis-Mortari A, Zhang J, Blazar BR. The host immune response is essential for the beneficial effect of adult stem cells after myocardial ischemia. *Experimental hematology*. Apr 2007;35(4):682-690.
40. van Amerongen M, Bou-Gharios G, Popa E, van Ark J, Petersen A, van Dam G, van Luyn M, Harmsen M. Bone marrow-derived myofibroblasts contribute functionally to scar formation after myocardial infarction. *The Journal of pathology*. Dec 19 2007.
41. van Amerongen MJ, Bou-Gharios G, Popa E, van Ark J, Petersen AH, van Dam GM, van Luyn MJ, Harmsen MC. Bone marrow-derived myofibroblasts contribute functionally to scar formation after myocardial infarction. *The Journal of pathology*. Feb 2008;214(3):377-386.
42. van der Bogt KE, Schrepfer S, Yu J, Sheikh AY, Hoyt G, Govaert JA, Velotta JB, Contag CH, Robbins RC, Wu JC. Comparison of transplantation of adipose tissue- and bone marrow-derived mesenchymal stem cells in the infarcted heart. *Transplantation*. Mar 15 2009;87(5):642-652.
43. van der Bogt KE, Sheikh AY, Schrepfer S, Hoyt G, Cao F, Ransohoff KJ, Swijnenburg RJ, Pearl J, Lee A, Fischbein M, Contag CH, Robbins RC, Wu JC. Comparison of different adult stem cell types for treatment of myocardial ischemia. *Circulation*. Sep 30 2008;118(14 Suppl):S121-129.
44. Wu JC, Inubushi M, Sundaresan G, Schelbert HR, Gambhir SS. Optical imaging of cardiac reporter gene expression in living rats. *Circulation*. April 2002;105(14):1631-1634.
45. Xie X, Cao F, Sheikh AY, Li Z, Connolly AJ, Pei X, Li RK, Robbins RC, Wu JC. Genetic modification of embryonic stem cells with VEGF enhances cell survival and improves cardiac function. *Cloning and stem cells*. Winter 2007;9(4):549-563.
46. Xiong YQ, Willard J, Kadurugamuwa JL, Yu J, Francis KP, Bayer AS. Real-time in vivo bioluminescent imaging for evaluating the efficacy of antibiotics in a rat *Staphylococcus aureus* endocarditis model. *Antimicrob Agents Chemother*. January 2005;49(1):380-387.
47. Yu J, Huang NF, Wilson KD, Velotta JB, Huang M, Li Z, Lee A, Robbins RC, Cooke JP, Wu JC. nAChRs mediate human embryonic stem cell-derived endothelial cells: proliferation, apoptosis, and angiogenesis. *PloS one*. 2009;4(9):e7040.
48. Zhang N, Ahsan MH, Zhu L, Sambucetti LC, Purchio AF, West DB. NF-kappaB and not the MAPK signaling pathway regulates GADD45beta expression during acute inflammation. *J Biol Chem*. June 2005;280(22):21400-21408.