

小动物活体光学成像技术在肿瘤研究中的应用


Revvity小动物活体光学成像技术已在生命科学基础研究、临床前医学研究及药物研发等领域得到广泛应用。在众多应用领域中，肿瘤研究是目前应用最为普遍的领域之一。常用于肿瘤活体成像的光学标记方法包括：1、利用萤火虫荧光素酶（Firefly Luciferase）或荧光蛋白作为报告基因，通过转基因技术体外标记肿瘤细胞而直接观测肿瘤的发展变化，或标记特定基因而研究肿瘤相关基因在肿瘤发展中的作用；2、通过外源注射功能性荧光探针，观测肿瘤发展过程中的分子事件，进而反映肿瘤的发展变化。宏观来说，应用小动物活体光学成像技术进行肿瘤研究主要集中于三个方面：1、长时间监测肿瘤生长及转移；2、抗肿瘤药物研发；3、癌症分子机理研究。下面结合一些具体实例进行阐述：

一.长时间监测肿瘤生长及转移


随着肿瘤研究的深入，应用传统方法（如卡尺测量肿瘤体积、肿瘤组织切片等）进行肿瘤研究已存在诸多限制。如进行组织切片观测前需要处死小鼠取出肿瘤组织，因此，在不同时间点或不同实验组都需要处死一批实验小鼠以获取足够的统计学数据，这样不但大大增加了实验成本，而且很难消除由于小鼠个体差异而产生的误差，无法获取可靠的重复性数据，同时，在制作切片时也无法保证实验的准确性，而利用活体光学成像技术可以对同一批小鼠进行不同时间点的长时间观测，进而大幅降低实验成本，并获取重复可靠的实验数据：




又如通过利用卡尺测量肿瘤体积的方法，只能等肿瘤发展至可以测量的程度才能开展实验，因此，无法进行肿瘤早期观测及微小转移灶的观测，而利用灵敏的生物发光成像技术在肿瘤发生早期即可进行有效观测，从而对肿瘤的整体发展过程进行全程监测，有力弥补了传统方法的缺陷。下面几个例子展示了应用生物发光成像技术长时间监测肿瘤的生长及转移。



上图：利用萤火虫荧光素酶标记 4T1 乳腺癌细胞建立皮下肿瘤模型，通过 IVIS 成像系统长期监测肿瘤的生长情况，在肿瘤细胞皮下注射的当天即可灵敏的观测到由 5 个被标记肿瘤细胞发出的光信号。定量分析显示，利用传统卡尺测量的方式（蓝色曲线）到 30 天左右才能明显看出肿瘤生长的差异，而利用生物发光成像技术（绿色曲线）在细胞皮下注射 7 天即可观测到肿瘤的生长，并对随后的发展变化进行长期观测。



上图：A.利用萤火虫荧光素酶标记 MDA-MB-231 乳腺癌细胞建立原位乳腺癌模型，通过 IVIS 成像系统长期监测肿瘤的生长情况，在注射后 27 天观测到其他部位的转移信号。B.体外组织成像进一步验证转移的发生。



上图：应用 IVIS 成像系统进行肿瘤信号的 3D 光学成像（橙色），同时应用 Quantum FX uCT 成像系统进行 3D 解剖学成像（骨架），并将 3D 光学功能性结果与 3D 解剖学结果进行影像融合，而确定肿瘤的骨转移。

二.抗肿瘤药物研发

Revvity的小动物活体光学成像技术已广泛应用于肿瘤治疗药物的临床前研发阶段，发挥越来越重要的作用。全球各大制药企业均已采用活体光学成像技术开展抗肿瘤新药的研发，其中已有6种药物获得FDA认证，另有8种药物处于临床测试阶段。




FDA Approved Drugs	Drugs currently in clinical trials
<ul style="list-style-type: none">Sutent (Pfizer): subcutaneous tumor xenograftDasatinib (Bristol-Myers Squibb): chronic myelogenous leukemiaTasigna (Novartis): leukemia/metastasis modelCubicin (Cubist Pharmaceuticals): bacterial peritonitis modelAflibercept (Sanofi-Aventis): orthotopic renal cancerVelcade (Millennium Pharmaceuticals): multiple myeloma	<ul style="list-style-type: none">ABT-888 (Abbott): multiple diverse tumor modelsRANKL Inhibitor (Amgen): denosumab and bone metastasis modelsPanzem (EntreMed Pharmaceuticals): orthotopic gliosarcomaAEE788 (Novartis): intraperitoneal tumor modelIT-101/CRLX 101 (Insert Therapeutics, Cerulean Pharma Inc.): Ewings sarcomaCHIR-258 (Novartis): orthotopic multiple myeloma modelNPL-0052 (Nereus Pharmaceuticals): subcutaneous tumor modelCG0070 (Cell Genesys): orthotopic bladder cancer (*status of trial unknown)

应用小动物活体光学成像技术进行新药研发，主要包括以下方面：1、在活体动物水平进行药效评价；2、观测药物在活体动物体内的靶向、分布及代谢。


1、药效评价

利用荧光素酶标记肿瘤细胞，并移植入动物体内建立肿瘤疾病动物模型，应用小动物活体光学成像技术观测给药后肿瘤光学信号的变化情况，进而评价不同药物、特定的给药途径、时间、剂量等给药策略对于肿瘤的治疗效果，如下图：




药效评价研究: 利用萤火虫荧光素酶标记 MDA-MB-435 乳腺癌细胞株, 将细胞注入小鼠肾包膜下构建肾包膜肿瘤模型, 进而对 3 种不同药物或不同剂量的治疗效果进行评价。(A,C) 应用 IVIS 成像系统长期观测 3 种药物对肾包膜乳腺癌移植瘤的治疗效果, 并进行定量分析, 结果显示 30 mg/kg Docetaxel (TXT) 对肿瘤的生长抑制效果最好; (B,D) 应用 IVIS 成像系统长期观测 3 种药物对肾包膜乳腺癌移植瘤肺部转移的抑制效果, 并进行定量分析, 结果显示 30 mg/kg Docetaxel (TXT) 对肿瘤转移的抑制效果最好。

相对于触诊、肿瘤体积测量等传统方法, 利用高灵敏度的生物发光成像技术进行药物评价, 可以更灵敏的发现残余病灶点或尽早发现肿瘤的复发, 从而更准确的对药物治疗效果进行判定, 如下例:




提高药效评价准确性: 利用萤火虫荧光素酶标记 PC-3M 人前列腺肿瘤细胞株, 建立肿瘤皮下移植模型进行药物评价。左图: 对照组, 右图: 治疗组。通过高灵敏度的生物发光成像技术可以准确检测出药物治疗后的残余病灶点, 从而进行正确的药效评价, 而如果通过触诊等传统方法则可能做出错误判断。


利用生物发光成像技术进行药效评价的另一独特优势在于，可以明确判断药物是否有效杀死肿瘤活细胞。这是由于生物发光的原理是基于活细胞环境的酶促反应，因此，能够发光的细胞必定是具有活性的。下图所示为辉瑞公司抗肿瘤药物 **Sutent** 的部分研究结果，研究人员首先利用卡尺测量肿瘤体积的方法观测该药物对于肿瘤的生长抑制情况，发现该药物能够延缓肿瘤的生长，但体积测量数据显示肿瘤并未变小，研究人员随后利用生物发光成像技术进行观测，发现给药一定时间后肿瘤光学信号显著降低，说明该药物对肿瘤活性细胞确实具有杀伤作用，同时说明单独依靠肿瘤体积测量的方式无法准确真实反映药物治疗效果。凭借活体光学成像的实验结果，**Sutent** 得以顺利通过 FDA 的认证。



除了应用生物发光技术研究药物对于肿瘤的治疗效果之外，荧光成像技术同样可以应用于此类研究，主要方式为通过外源注射功能性荧光试剂观测药物对于肿瘤某一方面的治疗情况，如利用反映血管生成的荧光试剂观测药物对于肿瘤血管新生的抑制情况，又如利用反映细胞凋亡的荧光试剂观测由于药物治疗而诱发的肿瘤细胞凋亡情况。



上图：应用 FMT 荧光断层成像系统结合 AngioSense750 荧光试剂研究 CT322 药物在抑制肿瘤血管新生方面的作用效果，结果表明 CT322 可以通过抑制肿瘤血管新生进而抑制肿瘤的生长。




上图：应用 FMT 荧光断层成像系统结合 Annexin-Vivo 750 荧光试剂观测药物处理后引发的肿瘤细胞凋亡，定量结果显示应用荧光成像技术能够灵敏观测到药物诱导的细胞凋亡发生，而肿瘤体积测量数据显示对照组和治疗组的肿瘤体积无明显差异，因而无法反映药物诱导的细胞凋亡情况。

2、观测药物靶向、分布及代谢

除了在药效学中的应用，小动物活体光学成像也广泛应用于药物在体内靶向、分布及代谢的研究。与药效学中的应用不同，此类应用是以药物为直接观测对象，因此标记方式通常是利用荧光探针直接标记药物本身，通过追踪荧光信号而反映药物在体内的分布情况。


例如在研究抗体或多肽类药物是否能够有效靶向肿瘤的实验，可以利用荧光染料

通过化学键的结合标记目标抗体或多肽，经尾静脉注射后，利用小动物活体光学成像系统观测上述标记对象的肿瘤靶向性，如下图：



药物的肿瘤靶向性研究：利用 VivoTag 645 荧光染料标记抗癌药物曲妥珠单抗(Trastuzumab)，尾静脉注入携带 HER2 阳性人卵巢癌 SKOV3 的 SCID 小鼠体内，通过荧光成像观测不同时间点药物对肿瘤的靶向情况（左上图后三列），肿瘤本身已被荧光素酶标记而通过生物发光成像（左上图最左列），右上图为荧光定量分析结果，标明曲妥珠单抗对 HER2 阳性的人卵巢癌 SKOV3 具有良好靶向性。


事实上，直接通过荧光染料标记药物注入体内的方式存在诸多问题，如低靶向性、低药效性、免疫排斥、药物毒性问题等。因此，构建新型药物载体也是目前药物研究的热点之一，应用小动物活体光学成像技术同样可以在这一领域发挥作用。如下图所示为研究人员通过小动物活体光学成像技术观测一种新型纳米共聚物给药载体在体内运载吗啉基反义寡核苷酸靶向治疗肿瘤的实验。此类给药载体包含多种组分，如起整体保护作用的外层多聚苹果酸骨架及聚乙二醇、起肿瘤靶向作用的抗体连接组分、药物释放组分、荧光染料结合位点等。结果表明，通过应用新型纳米给药载体可以有效提升药物的靶向治疗效果。



上图左：新型纳米共聚物给药载体结构示意图；上图右：应用 IVIS 成像系统观测不同时间点药物运载系统在体内的分布及对肿瘤的靶向。


在研究药物靶向的同时，了解不同时间点药物在动物体各个器官的分布及最终代谢情况同样必不可少。Revvity 的 FMT 小动物活体荧光断层成像系统可以很好的满足

此类应用需求。应用 FMT 成像系统，能够对荧光标记的药物在深层器官的分布进行断层扫描及三维重建，获得真实准确的三维定量数据，进而对药物的体内分布代谢情况作出正确分析。如下图：



上图左：利用 FMT 成像系统观测不同时间点经荧光染料 VivioTag 680 标记的 BSA 在小鼠不同器官的分布；上图右：不同时间点不同器官 BSA 分布的定量结果。

应用活体荧光成像技术进行药物分布的观测目前主要局限于对生物大分子药物的研究，而对天然或化学小分子药物的分布代谢研究主要依靠的是放射性核素标记成像（PET 或 SPECT），原因是用相对大分子量的荧光染料标记小分子药物，则荧光染料本身即会对小分子药物在体内的分布代谢产生影响，因此无法将活体荧光成像技术应用于此类研究。然而，科研人员最近发现应用高灵敏度的 IVIS 活体光学成像系统可以观测到放射性核素在小鼠体内发出的光学信号，其所依据的是切伦科夫辐射（Cherenkov Radiation）原理，由此拓展了活体光学成像技术的应用范畴，即利用活体光学成像系统观测经放射性核素标记的小分子药物在体内的分布代谢。



切伦科夫辐射成像：应用 IVIS 成像系统结合 $[^{18}\text{F}]\text{FDG}$ 放射性探针观测肿瘤。A.肿瘤生物发光成像结果；B.应用 IVIS 系统进行切伦科夫辐射成像，观测不同时间点 $[^{18}\text{F}]\text{FDG}$ 在肿瘤部位的代谢；C.应用 PET 系统成像，观测不同时间点 $[^{18}\text{F}]\text{FDG}$ 在肿瘤部位的代谢，结果与 B 一致。


三.癌症分子机理研究

对于癌症分子机理的研究之前一直局限于体外水平，体外研究的缺陷在于无法模拟肿瘤在动物体内真实的生理微环境，因此，单一的体外研究结果并不能完全反映癌症的发生发展机理。小动物活体光学成像技术使科研人员能够进一步将癌症分子机理的研究由体外拓展至体内，如在活体动物水平研究癌症相关基因在癌症发生发展进程中的作用、观测肿瘤发生发展过程中特异性分子事件的发生等。

1、应用生物发光技术研究癌症相关基因的作用


应用生物发光技术研究癌症相关基因的研究方法，主要利用荧光素酶标记特定基因，构建特定基因-荧光素酶的共表达载体，通过荧光素酶产生的生物发光信号反映该基因的表达情况，研究该基因的相关作用。以下几个例子是应用生物发光技术研究癌症相关基因在不同肿瘤模型中作用：

p53 是调节细胞正常生命活动的一种重要基因，控制着细胞周期的启动。p53 也同时被认为是一种重要的抑癌基因，在人类 50% 以上的肿瘤组织中均发现了 p53 基因的突变，这是肿瘤中最常见的遗传学改变，说明该基因的改变很可能是人类肿瘤产生的主要发病因素。下图所示是发表于 2007 年 Nature 杂志上的一篇研究 p53 抑癌作用的文献结果。研究者将一个同时带有致癌基因 ras、四环素反式激活因子 tTA ('tet-off')及四环素依赖性 p53 shRNA 的逆转录病毒表达载体与荧光素酶基因质粒共转染从小鼠胚胎提取的成肝细胞，并将转染细胞接种于小鼠肝脏，观测 p53 基因表达关闭或开启对于肝癌发生或消亡的作用。结果显示，当 p53 的表达被 shRNA 抑制时，通过生物发光观测到的肿瘤信号逐渐增强，而当给小鼠注入强力霉素 (doxycycline) 开启 p53 的表达后，肝癌细胞的生物发光信号逐渐减弱，说明 p53 的表达能够抑制肿瘤的生长并促使肿瘤消亡。




p53 基因的抑癌作用：上图 a，实验设计示意图；上图 b，关闭-开启 p53 的表达对肿瘤发展的影响；上图 c，关闭-开启-再关闭 p53 的表达对肿瘤发展的影响。

乳腺癌是女性最常见的恶性肿瘤之一，而乳腺癌转移是导致患者死亡的主要原因，乳腺癌转移的分子机理目前还不完全清楚。LSD1 是第一个被发现的组蛋白去甲基化酶，理论上对基因转录起广泛调控作用，但近期研究却表明 LSD1 只参与一些特异的细胞信号通路的调控而且与多种肿瘤的发生发展高度相关。北大医学部尚永丰教授课题组应用 IVIS 活体光学成像系统观测到 LSD1 能抑制乳腺癌的侵袭和转移，从而揭示了 LSD1 这一表观调控因子在抑制乳腺癌转移中的重要作用，并为乳腺癌转移的干预提供了新的分子靶点。上述结果发表于 2009 年 Cell 杂志。



LSD1 对乳腺癌转移的抑制作用: 将带有 LSD1 或 LSD1 siRNA 或相应对照的表达载体通过慢病毒转染经荧光素酶标记的人乳腺癌细胞株 MDA-MB-231-Luc-D3H2LN, 并原位接种于小鼠腹部乳腺脂肪垫, 观测 LSD1 的表达与抑制对乳腺癌细胞肺部转移的影响。上图左: 应用 IVIS 系统进行生物发光成像结果; 上图右: 光学定量结果。

肿瘤细胞生长、转移和复发的特点与干细胞的基本特性十分相似，因此，有学者提出肿瘤干细胞 (Cancer/Tumor Stem Cell, CSC/TSC) 的理论。肿瘤干细胞可以长时间处于休眠状态并具有多种耐药分子而对杀伤肿瘤细胞的外界理化因素不敏感，导致肿瘤往往在常规治疗方法消灭大部分普通肿瘤细胞后一段时间复发。因此，如果能够寻找到肿瘤干细胞特异性治疗靶标，将为肿瘤治疗开辟新的路径。OCT4 是参与调控胚胎干细胞自我更新和维持其全能性的最为重要的转录因子之一，因此，OCT4 可能是杀伤肿瘤干细胞的潜在靶点之一。下述实验是应用 IVIS 系统观测 OCT4 在乳腺癌细胞成瘤中的作用。研究者将带有 OCT4 基因的载体转染经荧光素酶标记的鼠源乳腺癌细胞株 4T1-luc, 并分选出高表达及低表达 OCT4 基因的细胞株，分别原位接种于小鼠胸部乳腺脂肪垫，观测结果显示 OCT4 高表达乳腺癌细胞株具有更高的成瘤性。




2、应用功能性荧光试剂观测肿瘤内部特异性分子事件的发生

肿瘤发展进程中伴随着诸多分子事件的发生，如某些蛋白酶的表达特异性升高、某些表面标识物的特异性表达、肿瘤周边血管的新生、肿瘤组织局部缺氧等，通过观测这些分子事件能够判断肿瘤的发展程度并作出预后。Revvity提供了一系列应用于肿瘤研究的功能性荧光试剂，使上述观测成为可能。以下为具体应用实例：

应用 PSA 750 FAST 酶激活类荧光试剂观测 PSA⁺前列腺肿瘤


前列腺特异抗原(Prostate Specific Antigen, PSA) 主要是由前列腺上皮细胞产生的蛋白分解酶，正常情况下被分泌入前列腺液或精液中以有活性的游离形式存在，血清中的 PSA 主要以结合形式存在。正常及良性前列腺增生的前列腺上皮均可分泌游离或结合形式的 PSA，但具有酶活性的 PSA 只存在于恶性前列腺肿瘤中。PSA 750 FAST 是一种只能被具有酶活性 PSA 特异性激活的荧光试剂，因此可用于监测前列腺肿瘤的恶性程度。



应用 IVIS Spectrum 及 FMT 成像系统结合 PSA 750 FAST 荧光试剂对裸鼠皮下接种的人前列腺肿瘤 PSA⁺ LNCaP 进行活体光学成像，成像时间为试剂尾静脉注射后 6h。


应用 IntegriSense 及 ProSense 荧光试剂观测肿瘤血管新生及组织蛋白酶的表达

肿瘤的发生发展伴随着诸多分子事件的共同发生，应用不同种类的荧光试剂，可以实现在一个实验中观测肿瘤内部的多个生物学进程。如下图所示，在一个实验中同时利用 IntegriSense 及 ProSense 两种荧光试剂对肿瘤进行观测：应用 IntegriSense750 靶向类荧光试剂可以靶向监测肿瘤血管上皮特异性表达的整联蛋白 $\alpha v\beta 3$ 而反映肿瘤的血管新生；而应用 ProSense680 酶激活类荧光试剂能够监测肿瘤细胞中组织蛋白酶的活性，进而揭示两种对象在肿瘤内部的不同分布及药物 (Avastin) 治疗后的不同变化。




应用 HypoxiSense 及 AngioSense 荧光试剂观测肿瘤微环境

肿瘤的发生和转移与肿瘤细胞所处的内外环境有着密切关系，了解肿瘤微环境对于肿瘤的诊断、防治和预后有着重要意义。下图所示为利用监测肿瘤组织缺氧的 HypoxiSense 靶向类荧光试剂与监测肿瘤血管生成的 AngioSense 血管生理类荧光试剂共同观测小鼠皮下接种的人宫颈癌肿瘤微环境。其中，HypoxiSense 能够特异性靶向缺氧肿瘤细胞表面上调表达的碳酸酐酶 9 (CAIX)，进而表征肿瘤的缺氧区域；而 AngioSense 通过富集于由于肿瘤血管新生而引发的血管渗漏区域，进而表征肿瘤的血管富集区域。定量结果显示缺氧部位主要位于肿瘤内部中心区。



应用 FMT 成像系统结合 HypoxiSense680 及 AngioSense750 荧光试剂观测裸鼠胸部脂肪垫接种的人宫颈癌肿瘤微环境，并进行定量分析。



IVIS APPLICATIONS IN ONCOLOGY RESEARCH

1. Abrams TJ, Murray LJ, Pesenti E, Holway VW, Colombo T, Lee LB, Cherrington JM, Pryer NK. Preclinical evaluation of the tyrosine kinase inhibitor SU11248 as a single agent and in combination with "standard of care" therapeutic agents for the treatment of breast cancer. *Molecular Cancer Therapy*. October 2003;2(10):1011-1021.
2. Acuff HB, Carter KJ, Fingleton B, Gorden DL, Matrisian LM. Matrix Metalloproteinase-9 from Bone Marrow-Derived Cells Contributes to Survival but not Growth of Tumor Cells in the Lung Microenvironment. *Cancer Res*. January 1, 2006 2006;66(1):259-266.
3. Adams JY, Johnson M, Sato M, Berger F, Gambhir SS, Carey M, Iruela-Arispe ML, Wu L. Visualization of advanced human prostate cancer lesions in living mice by a targeted gene transfer vector and optical imaging. *Nature medicine*. August 2002;8(8):891-897.
4. Agrawal VK, Copeland KM, Barbachano Y, Rahim A, Seth R, White CL, Hingorani M, Nutting CM, Kelly M, Harris P, Pandha H, Melcher AA, Vile RG, Porter C, Harrington KJ. Microvascular free tissue transfer for gene delivery: in vivo evaluation of different routes of plasmid and adenoviral delivery. *Gene therapy*. Jan 2009;16(1):78-92.
5. Ahmed N, Ratnayake M, Savoldo B, Perlaky L, Dotti G, Wels WS, Bhattacharjee MB, Gilbertson RJ, Shine HD, Weiss HL, Rooney CM, Heslop HE, Gottschalk S. Regression of experimental medulloblastoma following transfer of HER2-specific T cells. *Cancer Res*. Jun 15 2007;67(12):5957-5964.
6. Ahmed N, Salsman VS, Kew Y, Shaffer D, Powell S, Zhang YJ, Grossman RG, Heslop HE, Gottschalk S. HER2-specific T cells target primary glioblastoma stem cells and induce regression of autologous experimental tumors. *Clin Cancer Res*. Jan 15 2010;16(2):474-485.
7. Ahmed N, Salsman VS, Yvon E, Louis CU, Perlaky L, Wels WS, Dishop MK, Kleinerman EE, Pule M, Rooney CM, Heslop HE, Gottschalk S. Immunotherapy for osteosarcoma: genetic modification of T cells overcomes low levels of tumor antigen expression. *Mol Ther*. Oct 2009;17(10):1779-1787.
8. Akens MK, Yee AJ, Wilson BC, Burch S, Johnson CL, Lilge L, Bisland SK. Photodynamic therapy of vertebral metastases: evaluating tumor-to-neural tissue uptake of BPD-MA and ALA-PpIX in a murine model of metastatic human breast carcinoma. *Photochemistry and photobiology*. Sep-Oct 2007;83(5):1034-1039.
9. Alain T, Lun X, Martineau Y, Sean P, Pulendran B, Petroulakis E, Zemp FJ, Lemay CG, Roy D, Bell JC, Thomas G, Kozma SC, Forsyth PA, Costa-Mattioli M, Sonenberg N. Vesicular stomatitis virus oncolysis is potentiated by impairing mTORC1-dependent type I IFN production. *Proceedings of the National Academy of Sciences of the United States of America*. Jan 26 2010;107(4):1576-1581.
10. Alphas HH, Gambhira R, Karanam B, Roberts JN, Jagu S, Schiller JT, Zeng W, Jackson DC, Roden RB. Protection against heterologous human papillomavirus challenge by a synthetic lipopeptide vaccine containing a broadly

- cross-neutralizing epitope of L2. *Proceedings of the National Academy of Sciences of the United States of America*. Apr 15 2008;105(15):5850-5855.
11. Antczak C, Kloepping C, Radu C, Genski T, Muller-Kuhr L, Siems K, de Stanchina E, Abramson DH, Djaballah H. Revisiting Old Drugs as Novel Agents for Retinoblastoma: In vitro and In vivo Antitumor Activity of Cardenolides. *Investigative ophthalmology & visual science*. Jan 17 2009.
 12. Aoi A, Watanabe Y, Mori S, Takahashi M, Vassaux G, Kodama T. Herpes Simplex Virus Thymidine Kinase-Mediated Suicide Gene Therapy Using Nano/Microbubbles and Ultrasound. *Ultrasound in medicine & biology*. Dec 18 2007.
 13. Armstrong SA, Kung AL, Mabon ME, Silverman LB, Stam RW, Den Boer ML, Pieters R, Kersey JH, Sallan SE, Fletcher JA, Golub TR, Griffin JD, Korsmeyer SJ. Inhibition of FLT3 in MLL. Validation of a therapeutic target identified by gene expression based classification. *Cancer cell*. February 2003;3(2):173-183.
 14. Arnold RS, Sun CQ, Richards JC, Grigoriev G, Coleman IM, Nelson PS, Hsieh CL, Lee JK, Xu Z, Rogatko A, Osunkoya AO, Zayzafoon M, Chung L, Petros JA. Mitochondrial DNA mutation stimulates prostate cancer growth in bone stromal environment. *The Prostate*. Jan 1 2009;69(1):1-11.
 15. Arumugam T, Simeone DM, Van Golen K, Logsdon CD. S100P Promotes Pancreatic Cancer Growth, Survival, and Invasion. *Clin Cancer Res*. August 1, 2005 2005;11(15):5356-5364.
 16. Aung W, Hasegawa S, Koshikawa-Yano M, Obata T, Ikehira H, Furukawa T, Aoki I, Saga T. Visualization of in vivo electroporation-mediated transgene expression in experimental tumors by optical and magnetic resonance imaging. *Gene therapy*. Jul 2009;16(7):830-839.
 17. Aung W, Hasegawa S, Koshikawa-Yano M, Tsuji AB, Sogawa C, Sudo H, Sugyo A, Koizumi M, Furukawa T, Saga T. Noninvasive assessment of regulable transferred-p53 gene expression and evaluation of therapeutic response with FDG-PET in tumor model. *Gene therapy*. May 6 2010.
 18. Azab AK, Runnels JM, Pitsillides C, Moreau AS, Azab F, Leleu X, Jia X, Wright R, Ospina B, Carlson AL, Alt C, Burwick N, Roccaro AM, Ngo HT, Farag M, Melhem MR, Sacco A, Munshi NC, Hideshima T, Rollins BJ, Anderson KC, Kung AL, Lin CP, Ghobrial IM. CXCR4 inhibitor AMD3100 disrupts the interaction of multiple myeloma cells with the bone marrow microenvironment and enhances their sensitivity to therapy. *Blood*. Apr 30 2009;113(18):4341-4351.
 19. Bai X, Yan Y, Song YH, Seidensticker M, Rabinovich B, Metzeler R, Bankson JA, Vykoukal D, Alt E. Both cultured and freshly isolated adipose tissue-derived stem cells enhance cardiac function after acute myocardial infarction. *European heart journal*. Feb 2010;31(4):489-501.
 20. Baia GS, Dinca EB, Ozawa T, Kimura ET, McDermott MW, James CD, Vandenberg SR, Lal A. An Orthotopic Skull Base Model of Malignant Meningioma. *Brain Pathol*. Dec 17 2007.
 21. Bailey CL, Kelly P, Casey PJ. Activation of Rap1 promotes prostate cancer metastasis. *Cancer Res*. Jun 15 2009;69(12):4962-4968.

22. Balani P, Boulaire J, Zhao Y, Zeng J, Lin J, Wang S. High mobility group box2 promoter-controlled suicide gene expression enables targeted glioblastoma treatment. *Mol Ther*. Jun 2009;17(6):1003-1011.
23. Banaszynski LA, Sellmyer MA, Contag CH, Wandless TJ, Thorne SH. Chemical control of protein stability and function in living mice. *Nature medicine*. Oct 2008;14(10):1123-1127.
24. Bao L, Gorin MA, Zhang M, Ventura AC, Pomerantz WC, Merajver SD, Teknos TN, Mapp AK, Pan Q. Preclinical development of a bifunctional cancer cell homing, PKCepsilon inhibitory peptide for the treatment of head and neck cancer. *Cancer Res*. Jul 15 2009;69(14):5829-5834.
25. Bartlett DW, Davis ME. Insights into the kinetics of siRNA-mediated gene silencing from live-cell and live-animal bioluminescent imaging. *Nucl. Acids Res*. January 12, 2006 2006;34(1):322-333.
26. Basseres DS, Ebbs A, Levantini E, Baldwin AS. Requirement of the NF-kappaB subunit p65/RelA for K-Ras-induced lung tumorigenesis. *Cancer Res*. May 1 2010;70(9):3537-3546.
27. Bastianutto C, Mian A, Symes J, Mocanu J, Alajez N, Sleep G, Shi W, Keating A, Crump M, Gospodarowicz M, Medin J, Minden M, Liu FF. Local radiotherapy induces homing of hematopoietic stem cells to the irradiated bone marrow. *Cancer Res*. Nov 1 2007;67(21):10112-10116.
28. Becher OJ, Hambardzumyan D, Fomchenko EI, Momota H, Mainwaring L, Bleau AM, Katz AM, Edgar M, Kenney AM, Cordon-Cardo C, Blasberg RG, Holland EC. Gli activity correlates with tumor grade in platelet-derived growth factor-induced gliomas. *Cancer Res*. Apr 1 2008;68(7):2241-2249.
29. Beck BH, Kim HG, Kim H, Samuel S, Liu Z, Shrestha R, Haines H, Zinn K, Lopez RD. Adoptively transferred ex vivo expanded gammadelta-T cells mediate in vivo antitumor activity in preclinical mouse models of breast cancer. *Breast cancer research and treatment*. Sep 18 2009.
30. Becker CM, Wright RD, Satchi-Fainaro R, Funakoshi T, Folkman J, Kung AL, D'Amato RJ. A Novel Noninvasive Model of Endometriosis for Monitoring the Efficacy of Antiangiogenic Therapy. *Am J Pathol*. June 1, 2006 2006;168(6):2074-2084.
31. Bertino P, Piccardi F, Porta C, Favoni R, Cilli M, Mutti L, Gaudino G. Imatinib mesylate enhances therapeutic effects of gemcitabine in human malignant mesothelioma xenografts. *Clin Cancer Res*. Jan 15 2008;14(2):541-548.
32. Bhaumik S, Gambhir SS. Optical imaging of Renilla luciferase reporter gene expression in living mice. *Proceedings of the National Academy of Sciences of the United States of America*. January 2002;99(1):377-382.
33. Bjorndahl MA, Cao R, Burton JB, Brakenhielm E, Religa P, Galter D, Wu L, Cao Y. Vascular Endothelial Growth Factor-A Promotes Peritumoral Lymphangiogenesis and Lymphatic Metastasis. *Cancer Res*. October 15, 2005 2005;65(20):9261-9268.
34. Black KL, Yin D, Ong JM, Hu J, Konda BM, Wang X, Ko MK, Bayan JA, Sacapano MR, Espinoza A, Irvin DK, Shu Y. PDE5 inhibitors enhance tumor permeability and

- efficacy of chemotherapy in a rat brain tumor model. *Brain research*. Sep 16 2008;1230:290-302.
35. Blankenberg FG, Mandl S, Cao YA, O'Connell-Rodwell C, Contag C, Mari C, Gaynutdinov TI, Vanderheyden JL, Backer MV, Backer JM. Tumor imaging using a standardized radiolabeled adapter protein docked to vascular endothelial growth factor. *The Journal of Nuclear Medicine*. August 2004;45(8):1373-1380.
 36. Blasberg RG. In vivo molecular-genetic imaging: multi-modality nuclear and optical combinations. *Nuclear Medicine and Biology*. November 2003;30(8):879-888.
 37. Blum G, von Degenfeld G, Merchant MJ, Blau HM, Bogyo M. Noninvasive optical imaging of cysteine protease activity using fluorescently quenched activity-based probes. *Nat Chem Biol*. Oct 2007;3(10):668-677.
 38. Bondareva A, Downey CM, Ayres F, Liu W, Boyd SK, Hallgrimsson B, Jirik FR. The lysyl oxidase inhibitor, beta-aminopropionitrile, diminishes the metastatic colonization potential of circulating breast cancer cells. *PLoS ONE*. 2009;4(5):e5620.
 39. Bortolanza S, Bunuales M, Otano I, Gonzalez-Asequinolaza G, Ortiz-de-Solorzano C, Perez D, Prieto J, Hernandez-Alcoceba R. Treatment of pancreatic cancer with an oncolytic adenovirus expressing interleukin-12 in Syrian hamsters. *Mol Ther*. Apr 2009;17(4):614-622.
 40. Bos PD, Zhang XH, Nadal C, Shu W, Gomis RR, Nguyen DX, Minn AJ, van de Vijver MJ, Gerald WL, Foekens JA, Massague J. Genes that mediate breast cancer metastasis to the brain. *Nature*. Jun 18 2009;459(7249):1005-1009.
 41. Bradbury MS, Tomishima M, Panagiotakos G, Chan BK, Zanzonico P, Vider J, Ponomarev V, Studer L, Tabar V. Optical bioluminescence imaging of human ES cell progeny in the rodent CNS. *J Neurochem*. Jun 7 2007.
 42. Brader P, Riedl CC, Woo Y, Ponomarev V, Zanzonico P, Wen B, Cai S, Hricak H, Fong Y, Blasberg R, Serganova I. Imaging of hypoxia-driven gene expression in an orthotopic liver tumor model. *Molecular cancer therapeutics*. Nov 2007;6(11):2900-2908.
 43. Brader P, Stritzker J, Riedl CC, Zanzonico P, Cai S, Burnazi EM, Ghani ER, Hricak H, Szalay AA, Fong Y, Blasberg R. Escherichia coli Nissle 1917 Facilitates Tumor Detection by Positron Emission Tomography and Optical Imaging. *Clin Cancer Res*. Mar 27 2008.
 44. Breitbach CJ, Paterson JM, Lemay CG, Falls TJ, McGuire A, Parato KA, Stojdl DF, Daneshmand M, Speth K, Kirn D, McCart JA, Atkins H, Bell JC. Targeted inflammation during oncolytic virus therapy severely compromises tumor blood flow. *Mol Ther*. Sep 2007;15(9):1686-1693.
 45. Bremer E, van Dam GM, de Bruyn M, van Riezen M, Dijkstra M, Kamps G, Helfrich W, Haisma H. Potent systemic anticancer activity of adenovirally expressed EGFR-selective TRAIL fusion protein. *Mol Ther*. Dec 2008;16(12):1919-1926.
 46. Brentjens RJ, Santos E, Nikhamin Y, Yeh R, Matsushita M, La Perle K, Quintas-Cardama A, Larson SM, Sadelain M. Genetically targeted T cells

-
- eradicate systemic acute lymphoblastic leukemia xenografts. *Clin Cancer Res*. Sep 15 2007;13(18 Pt 1):5426-5435.
47. Briat A, Vassaux G. A new transgenic mouse line to image chemically induced p53 activation in vivo. *Cancer science*. Apr 2008;99(4):683-688.
 48. Brody JD, Goldstein MJ, Czerwinski DK, Levy R. Immunotransplantation preferentially expands T-effector cells over T-regulatory cells and cures large lymphoma tumors. *Blood*. Jan 1 2009;113(1):85-94.
 49. Bryant MJ, Chuah TL, Luff J, Lavin MF, Walker DG. A novel rat model for glioblastoma multiforme using a bioluminescent F98 cell line. *J Clin Neurosci*. May 2008;15(5):545-551.
 50. Burgos JS, Rosol M, Moats RA, Khankaldyyan V, Kohn DB, Nelson MD, Jr., Laug WE. Time course of bioluminescent signal in orthotopic and heterotopic brain tumors in nude mice. *BioTechniques*. June 2003;34(6):1184-1188.
 51. Burton JB, Johnson M, Sato M, Koh SB, Mulholland DJ, Stout D, Chatziioannou AF, Phelps ME, Wu H, Wu L. Adenovirus-mediated gene expression imaging to directly detect sentinel lymph node metastasis of prostate cancer. *Nature medicine*. Aug 2008;14(8):882-888.
 52. Butterworth KT, McCarthy HO, Devlin A, Ming L, Robson T, McKeown SR, Worthington J. Hypoxia selects for androgen independent LNCaP cells with a more malignant geno- and phenotype. *Int J Cancer*. Aug 15 2008;123(4):760-768.
 53. Cai W, Chen X. Multimodality imaging of vascular endothelial growth factor and vascular endothelial growth factor receptor expression. *Front Biosci*. 2007;12:4267-4279.
 54. Cai W, Chen X. Preparation of peptide-conjugated quantum dots for tumor vasculature-targeted imaging. *Nature protocols*. 2008;3(1):89-96.
 55. Calabrese C, Poppleton H, Kocak M, Hogg TL, Fuller C, Hamner B, Oh EY, Gaber MW, Finklestein D, Allen M, Frank A, Bayazitov IT, Zakharenko SS, Gajjar A, Davidoff A, Gilbertson RJ. A perivascular niche for brain tumor stem cells. *Cancer cell*. Jan 2007;11(1):69-82.
 56. Canon JR, Roudier M, Bryant R, Morony S, Stolina M, Kostenuik PJ, Dougall WC. Inhibition of RANKL blocks skeletal tumor progression and improves survival in a mouse model of breast cancer bone metastasis. *Clin Exp Metastasis*. Dec 5 2007.
 57. Casarez EV, Dunlap-Brown ME, Conaway MR, Amorino GP. Radiosensitization and modulation of p44/42 mitogen-activated protein kinase by 2-Methoxyestradiol in prostate cancer models. *Cancer Res*. Sep 1 2007;67(17):8316-8324.
 58. Casey G, Cashman JP, Morrissey D, Whelan MC, Larkin JO, Soden DM, Tangney M, O'Sullivan GC. Sonoporation mediated immunogene therapy of solid tumors. *Ultrasound in medicine & biology*. Mar 2010;36(3):430-440.
 59. Cawood R, Chen HH, Carroll F, Bazan-Peregrino M, van Rooijen N, Seymour LW. Use of tissue-specific microRNA to control pathology of wild-type adenovirus without attenuation of its ability to kill cancer cells. *PLoS pathogens*. May 2009;5(5):e1000440.
 60. Chan CT, Paulmurugan R, Gheysens OS, Kim J, Chiosis G, Gambhir SS. Molecular imaging of the efficacy of heat shock protein 90 inhibitors in living

-
- subjects. *Cancer Res.* Jan 1 2008;68(1):216-226.
61. Chan CT, Paulmurugan R, Reeves RE, Solow-Cordero D, Gambhir SS. Molecular imaging of phosphorylation events for drug development. *Mol Imaging Biol.* May-Jun 2009;11(3):144-158.
 62. Chan JK, Hamilton CA, Cheung MK, Karimi M, Baker J, Gall JM, Schulz S, Thorne SH, Teng NN, Contag CH, Lum LG, Negrin RS. Enhanced killing of primary ovarian cancer by retargeting autologous cytokine-induced killer cells with bispecific antibodies: a preclinical study. *Clin Cancer Res.* Mar 15 2006;12(6):1859-1867.
 63. Chanda D, Isayeva T, Kumar S, Hensel JA, Sawant A, Ramaswamy G, Siegal GP, Beatty MS, Ponnazhagan S. Therapeutic potential of adult bone marrow-derived mesenchymal stem cells in prostate cancer bone metastasis. *Clin Cancer Res.* Dec 1 2009;15(23):7175-7185.
 64. Chanda D, Isayeva T, Kumar S, Siegal GP, Szafran AA, Zinn KR, Reddy VV, Ponnazhagan S. Systemic osteoprotegerin gene therapy restores tumor-induced bone loss in a therapeutic model of breast cancer bone metastasis. *Mol Ther.* May 2008;16(5):871-878.
 65. Chandna P, Saad M, Wang Y, Ber E, Khandare J, Vetcher AA, Soldatenkov VA, Minko T. Targeted Proapoptotic Anticancer Drug Delivery System. *Molecular pharmaceuticals.* Oct 1 2007;4(5):668-678.
 66. Charafe-Jauffret E, Ginestier C, Iovino F, Tarpin C, Diebel M, Esterni B, Houvenaeghel G, Extra JM, Bertucci F, Jacquemier J, Xerri L, Dontu G, Stassi G, Xiao Y, Barsky SH, Birnbaum D, Viens P, Wicha MS. Aldehyde dehydrogenase 1-positive cancer stem cells mediate metastasis and poor clinical outcome in inflammatory breast cancer. *Clin Cancer Res.* Jan 1 2010;16(1):45-55.
 67. Charafe-Jauffret E, Ginestier C, Iovino F, Wicinski J, Cervera N, Finetti P, Hur MH, Diebel ME, Monville F, Dutcher J, Brown M, Viens P, Xerri L, Bertucci F, Stassi G, Dontu G, Birnbaum D, Wicha MS. Breast cancer cell lines contain functional cancer stem cells with metastatic capacity and a distinct molecular signature. *Cancer Res.* Feb 15 2009;69(4):1302-1313.
 68. Chen L, Yang S, Jakoncic J, Zhang JJ, Huang XY. Migrastatin analogues target fascin to block tumour metastasis. *Nature.* Apr 15 2010;464(7291):1062-1066.
 69. Chen X, Conti PS, Moats RA. In vivo near-infrared fluorescence imaging of integrin $\alpha v \beta 3$ in brain tumor xenografts. *Cancer research.* November 2004;64(21):8009-8014.
 70. Cheng CM, Lu YL, Chuang KH, Hung WC, Shiea J, Su YC, Kao CH, Chen BM, Roffler S, Cheng TL. Tumor-targeting prodrug-activating bacteria for cancer therapy. *Cancer gene therapy.* Jun 2008;15(6):393-401.
 71. Cheng JC, Kinjo K, Judelson D, Chang J, Wu WS, Schmid I, Shankar DB, Kasahara N, Stripecke R, Bhatia R, Landaw EM, Sakamoto KM. CREB is a critical regulator of normal hematopoiesis and leukemogenesis. *Blood.* Nov 1 2007.
 72. Cheng Z, Levi J, Xiong Z, Gheysens O, Keren S, Chen X, Gambhir SS. Near-Infrared Fluorescent Deoxyglucose Analogue for Tumor Optical Imaging in Cell Culture and Living Mice. *Bioconjugate Chem.* May 17, 2006

-
- 2006;17(3):662-669.
73. Cheng Z, Wu Y, Xiong Z, Gambhir SS, Chen X. Near-infrared fluorescent RGD peptides for optical imaging of integrin alphavbeta3 expression in living mice. *Bioconjugate chemistry*. Nov-Dec 2005;16(6):1433-1441.
 74. Chesler L, Goldenberg DD, Collins R, Grimmer M, Kim GE, Tihan T, Nguyen K, Yakovenko S, Matthay KK, Weiss WA. Chemotherapy-induced apoptosis in a transgenic model of neuroblastoma proceeds through p53 induction. *Neoplasia (New York, N.Y.)*. Nov 2008;10(11):1268-1274.
 75. Chiu CC, Lin CY, Lee LY, Chen YJ, Kuo TF, Chang JT, Liao CT, Wang HM, Yen TC, Shen CR, Liao SK, Cheng AJ. Glucose-regulated protein 78 regulates multiple malignant phenotypes in head and neck cancer and may serve as a molecular target of therapeutic intervention. *Molecular cancer therapeutics*. Sep 2008;7(9):2788-2797.
 76. Choi CY, Chan DA, Paulmurugan R, Sutphin PD, Le QT, Koong AC, Zundel W, Gambhir SS, Giaccia AJ. Molecular imaging of hypoxia-inducible factor 1 alpha and von Hippel-Lindau interaction in mice. *Mol Imaging*. May-Jun 2008;7(3):139-146.
 77. Choi YS, Hoory T, Monie A, Wu A, Connolly D, Hung CF. alpha-Galactosylceramide enhances the protective and therapeutic effects of tumor cell based vaccines for ovarian tumors. *Vaccine*. Oct 29 2008;26(46):5855-5863.
 78. Choy G, Choyke P, Libutti SK. Current advances in molecular imaging: noninvasive in vivo bioluminescent and fluorescent optical imaging in cancer research. *Molecular Imaging*. October 2003;2(4):303-312.
 79. Christian DA, Cai S, Garbuzenko OB, Harada T, Zajac AL, Minko T, Discher DE. Flexible Filaments for in Vivo Imaging and Delivery: Persistent Circulation of Filomicelles Opens the Dosage Window for Sustained Tumor Shrinkage. *Molecular pharmaceutics*. Mar 25 2009.
 80. Chu K, Cheng CJ, Ye X, Lee YC, Zurita AJ, Chen DT, Yu-Lee LY, Zhang S, Yeh ET, Hu MC, Logothetis CJ, Lin SH. Cadherin-11 promotes the metastasis of prostate cancer cells to bone. *Mol Cancer Res*. Aug 2008;6(8):1259-1267.
 81. Chuang CM, Monie A, Wu A, Mao CP, Hung CF. Treatment with LL-37 peptide enhances antitumor effects induced by CpG oligodeoxynucleotides against ovarian cancer. *Human gene therapy*. Apr 2009;20(4):303-313.
 82. Chuang CM, Monie A, Wu A, Pai SI, Hung CF. Combination of viral oncolysis and tumor-specific immunity to control established tumors. *Clin Cancer Res*. Jul 15 2009;15(14):4581-4588.
 83. Chung E, Yamashita H, Au P, Tannous BA, Fukumura D, Jain RK. Secreted Gaussia luciferase as a biomarker for monitoring tumor progression and treatment response of systemic metastases. *PLoS ONE*. 2009;4(12):e8316.
 84. Coffelt SB, Marini FC, Watson K, Zvezdaryk KJ, Dembinski JL, LaMarca HL, Tomchuck SL, Honer zu Bentrup K, Danka ES, Henkle SL, Scandurro AB. The pro-inflammatory peptide LL-37 promotes ovarian tumor progression through recruitment of multipotent mesenchymal stromal cells. *Proceedings of the*

-
- National Academy of Sciences of the United States of America*. Mar 10 2009;106(10):3806-3811.
85. Collisson EA, Kleer C, Wu M, De A, Gambhir SS, Merajver SD, Kolodney MS. Atorvastatin prevents RhoC isoprenylation, invasion, and metastasis in human melanoma cells. *Molecular Cancer Therapy*. October 2003;2(10):941-948.
 86. Contag CH, Jenkins D, Contag PR, Negrin RS. Use of reporter genes for optical measurements of neoplastic disease in vivo. *Neoplasia (New York, N.Y.)* Jan-April 2000;2(1-2):41-52.
 87. Contag CH, Ross BD. It's not just about anatomy: in vivo bioluminescence imaging as an eyepiece into biology. *Journal of Magnetic Resonance Imaging*. October 2002;16(4):378-387.
 88. Contag PR. Unraveling the complexity of oncogenesis through in vivo optical imaging. In: LaRochelle WJ, Shimkets RA, eds. *Cancer Drug Discovery and Development: the Oncogenomics Handbook*. Vol Chapter 17. Totowa, NJ: Humana Press, Inc.; 2005:263-274.
 89. Contero A, Richer E, Gondim A, Mason RP. High-throughput quantitative bioluminescence imaging for assessing tumor burden. *Methods in molecular biology (Clifton, N.J.)* 2009;574:37-45.
 90. Conway AE, Lindgren A, Galic Z, Pyle AD, Wu H, Zack JA, Pelligrini M, Teitell MA, Clark A. A Pluripotency and Self-Renewal Program Controls the Expansion of Genetically Unstable Cancer Stem Cells in Pluripotent Stem Cell-Derived Tumors. *Stem cells (Dayton, Ohio)*. Oct 2 2008.
 91. Cooper LJ, Al-Kadhimi Z, Serrano LM, Pfeiffer T, Olivares S, Castro A, Chang WC, Gonzalez S, Smith D, Forman SJ, Jensen MC. Enhanced antilymphoma efficacy of CD19-redirected influenza MP1-specific CTLs by cotransfer of T cells modified to present influenza MP1. *Blood*. February 2005;105(4):1622-1631.
 92. Corsello SM, Roti G, Ross KN, Chow KT, Galinsky I, DeAngelo DJ, Stone RM, Kung AL, Golub TR, Stegmaier K. Identification of AML1-ETO modulators by chemical genomics. *Blood*. Jun 11 2009;113(24):6193-6205.
 93. Cowey S, Szafran AA, Kappes J, Zinn KR, Siegal GP, Desmond RA, Kim H, Evans L, Hardy RW. Breast cancer metastasis to bone: evaluation of bioluminescent imaging and microSPECT/CT for detecting bone metastasis in immunodeficient mice. *Clin Exp Metastasis*. May 31 2007.
 94. Craft N, Bruhn KW, Nguyen BD, Prins R, Liau LM, Collisson EA, De A, Kolodney MS, Gambhir SS, Miller JF. Bioluminescent imaging of melanoma in live mice. *The Journal of investigative dermatology*. Jul 2005;125(1):159-165.
 95. Creamer BA, Triplett AA, Wagner KU. Longitudinal analysis of mammaryogenesis using a novel tetracycline-inducible mouse model and in vivo imaging. *Genesis*. Apr 2009;47(4):234-245.
 96. Dai Y, Liu M, Tang W, DeSano J, Burstein E, Davis M, Pienta K, Lawrence T, Xu L. Molecularly targeted radiosensitization of human prostate cancer by modulating inhibitor of apoptosis. *Clin Cancer Res*. Dec 1 2008;14(23):7701-7710.
 97. Dalgard CL, Van Quill KR, O'Brien JM. Evaluation of the in vitro and in vivo antitumor activity of histone deacetylase inhibitors for the therapy of

-
- retinoblastoma. *Clin Cancer Res.* May 15 2008;14(10):3113-3123.
98. Dallas NA, Gray MJ, Xia L, Fan F, van Buren G, 2nd, Gaur P, Samuel S, Lim SJ, Arumugam T, Ramachandran V, Wang H, Ellis LM. Neuropilin-2-mediated tumor growth and angiogenesis in pancreatic adenocarcinoma. *Clin Cancer Res.* Dec 15 2008;14(24):8052-8060.
99. Dangi-Garimella S, Yun J, Eves EM, Newman M, Erkeland SJ, Hammond SM, Minn AJ, Rosner MR. Raf kinase inhibitory protein suppresses a metastasis signalling cascade involving LIN28 and let-7. *The EMBO journal.* Feb 18 2009;28(4):347-358.
100. Danielsson A, Dzojic H, Nilsson B, Essand M. Increased therapeutic efficacy of the prostate-specific oncolytic adenovirus Ad[Δ l/PPT-E1A] by reduction of the insulator size and introduction of the full-length E3 region. *Cancer gene therapy.* Apr 2008;15(4):203-213.
101. Darragh MR, Schneider EL, Lou J, Phojanakong PJ, Farady CJ, Marks JD, Hann BC, Craik CS. Tumor detection by imaging proteolytic activity. *Cancer Res.* Feb 15 2010;70(4):1505-1512.
102. De A, Lewis XZ, Gambhir SS. Noninvasive imaging of lentiviral-mediated reporter gene expression in living mice. *Mol Ther.* May 2003;7(5 Pt 1):681-691.
103. Deane NG, Manning HC, Foutch AC, Washington MK, Aronow BJ, Bornhop DJ, Coffey RJ. Targeted imaging of colonic tumors in smad3^{-/-} mice discriminates cancer and inflammation. *Mol Cancer Res.* Apr 2007;5(4):341-349.
104. Deckers M, van Dinther M, Buijs J, Que I, Lowik C, van der Pluijm G, ten Dijke P. The Tumor Suppressor Smad4 Is Required for Transforming Growth Factor {beta}-Induced Epithelial to Mesenchymal Transition and Bone Metastasis of Breast Cancer Cells. *Cancer Res.* February 15, 2006 2006;66(4):2202-2209.
105. Dennis A, Kudo T, Kruidenier L, Girard F, Crepin VF, MacDonald TT, Frankel G, Wiles S. The p50 subunit of NF-kappaB is critical for in vivo clearance of the noninvasive enteric pathogen *Citrobacter rodentium*. *Infection and immunity.* Nov 2008;76(11):4978-4988.
106. Diagaradjane P, Orenstein-Cardona JM, N EC-C, Deorukhkar A, Shentu S, Kuno N, Schwartz DL, Gelovani JG, Krishnan S. Imaging epidermal growth factor receptor expression in vivo: pharmacokinetic and biodistribution characterization of a bioconjugated quantum dot nanoprobe. *Clin Cancer Res.* Feb 1 2008;14(3):731-741.
107. Dickson PV, Hamner JB, Burger RA, Garcia E, Ouma AA, Kim SU, Ng CY, Gray JT, Aboody KS, Danks MK, Davidoff AM. Intravascular administration of tumor tropic neural progenitor cells permits targeted delivery of interferon-beta and restricts tumor growth in a murine model of disseminated neuroblastoma. *J Pediatr Surg.* Jan 2007;42(1):48-53.
108. Dickson PV, Hamner JB, Cauthen LA, Ng CY, McCarville MB, Davidoff AM. Efficacy of zoledronate against neuroblastoma. *Surgery.* Aug 2006;140(2):227-235.
109. Dickson PV, Hamner JB, Streck CJ, Ng CY, McCarville MB, Calabrese C, Gilbertson RJ, Stewart CF, Wilson CM, Gaber MW, Pfeffer LM, Skapek SX,

-
- Nathwani AC, Davidoff AM. Continuous delivery of IFN-beta promotes sustained maturation of intratumoral vasculature. *Mol Cancer Res.* Jun 2007;5(6):531-542.
110. Difeo A, Huang F, Sangodkar J, Terzo EA, Leake D, Narla G, Martignetti JA. KLF6-SV1 is a novel antiapoptotic protein that targets the BH3-only protein NOXA for degradation and whose inhibition extends survival in an ovarian cancer model. *Cancer Res.* Jun 1 2009;69(11):4733-4741.
111. Dinca EB, Lu KV, Sarkaria JN, Pieper RO, Prados MD, Haas-Kogan DA, Vandenberg SR, Berger MS, James CD. p53 Small-molecule inhibitor enhances temozolomide cytotoxic activity against intracranial glioblastoma xenografts. *Cancer Res.* Dec 15 2008;68(24):10034-10039.
112. Dobrenkov K, Olszewska M, Likar Y, Shenker L, Gunset G, Cai S, Pillarsetty N, Hricak H, Sadelain M, Ponomarev V. Monitoring the efficacy of adoptively transferred prostate cancer-targeted human T lymphocytes with PET and bioluminescence imaging. *J Nucl Med.* Jul 2008;49(7):1162-1170.
113. Dong M, How T, Kirkbride KC, Gordon KJ, Lee JD, Hempel N, Kelly P, Moeller BJ, Marks JR, Blobe GC. The type III TGF-beta receptor suppresses breast cancer progression. *The Journal of clinical investigation.* Jan 2007;117(1):206-217.
114. Dong Z, Neiva KG, Jin T, Zhang Z, Hall DE, Mooney DJ, Polverini PJ, Nor JE. Quantification of human angiogenesis in immunodeficient mice using a photon counting-based method. *BioTechniques.* Jul 2007;43(1):73-77.
115. Dorn DC, Kou CA, Png KJ, Moore MA. The effect of cantharidins on leukemic stem cells. *Int J Cancer.* May 1 2009;124(9):2186-2199.
116. Doubrovin M, Serganova I, Mayer-Kuckuk P, Ponomarev V, Blasberg RG. Multimodality in vivo molecular-genetic imaging. *Bioconjugate chemistry.* Nov-December 2004;15(6):1376-1388.
117. Douma S, Van Laar T, Zevenhoven J, Meuwissen R, Van Garderen E, Peeper DS. Suppression of anoikis and induction of metastasis by the neurotrophic receptor TrkB. *Nature.* August 2004;430(7003):1034-1039.
118. Drake JM, Gabriel C, L. , Henry M, D. . Assessing Tumor Growth and Distribution in a Model of Prostate Cancer Metastasis using Bioluminescence Imaging. *Clinical and Experimental Metastasis.* 2005;V22(8):674-684.
119. Drake JM, Strohhahn G, Bair TB, Moreland JG, Henry MD. ZEB1 enhances transendothelial migration and represses the epithelial phenotype of prostate cancer cells. *Molecular biology of the cell.* Apr 2009;20(8):2207-2217.
120. Du J, Zhou Y, Fu QX, Gong WL, Zhao F, Peng JC, Zhan LS. Bioluminescence imaging of hepatitis B virus enhancer and promoter activities in mice. *FEBS letters.* Oct 15 2008;582(23-24):3552-3556.
121. Dubernard G, Oster M, Chareyre F, Antoine M, Rouzier R, Uzan S, Aractingi S, Khosrotehrani K. Increased fetal cell microchimerism in high grade breast carcinomas occurring during pregnancy. *Int J Cancer.* Mar 1 2009;124(5):1054-1059.
122. Dubikovskaya EA, Thorne SH, Pillow TH, Contag CH, Wender PA. Overcoming multidrug resistance of small-molecule therapeutics through conjugation with releasable octaarginine transporters. *Proceedings of the National Academy of*

-
- Sciences of the United States of America*. Aug 26 2008;105(34):12128-12133.
123. Duiker EW, de Vries EG, Mahalingam D, Meersma GJ, Boersma-van Ek W, Hollema H, Lub-de Hooge MN, van Dam GM, Cool RH, Quax WJ, Samali A, van der Zee AG, de Jong S. Enhanced antitumor efficacy of a DR5-specific TRAIL variant over recombinant human TRAIL in a bioluminescent ovarian cancer xenograft model. *Clin Cancer Res*. Mar 15 2009;15(6):2048-2057.
 124. Dumitru CD, Antonysamy MA, Gorski KS, Johnson DD, Reddy LG, Lutterman JL, Piri MM, Proksch J, McGurran SM, Egging EA, Cochran FR, Lipson KE, Tomai MA, Gullikson GW. NK1.1+ cells mediate the antitumor effects of a dual Toll-like receptor 7/8 agonist in the disseminated B16-F10 melanoma model. *Cancer Immunol Immunother*. Apr 2009;58(4):575-587.
 125. Durai M, Krueger C, Ye Z, Cheng L, Mackensen A, Oelke M, Schneck JP. In vivo functional efficacy of tumor-specific T cells expanded using HLA-Ig based artificial antigen presenting cells (aAPC). *Cancer Immunol Immunother*. Feb 2009;58(2):209-220.
 126. Dydensborg AB, Rose AA, Wilson BJ, Grote D, Paquet M, Giguere V, Siegel PM, Bouchard M. GATA3 inhibits breast cancer growth and pulmonary breast cancer metastasis. *Oncogene*. Jul 23 2009;28(29):2634-2642.
 127. Ebos JM, Lee CR, Bogdanovic E, Alami J, Van Slyke P, Francia G, Xu P, Mutsaers AJ, Dumont DJ, Kerbel RS. Vascular endothelial growth factor-mediated decrease in plasma soluble vascular endothelial growth factor receptor-2 levels as a surrogate biomarker for tumor growth. *Cancer Res*. Jan 15 2008;68(2):521-529.
 128. Echchannaoui H, Frei K, Schnell C, Leib SL, Zimmerli W, Landmann R. Toll-like receptor 2-deficient mice are highly susceptible to *Streptococcus pneumoniae* meningitis because of reduced bacterial clearing and enhanced inflammation. *Journal of Infectious Disease*. September 2002;186(6):798-806.
 129. Edgington LE, Berger AB, Blum G, Albrow VE, Paulick MG, Lineberry N, Bogoyo M. Noninvasive optical imaging of apoptosis by caspase-targeted activity-based probes. *Nature medicine*. Aug 2009;15(8):967-973.
 130. Edinger M, Cao YA, Hornig YS, Jenkins DE, Verneris MR, Bachmann MH, Negrin RS, Contag CH. Advancing animal models of neoplasia through in vivo bioluminescence imaging. *Eur J Cancer*. November 2002;38(16):2128-2136.
 131. Edinger M, Cao YA, Verneris MR, Bachmann MH, Contag CH, Negrin RS. Revealing lymphoma growth and the efficacy of immune cell therapies using in vivo bioluminescence imaging. *Blood*. January 2003;101(2):640-648.
 132. Edinger M, Sweeney TJ, Tucker AA, Olomu AB, Negrin RS, Contag CH. Noninvasive assessment of tumor cell proliferation in animal models. *Neoplasia (New York, N.Y.)*. October 1999;1(4):303-310.
 133. Eefting D, Grimbergen JM, De Vries MR, Van Weel V, Kaijzel EL, Que I, Moon RT, Lowik CW, Van Bockel JH, Quax PH. Prolonged In Vivo Gene Silencing by Electroporation-Mediated Plasmid Delivery of Small Interfering RNA. *Human gene therapy*. Sep 12 2007.
 134. Ehrbar M, Zeisberger SM, Raeber GP, Hubbell JA, Schnell C, Zisch AH. The role of actively released fibrin-conjugated VEGF for VEGF receptor 2 gene activation

-
- and the enhancement of angiogenesis. *Biomaterials*. Dec 21 2007.
135. Elankumaran S, Chavan V, Qiao D, Shobana R, Moorkanat G, Biswas M, Samal SK. Type I Interferon Sensitive Recombinant Newcastle Disease Virus For Oncolytic Virotherapy. *Journal of virology*. Feb 10 2010.
 136. Elia L, Aurisicchio L, Facciabene A, Giannetti P, Ciliberto G, La Monica N, Palombo F. CD4+CD25+ regulatory T-cell-inactivation in combination with adenovirus vaccines enhances T-cell responses and protects mice from tumor challenge. *Cancer gene therapy*. Feb 2007;14(2):201-210.
 137. Ellwood-Yen K, Wongvipat J, Sawyers C. Transgenic Mouse Model for Rapid Pharmacodynamic Evaluation of Antiandrogens. *Cancer Res*. November 1, 2006 2006;66(21):10513-10516.
 138. Evgin L, Vaha-Koskela M, Rintoul J, Falls T, Le Boeuf F, Barrett JW, Bell JC, Stanford MM. Potent Oncolytic Activity of Raccoonpox Virus in the Absence of Natural Pathogenicity. *Mol Ther*. Feb 16 2010.
 139. Ezhilarasan R, Jadhav U, Mohanam I, Rao JS, Gujrati M, Mohanam S. The hemopexin domain of MMP-9 inhibits angiogenesis and retards the growth of intracranial glioblastoma xenograft in nude mice. *Int J Cancer*. Jan 15 2009;124(2):306-315.
 140. Fei P, Wang W, Kim SH, Wang S, Burns TF, Sax JK, Buzzai M, Dicker DT, McKenna WG, Bernhard EJ, El-Deiry WS. Bnip3L is induced by p53 under hypoxia, and its knockdown promotes tumor growth. *Cancer cell*. December 2004;6(6):597-609.
 141. Fels DR, Ye J, Segan AT, Kridel SJ, Spiotto M, Olson M, Koong AC, Koumenis C. Preferential cytotoxicity of bortezomib toward hypoxic tumor cells via overactivation of endoplasmic reticulum stress pathways. *Cancer Res*. Nov 15 2008;68(22):9323-9330.
 142. Fleshker S, Preise D, Kalchenko V, Scherz A, Salomon Y. Prompt Assessment of WST11-VTP Outcome Using Luciferase Transfected Tumors Enables Second Treatment and Increase in Overall Therapeutic Rate. *Photochemistry and photobiology*. Apr 8 2008.
 143. Fogal V, Richardson AD, Karmali PP, Scheffler IE, Smith JW, Ruoslahti E. Mitochondrial p32 Protein Is a Critical Regulator of Tumor Metabolism via Maintenance of Oxidative Phosphorylation. *Molecular and cellular biology*. Mar 2010;30(6):1303-1318.
 144. Fogar P, Navaglia F, Basso D, Zambon CF, Moserle L, Indraccolo S, Stranges A, Greco E, Fadi E, Padoan A, Pantano G, Sanzari MC, Pedrazzoli S, Montecucco C, Plebani M. Heat-induced transcription of diphtheria toxin A or its variants, CRM176 and CRM197: implications for pancreatic cancer gene therapy. *Cancer gene therapy*. Jul 17 2009.
 145. Foster AE, Dotti G, Lu A, Khalil M, Brenner MK, Heslop HE, Rooney CM, Bollard CM. Antitumor activity of EBV-specific T lymphocytes transduced with a dominant negative TGF-beta receptor. *J Immunother*. Jun 2008;31(5):500-505.
 146. Francia G, Man S, Lee CJ, Lee CR, Xu P, Mossoba ME, Emmenegger U, Medin JA, Kerbel RS. Comparative impact of trastuzumab and cyclophosphamide on

-
- HER-2-positive human breast cancer xenografts. *Clin Cancer Res.* Oct 15 2009;15(20):6358-6366.
147. Fukukawa C, Hanaoka H, Nagayama S, Tsunoda T, Toguchida J, Endo K, Nakamura Y, Katagiri T. Radioimmunotherapy of human synovial sarcoma using a monoclonal antibody against FZD10. *Cancer science.* Feb 2008;99(2):432-440.
148. Fushiki H, Kanoh-Azuma T, Katoh M, Kawabata K, Jiang J, Tsuchiya N, Satow A, Tamai Y, Hayakawa Y. Quantification of mouse pulmonary cancer models by microcomputed tomography imaging. *Cancer science.* Aug 2009;100(8):1544-1549.
149. Gade TPF, Hassen W, Santos E, Gunset G, Saudemont A, Gong MC, Brentjens R, Zhong X-S, Stephan M, Stefanski J, Lyddane C, Osborne JR, Buchanan IM, Hall SJ, Heston WD, Riviere I, Larson SM, Koutcher JA, Sadelain M. Targeted Elimination of Prostate Cancer by Genetically Directed Human T Lymphocytes. *Cancer Res.* October 1, 2005 2005;65(19):9080-9088.
150. Ganesh S, Gonzalez-Edick M, Gibbons D, Van Roey M, Jooss K. Intratumoral coadministration of hyaluronidase enzyme and oncolytic adenoviruses enhances virus potency in metastatic tumor models. *Clin Cancer Res.* Jun 15 2008;14(12):3933-3941.
151. Gao Y, Ng SS, Chau DH, Yao H, Yang C, Man K, Huang PT, Huang C, Huang JJ, Kung HF, Lin MC. Development of recombinant adeno-associated virus and adenovirus cocktail system for efficient hTERTC27 polypeptide-mediated cancer gene therapy. *Cancer gene therapy.* Nov 2008;15(11):723-732.
152. Gao YJ, Cheng JK, Zeng Q, Xu ZZ, Decosterd I, Xu X, Ji RR. Selective inhibition of JNK with a peptide inhibitor attenuates pain hypersensitivity and tumor growth in a mouse skin cancer pain model. *Experimental neurology.* Sep 2009;219(1):146-155.
153. Garbuzenko OB, Saad M, Betigeri S, Zhang M, Vetcher AA, Soldatenkov VA, Reimer DC, Pozharov VP, Minko T. Intratracheal versus intravenous liposomal delivery of siRNA, antisense oligonucleotides and anticancer drug. *Pharmaceutical research.* Feb 2009;26(2):382-394.
154. Garzia L, Andolfo I, Cusanelli E, Marino N, Petrosino G, De Martino D, Esposito V, Galeone A, Navas L, Esposito S, Gargiulo S, Fattet S, Donofrio V, Cinalli G, Brunetti A, Vecchio LD, Northcott PA, Delattre O, Taylor MD, Iolascon A, Zollo M. MicroRNA-199b-5p impairs cancer stem cells through negative regulation of HES1 in medulloblastoma. *PLoS ONE.* 2009;4(3):e4998.
155. Ghods AJ, Irvin D, Liu G, Yuan X, Abdulkadir IR, Tunici P, Konda B, Wachsmann-Hogiu S, Black KL, Yu JS. Spheres isolated from 9L gliosarcoma rat cell line possess chemoresistant and aggressive cancer stem-like cells. *Stem cells (Dayton, Ohio).* Jul 2007;25(7):1645-1653.
156. Ghoroghchian PP, Frail PR, Susumu K, Blessington D, Brannan AK, Bates FS, Chance B, Hammer DA, Therien MJ. Near-infrared-emissive polymersomes: Self-assembled soft matter for in vivo optical imaging. *PNAS.* February 22, 2005 2005;102(8):2922-2927.
157. Gibson DL, Ma C, Bergstrom KS, Huang JT, Man C, Vallance BA. MyD88

-
- signalling plays a critical role in host defence by controlling pathogen burden and promoting epithelial cell homeostasis during *Citrobacter rodentium*-induced colitis. *Cellular microbiology*. Mar 2008;10(3):618-631.
158. Giri RS, Thaker HM, Giordano T, Williams J, Rogers D, Vasu KK, Sudarsanam V. Design, synthesis and evaluation of novel 2-thiophen-5-yl-3H-quinazolin-4-one analogues as inhibitors of transcription factors NF-small ka, CyrillicB and AP-1 mediated transcriptional activation: Their possible utilization as anti-inflammatory and anti-cancer agents. *Bioorganic & medicinal chemistry*. Jan 11 2010.
159. Goldhoff P, Warrington NM, Limbrick DD, Jr., Hope A, Woerner BM, Jackson E, Perry A, Piwnica-Worms D, Rubin JB. Targeted inhibition of cyclic AMP phosphodiesterase-4 promotes brain tumor regression. *Clin Cancer Res*. Dec 1 2008;14(23):7717-7725.
160. Gondi CS, Veeravalli KK, Gorantla B, Dinh DH, Fassett D, Klopfenstein JD, Gujrati M, Rao JS. Human umbilical cord blood stem cells show PDGF-D-dependent glioma cell tropism in vitro and in vivo. *Neuro-oncology*. May 2010;12(5):453-465.
161. Gopal AK, Press OW, Wilbur SM, Maloney DG, Pagel JM. Rituximab blocks binding of radiolabeled anti-CD20 antibodies (Ab) but not radiolabeled anti-CD45 Ab. *Blood*. Aug 1 2008;112(3):830-835.
162. Goun EA, Shinde R, Dehnert KW, Adams-Bond A, Wender PA, Contag CH, Franc BL. Intracellular cargo delivery by an octaarginine transporter adapted to target prostate cancer cells through cell surface protease activation. *Bioconjugate chemistry*. May-Jun 2006;17(3):787-796.
163. Granot D, Addadi Y, Kalchenko V, Harmelin A, Kunz-Schughart LA, Neeman M. In vivo imaging of the systemic recruitment of fibroblasts to the angiogenic rim of ovarian carcinoma tumors. *Cancer Res*. Oct 1 2007;67(19):9180-9189.
164. Gray MJ, Van Buren G, Dallas NA, Xia L, Wang X, Yang AD, Somcio RJ, Lin YG, Lim S, Fan F, Mangala LS, Arumugam T, Logsdon CD, Lopez-Berestein G, Sood AK, Ellis LM. Therapeutic targeting of neuropilin-2 on colorectal carcinoma cells implanted in the murine liver. *Journal of the National Cancer Institute*. Jan 16 2008;100(2):109-120.
165. Greig JA, Buckley SM, Waddington SN, Parker AL, Bhella D, Pink R, Rahim AA, Morita T, Nicklin SA, McVey JH, Baker AH. Influence of Coagulation Factor X on In Vitro and In Vivo Gene Delivery by Adenovirus (Ad) 5, Ad35, and Chimeric Ad5/Ad35 Vectors. *Mol Ther*. Jul 14 2009.
166. Griesenbach U, Meng C, Farley R, Gardner A, Brake MA, Frankel GM, Gruenert DC, Cheng SH, Scheule RK, Alton EW. The role of doxorubicin in non-viral gene transfer in the lung. *Biomaterials*. Apr 2009;30(10):1971-1977.
167. Gross S, Abraham U, Prior JL, Herzog ED, Piwnica-Worms D. Continuous Delivery of D-Luciferin by Implanted Micro-osmotic Pumps Enables True Real-Time Bioluminescence Imaging of Luciferase Activity In Vivo. *Mol Imaging*. Apr-Jun 2007;6(2):121-130.
168. Gross S, Gammon ST, Moss BL, Rauch D, Harding J, Heinecke JW, Ratner L, Piwnica-Worms D. Bioluminescence imaging of myeloperoxidase activity in vivo. *Nature medicine*. Apr 2009;15(4):455-461.

169. Gross S, Piwnica-Worms D. Spying on cancer: molecular imaging in vivo with genetically encoded reporters. *Cancer cell*. January 2005;7(1):5-15.
170. Gu L, Tsark WM, Brown DA, Blanchard S, Synold TW, Kane SE. A new model for studying tissue-specific *mdr1a* gene expression in vivo by live imaging. *Proceedings of the National Academy of Sciences of the United States of America*. Mar 31 2009;106(13):5394-5399.
171. Gueven N, Fukao T, Luff J, Paterson C, Kay G, Kondo N, Lavin MF. Regulation of the *Atm* promoter in vivo. *Genes, Chromosomes and Cancer*. 2006;45(1):61-71.
172. Gumireddy K, Sun F, Klein-Szanto AJ, Gibbins JM, Gimotty PA, Saunders AJ, Schultz PG, Huang Q. In vivo selection for metastasis promoting genes in the mouse. *Proceedings of the National Academy of Sciences of the United States of America*. Apr 17 2007;104(16):6696-6701.
173. Guse K, Dias JD, Bauerschmitz GJ, Hakkarainen T, Aavik E, Ranki T, Pisto T, Sarkioja M, Desmond RA, Kanerva A, Hemminki A. Luciferase imaging for evaluation of oncolytic adenovirus replication in vivo. *Gene therapy*. Jun 2007;14(11):902-911.
174. Guse K, Ranki T, Ala-Opas M, Bono P, Sarkioja M, Rajeci M, Kanerva A, Hakkarainen T, Hemminki A. Treatment of metastatic renal cancer with capsid-modified oncolytic adenoviruses. *Molecular cancer therapeutics*. Oct 2007;6(10):2728-2736.
175. Guse K, Sloniecka M, Diaconu I, Ottolino-Perry K, Tang N, Ng C, Le Boeuf F, Bell JC, McCart JA, Ristimaki A, Pesonen S, Cerullo V, Hemminki A. Antiangiogenic arming of an oncolytic vaccinia virus enhances antitumor efficacy in renal cell cancer models. *Journal of virology*. Jan 2010;84(2):856-866.
176. Hadaschik BA, Adomat H, Fazli L, Fradet Y, Andersen RJ, Gleave ME, So AI. Intravesical Chemotherapy of High-Grade Bladder Cancer with HTI-286, A Synthetic Analogue of the Marine Sponge Product Hemiasterlin. *Clin Cancer Res*. Mar 1 2008;14(5):1510-1518.
177. Hadaschik BA, Black PC, Sea JC, Metwalli AR, Fazli L, Dinney CP, Gleave ME, So AI. A validated mouse model for orthotopic bladder cancer using transurethral tumour inoculation and bioluminescence imaging. *BJU international*. Sep 11 2007.
178. Hadaschik BA, Jackson J, Fazli L, Zoubeydi A, Burt HM, Gleave ME, So AI. Intravesically administered antisense oligonucleotides targeting heat-shock protein-27 inhibit the growth of non-muscle-invasive bladder cancer. *BJU international*. Aug 5 2008;102(5):610-616.
179. Hadaschik BA, ter Borg MG, Jackson J, Sowery RD, So AI, Burt HM, Gleave ME. Paclitaxel and cisplatin as intravesical agents against non-muscle-invasive bladder cancer. *BJU international*. Jun 2008;101(11):1347-1355.
180. Hajitou A, Lev DC, Hannay JA, Korchin B, Staquicini FI, Soghomonyan S, Alauddin MM, Benjamin RS, Pollock RE, Gelovani JG, Pasqualini R, Arap W. A preclinical model for predicting drug response in soft-tissue sarcoma with targeted AAVP molecular imaging. *Proceedings of the National Academy of Sciences of the United States of America*. Mar 18 2008;105(11):4471-4476.
181. Halder J, Lin YG, Merritt WM, Spannuth WA, Nick AM, Honda T, Kamat AA, Han

-
- LY, Kim TJ, Lu C, Tari AM, Bornmann W, Fernandez A, Lopez-Berestein G, Sood AK. Therapeutic efficacy of a novel focal adhesion kinase inhibitor TAE226 in ovarian carcinoma. *Cancer Res.* Nov 15 2007;67(22):10976-10983.
182. Haley ES, Au GG, Carlton BR, Barry RD, Shafren DR. Regional administration of oncolytic Echovirus 1 as a novel therapy for the peritoneal dissemination of gastric cancer. *Journal of molecular medicine (Berlin, Germany).* Apr 2009;87(4):385-399.
183. Hambardzumyan D, Amankulor NM, Helmy KY, Becher OJ, Holland EC. Modeling Adult Gliomas Using RCAS/t-va Technology. *Translational oncology.* May 2009;2(2):89-95.
184. Hamner JB, Dickson PV, Sims TL, Zhou J, Spence Y, Ng CY, Davidoff AM. Bortezomib inhibits angiogenesis and reduces tumor burden in a murine model of neuroblastoma. *Surgery.* Aug 2007;142(2):185-191.
185. Hamstra DA, Bhojani MS, Griffin LB, Laxman B, Ross BD, Rehemtulla A. Real-time Evaluation of p53 Oscillatory Behavior In vivo Using Bioluminescent Imaging. *Cancer Res.* August 1, 2006 2006;66(15):7482-7489.
186. Han Z, Fu A, Wang H, Diaz R, Geng L, Onishko H, Hallahan DE. Noninvasive assessment of cancer response to therapy. *Nature medicine.* Mar 2008;14(3):343-349.
187. Hara M, Murakami T, Kobayashi E. In vivo bioimaging using photogenic rats: fate of injected bone marrow-derived mesenchymal stromal cells. *Journal of autoimmunity.* May 2008;30(3):163-171.
188. Harada H, Itasaka S, Kizaka-Kondoh S, Shibuya K, Morinibu A, Shinomiya K, Hiraoka M. The Akt/mTOR pathway assures the synthesis of HIF-1alpha protein in a glucose- and reoxygenation-dependent manner in irradiated tumors. *The Journal of biological chemistry.* Feb 20 2009;284(8):5332-5342.
189. Harada H, Itasaka S, Zhu Y, Zeng L, Xie X, Morinibu A, Shinomiya K, Hiraoka M. Treatment regimen determines whether an HIF-1 inhibitor enhances or inhibits the effect of radiation therapy. *British journal of cancer.* Mar 10 2009;100(5):747-757.
190. Harada H, Kizaka-Kondoh S, Hiraoka M. Optical Imaging of Tumor Hypoxia and Evaluation of Efficacy of a Hypoxia-Targeting Drug in Living Animals. *Molecular Imaging.* 2005;4(3):182-193.
191. Harada H, Kizaka-Kondoh S, Li G, Itasaka S, Shibuya K, Inoue M, Hiraoka M. Significance of HIF-1-active cells in angiogenesis and radioresistance. *Oncogene.* Jun 11 2007.
192. Hariri G, Zhang Y, Fu A, Han Z, Brechbiel M, Tantawy MN, Peterson TE, Mernaugh R, Hallahan D. Radiation-guided P-selectin antibody targeted to lung cancer. *Annals of biomedical engineering.* May 2008;36(5):821-830.
193. Harrell MI, Iritani BM, Ruddell A. Tumor-induced sentinel lymph node lymphangiogenesis and increased lymph flow precede melanoma metastasis. *Am J Pathol.* Feb 2007;170(2):774-786.
194. Hart LS, El-Deiry WS. Invincible, but not invisible: imaging approaches toward in vivo detection of cancer stem cells. *J Clin Oncol.* Jun 10 2008;26(17):2901-2910.
195. Hashizume R, Gupta N, Berger MS, Banerjee A, Prados MD, Ayers-Ringler J,

-
- James CD, VandenBerg SR. Morphologic and molecular characterization of ATRT xenografts adapted for orthotopic therapeutic testing. *Neuro-oncology*. Apr 2010;12(4):366-376.
196. Hata N, Shinojima N, Gumin J, Yong R, Marini F, Andreeff M, Lang FF. Platelet-derived growth factor BB mediates the tropism of human mesenchymal stem cells for malignant gliomas. *Neurosurgery*. Jan 2010;66(1):144-156; discussion 156-147.
197. Hatton BA, Villavicencio EH, Tsuchiya KD, Pritchard JI, Ditzler S, Pullar B, Hansen S, Knoblaugh SE, Lee D, Eberhart CG, Hallahan AR, Olson JM. The Smo/Smo model: hedgehog-induced medulloblastoma with 90% incidence and leptomeningeal spread. *Cancer Res*. Mar 15 2008;68(6):1768-1776.
198. Havens AM, Pedersen EA, Shiozawa Y, Ying C, Jung Y, Sun Y, Neeley C, Wang J, Mehra R, Keller ET, McCauley LK, Loberg RD, Pienta KJ, Taichman RS. An in vivo mouse model for human prostate cancer metastasis. *Neoplasia (New York, N.Y.)*. Apr 2008;10(4):371-380.
199. Hawes JJ, Reilly KM. Bioluminescent approaches for measuring tumor growth in a mouse model of neurofibromatosis. *Toxicologic pathology*. 2010;38(1):123-130.
200. He J, Rusckowski M, Wang Y, Dou S, Liu X, Zhang S, Liu G, Hnatowich DJ. Optical pretargeting of tumor with fluorescent MORF oligomers. *Mol Imaging Biol*. Jan-Feb 2007;9(1):17-23.
201. He Y, Rajantie I, Pajusola K, Jeltsch M, Holopainen T, Yla-Herttuala S, Harding T, Jooss K, Takahashi T, Alitalo K. Vascular Endothelial Cell Growth Factor Receptor 3-Mediated Activation of Lymphatic Endothelium Is Crucial for Tumor Cell Entry and Spread via Lymphatic Vessels. *Cancer Res*. June 1, 2005 2005;65(11):4739-4746.
202. Heishi T, Hosaka T, Suzuki Y, Miyashita H, Oike Y, Takahashi T, Nakamura T, Arioka S, Mitsuda Y, Takakura T, Hojo K, Matsumoto M, Yamauchi C, Ohta H, Sonoda H, Sato Y. Endogenous Angiogenesis Inhibitor Vasohibin1 Exhibits Broad-Spectrum Antilymphangiogenic Activity and Suppresses Lymph Node Metastasis. *Am J Pathol*. Feb 4 2010.
203. Hemrajani C, Marches O, Wiles S, Girard F, Dennis A, Dziva F, Best A, Phillips AD, Berger CN, Mousnier A, Crepin VF, Kruidenier L, Woodward MJ, Stevens MP, La Ragione RM, MacDonald TT, Frankel G. Role of NleH, a type III secreted effector from attaching and effacing pathogens, in colonization of the bovine, ovine, and murine gut. *Infection and immunity*. Nov 2008;76(11):4804-4813.
204. Henke E, Perk J, Vider J, de Candia P, Chin Y, Solit DB, Ponomarev V, Cartegni L, Manova K, Rosen N, Benzra R. Peptide-conjugated antisense oligonucleotides for targeted inhibition of a transcriptional regulator in vivo. *Nature biotechnology*. Jan 2008;26(1):91-100.
205. Henriquez NV, van Overveld PG, Que I, Buijs JT, Bachelier R, Kaijzel EL, Lowik CW, Clezardin P, van der Pluijm G. Advances in optical imaging and novel model systems for cancer metastasis research. *Clin Exp Metastasis*. 2007;24(8):699-705.
206. Henshaw J, Mossop B, Yuan F. Relaxin treatment of solid tumors: effects on

-
- electric field-mediated gene delivery. *Molecular cancer therapeutics*. Aug 2008;7(8):2566-2573.
- 207.** Hesselink JW, van Dam GM. [Bioluminescence: the potential of a non-invasive bio-optical imaging technique and improvement of animal research]. *Tijdschr Diergeneeskd*. Apr 1 2007;132(7):244-250.
- 208.** Hickerson RP, Smith FJ, Reeves RE, Contag CH, Leake D, Leachman SA, Milstone LM, McLean WH, Kaspar RL. Single-nucleotide-specific siRNA targeting in a dominant-negative skin model. *The Journal of investigative dermatology*. Mar 2008;128(3):594-605.
- 209.** Hickson J, Ackler S, Klaubert D, Bouska J, Ellis P, Foster K, Oleksijew A, Rodriguez L, Schlessinger S, Wang B, Frost D. Noninvasive molecular imaging of apoptosis in vivo using a modified firefly luciferase substrate, Z-DEVD-aminoluciferin. *Cell death and differentiation*. Jan 8 2010.
- 210.** Hildebrandt IJ, Iyer M, Wagner E, Gambhir SS. Optical imaging of transferrin targeted PEI/DNA complexes in living subjects. *Gene therapy*. May 2003;10(9):758-764.
- 211.** Hingorani M, White CL, Merron A, Peerlinck I, Gore ME, Slade A, Scott SD, Nutting CM, Pandha HS, Melcher AA, Vile RG, Vassaux G, Harrington KJ. Inhibition of repair of radiation-induced DNA damage enhances gene expression from replication-defective adenoviral vectors. *Cancer Res*. Dec 1 2008;68(23):9771-9778.
- 212.** Hingorani M, White CL, Zaidi S, Merron A, Peerlinck I, Gore ME, Nutting CM, Pandha HS, Melcher AA, Vile RG, Vassaux G, Harrington KJ. Radiation-mediated up-regulation of gene expression from replication-defective adenoviral vectors: implications for sodium iodide symporter gene therapy. *Clin Cancer Res*. Aug 1 2008;14(15):4915-4924.
- 213.** Hirbe AC, Roelofs AJ, Floyd DH, Deng H, Becker SN, Lanigan LG, Apicelli AJ, Xu Z, Prior JL, Eagleton MC, Piwnica-Worms D, Rogers MJ, Weilbaecher K. The bisphosphonate zoledronic acid decreases tumor growth in bone in mice with defective osteoclasts. *Bone*. May 2009;44(5):908-916.
- 214.** Hokaiwado N, Takeshita F, Naiki-Ito A, Asamoto M, Ochiya T, Shirai T. Glutathione S-transferase Pi mediates proliferation of androgen-independent prostate cancer cells. *Carcinogenesis*. Jun 2008;29(6):1134-1138.
- 215.** Holland SJ, Pan A, Franci C, Hu Y, Chang B, Li W, Duan M, Torneros A, Yu J, Heckrodt TJ, Zhang J, Ding P, Apatira A, Chua J, Brandt R, Pine P, Goff D, Singh R, Payan DG, Hitoshi Y. R428, a selective small molecule inhibitor of Axl kinase, blocks tumor spread and prolongs survival in models of metastatic breast cancer. *Cancer Res*. Feb 15 2010;70(4):1544-1554.
- 216.** Hollingshead MG, Bonomi CA, Borgel SD, Carter JP, Shoemaker R, Melillo G, Sausville EA. A potential role for imaging technology in anticancer efficacy evaluations. *Eur J Cancer*. April 2004;40(6):890-898.
- 217.** Holopainen T, Huang H, Chen C, Kim KE, Zhang L, Zhou F, Han W, Li C, Yu J, Wu J, Koh GY, Alitalo K, He Y. Angiopoietin-1 overexpression modulates vascular endothelium to facilitate tumor cell dissemination and metastasis establishment.