


小动物活体光学成像技术在基因和细胞治疗中的应用

Revivity小动物活体光学成像技术已在生命科学基础研究、临床前医学研究及药物研发等领域得到广泛应用。在众多应用领域中，基因和细胞治疗是活体光学成像技术的热点之一。在活体光学成像实验中，常用于基因和细胞治疗的光学标记方法包括：1. 利用萤火虫荧光素酶（Firefly Luciferase）或者荧光蛋白作为报告基因，通过转基因技术标记细胞，病毒，细菌或者基因，进而观测基因和细胞治疗效果。2. 通过外源注射功能性荧光探针，观测基因和细胞治疗效果。

一. 实时监测非侵入性基因递送和治疗

1. 不同载体的基因传递


光学成像技术可以用来监测基因传递过程。进行转基因的载体包括病毒载体比如腺病毒（Adenovirus, Ad），疱疹性口腔炎病毒（Vesicular Stomatitis Virus, VSV）和单纯疱疹病毒（Herpes Simplex Virus, HSV）。其他的载体包括化学转染试剂比如脂质转染（Lipofection）和高分子聚合物（Polymer reagents）。下图所示，研究人员使用荧光虫荧光素酶标记的 pMVC-luc 质粒进行多种基因传递方式进行转染，如超声（US-Ultrasound），电穿孔（EP-Electroporation），脂质转染（Lipo-Lipofectamine 2000），腺病毒（Ad-Adenovirus 5）和腺相关病毒（AAV-Adeno-associated virus 2）转导。利用 Revvity 的 IVIS 成像系统可显示出腺病毒和腺相关病毒转导后所显示生物发光最强，说明这两种方法的效率最高。



(Tangney and Francis, Current Gene Therapy, 2012)

上图：在小鼠皮下 JBS 纤维肉瘤移植后，使用超声，电穿孔，脂质转染，腺病毒和腺相关病毒转导进行质粒基因传递。

研究人员使用电穿孔的方法在小鼠的不同部位如肝脏和肌肉部位进行萤火虫荧光素酶标记的 pMVC-luc 质粒转染，利用 Revvity 的 IVIS 成像系统观测转染后的生物发光强度，发现质粒在肌肉中的转染效率更高（下图）。




pMVC-luc 质粒转染到小鼠的肝脏和肌肉部位后，观测体内的生物发光强度。

(Morrisey et al, J Biomed Biotech, 2012)

2. 传递 DNA 治疗疾病


基因传递细胞因子基因比如 IFN 质粒能够调节天生免疫系统和适应性免疫系统。研究人员在回结肠动脉注射萤火虫荧光素标记结肠癌细胞株 (luc-Conlon26 cell) 使肿瘤转移到肝脏，然后静脉注射 IFN-λ 质粒。Revvity 的 IVIS 系统成像结果显示注射 IFN-λ 质粒的小鼠结肠癌生物发光强度明显下降，而且结肠癌小鼠的存活率明显提高。这说明 IFN-λ 质粒传递能够显著的阻止结肠肿瘤细胞的生长，从而帮助治疗人类恶性肿瘤（下图）。



(Sato et al, J. Immunology, 2006)

3. 使用 RNA 治疗疾病


细胞衰老在癌症的发展中起到了重要作用，microRNAs 可以调节细胞衰老。研究人员在萤火虫荧光素酶标记 MDA-MB-231-luc 乳腺癌细胞建立的肿瘤模型中，注射 microRNA（衰老相关 miR22）。Revvity 的 IVIS 成像系统结果显示 miR-22 在小鼠体内能够通过诱导细胞衰老从而抑制肿瘤的生长和转移（下图）。



在小鼠乳腺癌部位注射 miR-22, 46 天后, 小鼠体内的乳腺癌细胞生物发光强度明显下降 (A)。取出不同器官包括肝脏 (L, Liver), 肾脏 (K, Kidney), 脾脏 (Sp, Spleen), 胃 (St, Stomach) 和小肠 (SI, Small Intestine)。这些器官的生物发光强度也明显下降 (B)。


(Xu et al, J Cell Biol, 2011)

神经细胞瘤是从交感神经系统的前体细胞引发的癌症。在 In vitro 研究中, microRNA (miR-184) 能够抑制神经细胞瘤细胞在细胞培养液中的生长。因此, 研究 miR-184 在活体中的作用是非常有意义的。研究人员在小鼠腹膜后注射 miR-184 高表达的萤火虫荧光素标记人神经细胞瘤细胞 (SKN-AS^{luc} or NB1691^{luc}), Revvity 的 IVIS 成像系统结果显示 miR-184 能够抑制神经细胞瘤细胞的发光强度, 而且提高肿瘤小鼠的存活率 (下图)。研究结果说明 miR-184 能够抑制神经细胞瘤细胞在原位异种移植瘤模型内的增殖, 应用 miR-184 是治疗神经胶质瘤的潜在疗法。



(Tivnan et al, BMC Cancer, 2011)

shRNA (short hairpin RNA) 能够用来干扰 RNA 使基因沉默。研究人员使用 shRNA 敲除阿片受体基因的肺癌小鼠，注射 Revvity公司的 Prosensen680 探针探测肿瘤的组织蛋白酶活性从而了解 shRNA 对于肺癌小鼠的治疗效果。使用 Revvity公司的FMT 小动物活体荧光断层成像系统观测 Prosense680 荧光强度，结果显示使用 shRNA 敲除阿片受体基因能够明显的抑制组织蛋白酶活性和肺癌肿瘤细胞的生成（下图）。




(Mathew et al, Int Anaesthesia Res, 2011)

使用野生型 (Wildtype)和 shRNA 敲除阿片受体基因 (MOR Knockout) 的小鼠皮下移植荧光蛋白 (GFP/RFP) 标记的 Lewis 肺癌细胞 (Lewis lung carcinoma, LLC)。尾静脉注射 Prosense680, 使用 FMT 观测和定量荧光强度 (A&B), 图像分析肿瘤细胞 (C)。

二. 细胞治疗


人类胚胎干细胞能够运用于细胞移植治疗，然而在细胞移植后会发生免疫排斥反应，因此研究低抗原性的胚胎干细胞并且提高移植后细胞活性是干细胞移植应用的热点。研究人员移植萤火虫荧光素酶标记人白细胞抗原敲除的胚胎干细胞到具有免疫活性的急性心肌梗塞小鼠模型，免疫排斥反应减轻而且胚胎干细胞的存活率提高（下图）。结果说明敲除白细胞抗原的胚胎干细胞能够激发 T 细胞忽视和减轻异种免疫排斥，低抗原性的胚胎干细胞能够提高移植后细胞的活性。



(Deuse et al, Circulation, 2011)

(A&B) 注射人胚胎干细胞到小鼠心肌梗塞部位 (C) 使用 Revvity公司的 IVIS 系统成像显示胚胎干细胞在免疫缺陷小鼠 (SCID-beige) 中的生物发光时间明显长于具有免疫活性 (Balb/c) 的小鼠 (D) 敲除人白细胞抗原的胚胎干细胞 (siRNA+IB) 在免疫缺陷小鼠 (SCID-beige) 的存活时间明显提高，减低了免疫排斥反应。


人体胚胎素在间充质干细胞和造血干细胞的细胞治疗中有非常重要的作用，包括提高造血干细胞的移植成活率，炎症调节，骨修复和癌症。人胚胎源贴壁细胞是从人产后胚胎分离出的间充质样干细胞。研究人员骨内移植 H929 骨髓瘤细胞到严重联合免疫缺陷 (SCID) 小鼠，静脉注射萤火虫荧光素酶和增强型绿色荧光蛋白 (EGFP) 标记人胚胎源贴壁细胞，Revvity 的 IVIS 系统显示人胚胎源贴壁细胞能够位移到骨髓瘤骨内从而促进细胞生长 (下图)。



静脉注射萤火虫荧光素酶和增强型绿色荧光蛋白标记人胚胎贴壁细胞到骨髓瘤小鼠模型 (A) 探测人胚胎贴壁细胞在小鼠体内生物发光强度 (B) 使用组织切片样品，免疫组化监测绿色荧光蛋白 (GFP) 发光强度。

(Li et al, Stem Cells, 2011)


肝癌是全球第五大癌症，治疗肝癌的首选方法是肝脏移植或者手术切除。肿瘤晚期病人手术切除治疗后有高复发率，因此使用化学疗法是治疗晚期肿瘤患者的第二选择。然而总体来说治疗效果并不理想，主要原因是肿瘤细胞的高化疗抗性和化疗治疗剂本身的毒性，因此研究肝癌的化疗抗性对提高治疗的效率和患者康复有重要作用。研究发现肿瘤干细胞对化学治疗有强的抵抗力，因此肿瘤干细胞对肿瘤在常规治疗后的复发起关键作用。研究者使用异种移植化疗抗性肝癌小鼠，顺铂 (cisplatin) 化学治疗后，CD24 (一种粘蛋白样细胞表面糖蛋白) 表达上调。因此为了研究 CD24 在肿瘤干细胞高化疗抗性的作用，静脉注射 CD24⁺ 人肝癌细胞到免疫缺陷/非肥胖糖尿病 (NOD/SCID) 小鼠，Revvity 的 IVIS 系统成像显示肿瘤在小鼠体内形成，而且能够自我维持，分化和转移，然而敲除 CD24 的人肝癌细胞抑制肿瘤干细胞形态形成 (下图)。结果说明 CD24 信号通路是肝癌患者的治疗靶位。



(A&B) *In vitro* 肿瘤细胞迁移和侵袭实验，CD24⁺ 人肝癌细胞相对于 CD⁻ 人肝癌细胞显示更高的移动和侵入效率。小鼠活体内 (C) 和肝脏组织 (D) 生物发光成像显示，注射 CD24⁺ 人肝癌细胞的小鼠形成肿瘤的比率高。

(Lee et al, Cell, 2011)


使用细菌能够用于实时监测传染性疾病和基因治疗研究的效果。皮下注射萤火虫荧光素酶标记结肠癌细胞株（HCT 116）形成肿瘤小鼠模型，然后尾静脉注射细菌荧光素酶标记的 *B.breve* UCC2003。Revvity的 IVIS 系统三维成像显示细菌荧光素酶标记 *B.breve* 和生物发光标记 HCT116 共定位（下图）。结果说明使用非病原性细菌可以作为载体用于癌症靶位治疗。



(a) 小鼠体内的多模式成像：标记细菌（*B.breve*）生物发光，结肠癌细胞（HCT116）生物发光和 μ CT 成像。（细菌生物发光：橘色，细胞癌生物发光：绿色）。
(b) 小鼠皮下肿瘤内放大成像。活肿瘤细胞（绿色/蓝色），脉管系统（造影剂-红色）和细菌（橘色/黄色）。

(Cronin et al, Plos, 2012)


人表皮生长因子受体 2（HER2）的高表达，点突变和基因删除是和多种癌症（乳腺癌，卵巢癌，胰腺癌，胃癌和直肠癌）紧密相联系的。单核细胞增多性李斯特菌（Lm, *Listeria monocytogenes*）是一种兼性胞内寄生菌，而且可以用于疫苗载体。给予 Balb/c 小鼠新型单核细胞增多性李斯特菌-人表皮生长因子受体 2 嵌合疫苗免疫注射，然后静脉注射萤火虫荧光素酶标记小鼠乳腺癌细胞（4T1-luc）。Revvity的 IVIS 系统成像显示嵌合疫苗免疫乳腺癌肿瘤小鼠的生物发光强度明显弱于单核细胞增多性李斯特菌免疫的对照组（下图）。结果说明细菌受体嵌合疫苗能够明显推迟肺癌肿瘤负载，转移和提高小鼠存活率。



Balb/c 小鼠免疫单核细胞增多性李斯特菌（control Lm），或者单核细胞增多性李斯特菌-人表皮生长因子受体 2 嵌合疫苗（Lm-hHer2/neu chimera）。然后静脉注射萤火虫荧光素酶标记小鼠乳腺癌细胞（4T1-luc）。


(Seavey et al, Clin. Cancer Res, 2009)

利用胰岛细胞移植治疗一型糖尿病的障碍是患者慢性免疫抑制的问题。然而服用免疫抑制药物不仅能够增加传染性疾病和恶性肿瘤发病率而且产生 β 细胞再生的紊乱。研究如何包膜细胞进行移植，降低和消除移植中的免疫抑制是现在研究的热点。研究者使用 TheraCyte 大包裹技术包裹萤火虫荧光素酶标记 β 细胞进行移植，Revivity 的 IVIS 系统成像显示 TheraCyte 大包裹包膜的 β 细胞在异种移植的小鼠中生物发光强度提高，结果说明 TheraCyte 大包裹能够提高异种移植的 β 细胞成活率（下图）。利用 TheraCyte 大包裹技术可以进行细胞治疗糖尿病。



上图：（A）实验原理的示意图。（B）萤火虫荧光素酶标记的新生 FVB 小鼠胰岛细胞包膜，移植到 wtFVB 小鼠。（C）生物发光定量分析结果。（D）无包膜的 FVB 小鼠胰岛细胞通过肾脏间室（Renal compartment, RC）或者皮下（contra-lateral subcutaneous, RC/SC）注射移植到 ICR 小鼠（Control-ICR）。TheraCyte 大包裹包膜的 FVB 小鼠胰岛细胞异种移植到 ICR 小鼠 (Device-ICR)。（E）生物发光定量分析利用 TheraCyte 大包裹技术进行移植和对照结果。

多形性成胶质细胞瘤是最具有侵越性的原发性脑瘤。目前来说，手术切除结合放射和化学疗法是治疗多形性成胶质细胞瘤的最好方法。然而对于中老年患者，放射和化学治疗只能减缓而不能阻止肿瘤细胞生长。CD133+肿瘤干细胞的化学和放射疗法耐药性能够解释传统疗法低效率的原因。正位异种移植萤火虫荧光素酶标记的多形性成胶质细胞瘤细胞到免疫缺陷小鼠（SCID）前额叶，建立多形性成胶质细胞瘤小鼠模型。肿瘤内注射 HER2 靶向的嵌合性 T 细胞，Revivity 的 IVIS 系统成像显示经过特异性 HER2+ T 细胞治疗小鼠的生物发光明显下降（下图）。结果说明特异性 HER2+ T 细胞能够抑制多形性成胶质细胞瘤的生长，从而进行有效治疗多形性成胶质细胞瘤患者。



上图：首行（没有治疗的多形性成胶质细胞瘤小鼠）中行（T 细胞治疗多形性成胶质细胞瘤小鼠）末行（特异性 HER2+ T 细胞治疗多形性成胶质细胞瘤小鼠）。

GENE THERAPY AND DELIVERY PAPERS IMAGED ON THE IVIS

1. Agrawal VK, Copeland KM, Barbachano Y, Rahim A, Seth R, White CL, Hingorani M, Nutting CM, Kelly M, Harris P, Pandha H, Melcher AA, Vile RG, Porter C, Harrington KJ. Microvascular free tissue transfer for gene delivery: in vivo evaluation of different routes of plasmid and adenoviral delivery. *Gene therapy*. Jan 2009;16(1):78-92.
2. Ahmed N, Ratnayake M, Savoldo B, Perlaky L, Dotti G, Wels WS, Bhattacharjee MB, Gilbertson RJ, Shine HD, Weiss HL, Rooney CM, Heslop HE, Gottschalk S. Regression of experimental medulloblastoma following transfer of HER2-specific T cells. *Cancer Res*. Jun 15 2007;67(12):5957-5964.
3. Ahmed N, Salsman VS, Kew Y, Shaffer D, Powell S, Zhang YJ, Grossman RG, Heslop HE, Gottschalk S. HER2-specific T cells target primary glioblastoma stem cells and induce regression of autologous experimental tumors. *Clin Cancer Res*. Jan 15 2010;16(2):474-485.
4. Ahmed N, Salsman VS, Yvon E, Louis CU, Perlaky L, Wels WS, Dishop MK, Kleinerman EE, Pule M, Rooney CM, Heslop HE, Gottschalk S. Immunotherapy for osteosarcoma: genetic modification of T cells overcomes low levels of tumor antigen expression. *Mol Ther*. Oct 2009;17(10):1779-1787.
5. Akimoto T, Sorg BS, Yan Z. Real-time imaging of peroxisome proliferator-activated receptor-gamma coactivator-1alpha promoter activity in skeletal muscles of living mice. *American Journal of Physiology, Cell Physiology*. September 2004;287(3):C790-796.
6. Alter J, Sennoga CA, Lopes DM, Eckersley RJ, Wells DJ. Microbubble stability is a major determinant of the efficiency of ultrasound and microbubble mediated in vivo gene transfer. *Ultrasound Med Biol*. Jun 2009;35(6):976-984.
7. Aoi A, Watanabe Y, Mori S, Takahashi M, Vassaux G, Kodama T. Herpes Simplex Virus Thymidine Kinase-Mediated Suicide Gene Therapy Using Nano/Microbubbles and Ultrasound. *Ultrasound Med Biol*. Dec 18 2007.
8. Arenas F, Hervias I, Uriz M, Joplin R, Prieto J, Medina JF. Combination of ursodeoxycholic acid and glucocorticoids upregulates the AE2 alternate promoter in human liver cells. *J Clin Invest*. Feb 2008;118(2):695-709.
9. Asokan A, Johnson JS, Li C, Samulski RJ. Bioluminescent virion shells: new tools for quantitation of AAV vector dynamics in cells and live animals. *Gene therapy*. Dec 2008;15(24):1618-1622.
10. Aung W, Hasegawa S, Koshikawa-Yano M, Obata T, Ikehira H, Furukawa T, Aoki I, Saga T. Visualization of in vivo electroporation-mediated transgene expression in experimental tumors by optical and magnetic resonance imaging. *Gene therapy*. Jul 2009;16(7):830-839.
11. Aung W, Hasegawa S, Koshikawa-Yano M, Tsuji AB, Sogawa C, Sudo H, Sugyo A, Koizumi M, Furukawa T, Saga T. Noninvasive assessment of regulable transferred-p53 gene expression and evaluation of therapeutic response with FDG-PET in tumor model. *Gene therapy*. May 6 2010.
12. Balani P, Boulaire J, Zhao Y, Zeng J, Lin J, Wang S. High mobility group box2 promoter-controlled suicide gene expression enables targeted glioblastoma treatment. *Mol Ther*. Jun 2009;17(6):1003-1011.
13. Barth AS, Kizana E, Smith RR, Terrovitis J, Dong P, Leppo MK, Zhang Y, Miake J, Olson EN, Schneider JW, Abraham MR, Marban E. Lentiviral vectors bearing the cardiac promoter of the Na⁺-Ca²⁺ exchanger report cardiogenic differentiation in stem cells. *Mol Ther*. May 2008;16(5):957-964.
14. Basile P, Dadali T, Jacobson J, Hasslund S, Ulrich-Vinther M, Soballe K, Nishio Y, Drissi MH, Langstein HN, Mitten DJ, O'Keefe RJ, Schwarz EM, Awad HA. Freeze-dried tendon allografts as tissue-engineering scaffolds for Gdf5 gene delivery. *Mol Ther*. Mar 2008;16(3):466-473.
15. Bayer M, Kantor B, Cockrell A, Ma H, Zeithaml B, Li X, McCown T, Kafri T. A large U3 deletion causes increased in vivo expression from a nonintegrating lentiviral vector. *Mol Ther*. Dec 2008;16(12):1968-1976.
16. Bell JB, Aronovich EL, Schreifels JM, Beadnell TC, Hackett PB. Duration of

- Expression and Activity of Sleeping Beauty Transposase in mouse liver following hydrodynamic DNA delivery. *Mol Ther.* Jul 13 2010.
17. Belur LR, Mclvor RS, Wilber A. Liver-directed gene therapy using the sleeping beauty transposon system. *Methods in molecular biology (Clifton, N.J.)* 2008;434:267-276.
 18. Belur LR, Podetz-Pedersen K, Frandsen J, Mclvor RS. Lung-directed gene therapy in mice using the nonviral Sleeping Beauty transposon system. *Nature protocols.* 2007;2(12):3146-3152.
 19. Bertoni C, Jarrachian S, Wheeler TM, Li Y, Olivares EC, Calos MP, Rando TA. Enhancement of plasmid-mediated gene therapy for muscular dystrophy by directed plasmid integration. *PNAS.* January 10, 2006 2006;103(2):419-424.
 20. Bhaumik S, Gambhir SS. Optical imaging of Renilla luciferase reporter gene expression in living mice. *Proceedings of the National Academy of Sciences of the United States of America.* January 2002;99(1):377-382.
 21. Bisanz K, Yu J, Edlund M, Spohn B, Hung MC, Chung LW, Hsieh CL. Targeting ECM-integrin interaction with liposome-encapsulated small interfering RNAs inhibits the growth of human prostate cancer in a bone xenograft imaging model. *Mol Ther.* Oct 2005;12(4):634-643.
 22. Bortolanza S, Bunuales M, Otano I, Gonzalez-Aseguinolaza G, Ortiz-de-Solorzano C, Perez D, Prieto J, Hernandez-Alcoceba R. Treatment of pancreatic cancer with an oncolytic adenovirus expressing interleukin-12 in Syrian hamsters. *Mol Ther.* Apr 2009;17(4):614-622.
 23. Bradbury MS, Tomishima M, Panagiotakos G, Chan BK, Zanzonico P, Vider J, Ponomarev V, Studer L, Tabar V. Optical bioluminescence imaging of human ES cell progeny in the rodent CNS. *J Neurochem.* Jun 7 2007.
 24. Breitbach CJ, Paterson JM, Lemay CG, Falls TJ, McGuire A, Parato KA, Stojdl DF, Daneshmand M, Speth K, Kirn D, McCart JA, Atkins H, Bell JC. Targeted inflammation during oncolytic virus therapy severely compromises tumor blood flow. *Mol Ther.* Sep 2007;15(9):1686-1693.
 25. Bremer E, van Dam GM, de Bruyn M, van Riezen M, Dijkstra M, Kamps G, Helfrich W, Haisma H. Potent systemic anticancer activity of adenovirally expressed EGFR-selective TRAIL fusion protein. *Mol Ther.* Dec 2008;16(12):1919-1926.
 26. Brian T. Feeley AHCOSLKISYCJRL. In vivo molecular imaging of adenoviral versus lentiviral gene therapy in two bone formation models. *Journal of Orthopaedic Research.* 2006;24(8):1709-1721.
 27. Buckley SM, Howe SJ, Sheard V, Ward NJ, Coutelle C, Thrasher AJ, Waddington SN, McKay TR. Lentiviral transduction of the murine lung provides efficient pseudotype and developmental stage-dependent cell-specific transgene expression. *Gene therapy.* Aug 2008;15(16):1167-1175.
 28. Burton JB, Johnson M, Sato M, Koh SB, Mulholland DJ, Stout D, Chatzioannou AF, Phelps ME, Wu H, Wu L. Adenovirus-mediated gene expression imaging to directly detect sentinel lymph node metastasis of prostate cancer. *Nature medicine.* Aug 2008;14(8):882-888.
 29. Cao F, Drukker M, Lin S, Sheikh AY, Xie X, Li Z, Connolly AJ, Weissman IL, Wu JC. Molecular imaging of embryonic stem cell misbehavior and suicide gene ablation. *Cloning and stem cells.* Spring 2007;9(1):107-117.
 30. Casey G, Cashman JP, Morrissey D, Whelan MC, Larkin JO, Soden DM, Tangney M, O'Sullivan GC. Sonoporation mediated immunogene therapy of solid tumors. *Ultrasound Med Biol.* Mar 2010;36(3):430-440.
 31. Cashman JP, Larkin JO, Casey G, Whelan MC, Collins C, Aarons S, Tangney M, O'Sullivan GC. Immune gene therapy as a neoadjuvant to surgical excision to control metastatic cancers. *Cancer letters.* Jan 5 2008.
 32. Chalberg TW, Genise HL, Vollrath D, Calos MP. phiC31 integrase confers genomic integration and long-term transgene expression in rat retina. *Investigative ophthalmology & visual science.* Jun 2005;46(6):2140-2146.
 33. Chalberg TW, Vankov A, Molnar FE, Butterwick AF, Huie P, Calos MP, Palanker DV. Gene transfer to rabbit retina with electron avalanche transfection. *Investigative ophthalmology & visual science.* Sep 2006;47(9):4083-4090.
 34. Chan J, Kumar S, Fisk NM. First trimester embryo-fetoscopic and ultrasound-guided

- fetal blood sampling for ex vivo viral transduction of cultured human fetal mesenchymal stem cells. *Human reproduction (Oxford, England)*. Nov 2008;23(11):2427-2437.
35. Chanda D, Isayeva T, Kumar S, Siegal GP, Szafran AA, Zinn KR, Reddy VV, Ponnazhagan S. Systemic osteoprotegerin gene therapy restores tumor-induced bone loss in a therapeutic model of breast cancer bone metastasis. *Mol Ther*. May 2008;16(5):871-878.
 36. Chang CL, Wu TC, Hung CF. Control of human mesothelin-expressing tumors by DNA vaccines. *Gene therapy*. Aug 2007;14(16):1189-1198.
 37. Chen CP, Kim JS, Liu D, Rettig GR, McAnuff MA, Martin ME, Rice KG. Synthetic PEGylated Glycoproteins and Their Utility in Gene Delivery. *Bioconjugate chemistry*. Mar-Apr 2007;18(2):371-378.
 38. Cheng CM, Lu YL, Chuang KH, Hung WC, Shiea J, Su YC, Kao CH, Chen BM, Roffler S, Cheng TL. Tumor-targeting prodrug-activating bacteria for cancer therapy. *Cancer Gene Ther*. Jun 2008;15(6):393-401.
 39. Chou CK, Hung JY, Liu JC, Chen CT, Hung MC. An attenuated Salmonella oral DNA vaccine prevents the growth of hepatocellular carcinoma and colon cancer that express alpha-fetoprotein. *Cancer Gene Ther*. Aug 2006;13(8):746-752.
 40. Choy G, Choyke P, Libutti SK. Current advances in molecular imaging: noninvasive in vivo bioluminescent and fluorescent optical imaging in cancer research. *Molecular Imaging*. October 2003;2(4):303-312.
 41. Chuang CM, Monie A, Wu A, Mao CP, Hung CF. Treatment with LL-37 peptide enhances antitumor effects induced by CpG oligodeoxynucleotides against ovarian cancer. *Hum Gene Ther*. Apr 2009;20(4):303-313.
 42. Contag CH, Bachmann MH. Advances in in vivo bioluminescence imaging of gene expression. *Annual Review of Biomedical Engineering*. 2002;4:235-260.
 43. Contag CH, Ross BD. It's not just about anatomy: in vivo bioluminescence imaging as an eyepiece into biology. *Journal of Magnetic Resonance Imaging*. October 2002;16(4):378-387.
 44. Cooney DS, Chambers CB, Suchy H, Hornick MC, Reichensperger J, Neumeister MW. Transgene expression in a model of composite tissue allotransplantation. *Plastic and reconstructive surgery*. Mar 2010;125(3):837-845.
 45. Cronin M, Morrissey D, Rajendran S, El Mashad SM, van Sinderen D, O'Sullivan GC, Tangney M. Orally administered bifidobacteria as vehicles for delivery of agents to systemic tumors. *Mol Ther*. Jul 2010;18(7):1397-1407.
 46. Danielsson A, Dzojic H, Nilsson B, Essand M. Increased therapeutic efficacy of the prostate-specific oncolytic adenovirus Ad[1/PPT-E1A] by reduction of the insulator size and introduction of the full-length E3 region. *Cancer Gene Ther*. Apr 2008;15(4):203-213.
 47. Davila HH, Magee TR, Vernet D, Rajfer J, Gonzalez-Cadavid NF. Gene transfer of inducible nitric oxide synthase complementary DNA regresses the fibrotic plaque in an animal model of Peyronie's disease. *Biology of Reproduction*. November 2004;71(5):1568-1577.
 48. De A, Lewis XZ, Gambhir SS. Noninvasive imaging of lentiviral-mediated reporter gene expression in living mice. *Mol Ther*. May 2003;7(5 Pt 1):681-691.
 49. De A, Yaghoubi SS, Gambhir SS. Applications of lentiviral vectors in noninvasive molecular imaging. *Methods in molecular biology (Clifton, N.J.)*. 2008;433:177-202.
 50. De Laporte L, Yang Y, Zelivyanskaya ML, Cummings BJ, Anderson AJ, Shea LD. Plasmid releasing multiple channel bridges for transgene expression after spinal cord injury. *Mol Ther*. Feb 2009;17(2):318-326.
 51. de Mare A, Lambeck AJ, Regts J, van Dam GM, Nijman HW, Snippe H, Wilschut J, Daemen T. Viral vector-based prime-boost immunization regimens: a possible involvement of T-cell competition. *Gene therapy*. Mar 2008;15(6):393-403.
 52. Deroose CM, Reumers V, Gijssbers R, Bormans G, Debyser Z, Mortelmans L, Baekelandt V. Noninvasive Monitoring of Long-Term Lentiviral Vector-Mediated Gene Expression in Rodent Brain with Bioluminescence Imaging. *Molecular Therapy*. 2006;In Press, Corrected Proof.
 53. Deuse T, Peter C, Fedak PW, Doyle T, Reichenspurner H, Zimmermann WH,

- Eschenhagen T, Stein W, Wu JC, Robbins RC, Schrepfer S. Hepatocyte growth factor or vascular endothelial growth factor gene transfer maximizes mesenchymal stem cell-based myocardial salvage after acute myocardial infarction. *Circulation*. Sep 15 2009;120(11 Suppl):S247-254.
54. Di Stasi A, De Angelis B, Rooney CM, Zhang L, Mahendravada A, Foster AE, Heslop HE, Brenner MK, Dotti G, Savoldo B. T lymphocytes coexpressing CCR4 and a chimeric antigen receptor targeting CD30 have improved homing and antitumor activity in a Hodgkin tumor model. *Blood*. Jun 18 2009;113(25):6392-6402.
 55. Dickson PV, Hamner JB, Burger RA, Garcia E, Ouma AA, Kim SU, Ng CY, Gray JT, Aboody KS, Danks MK, Davidoff AM. Intravascular administration of tumor tropic neural progenitor cells permits targeted delivery of interferon-beta and restricts tumor growth in a murine model of disseminated neuroblastoma. *Journal of pediatric surgery*. Jan 2007;42(1):48-53.
 56. Doubrovin M, Serganova I, Mayer-Kuckuk P, Ponomarev V, Blasberg RG. Multimodality in vivo molecular-genetic imaging. *Bioconjugate chemistry*. Nov-December 2004;15(6):1376-1388.
 57. Duda J, Karimi M, Negrin RS, Contag CH. Methods for imaging cell fates in hematopoiesis. *Methods in molecular medicine*. 2007;134:17-34.
 58. Dunn CA, Jin Q, Taba M, Jr., Franceschi RT, Bruce Rutherford R, Giannobile WV. BMP gene delivery for alveolar bone engineering at dental implant defects. *Molecular Therapy*. February 2005;11(2):294-299.
 59. Ebos JM, Lee CR, Bogdanovic E, Alami J, Van Slyke P, Francia G, Xu P, Mutsaers AJ, Dumont DJ, Kerbel RS. Vascular endothelial growth factor-mediated decrease in plasma soluble vascular endothelial growth factor receptor-2 levels as a surrogate biomarker for tumor growth. *Cancer Res*. Jan 15 2008;68(2):521-529.
 60. Eefting D, Grimbergen JM, De Vries MR, Van Weel V, Kaijzel EL, Que I, Moon RT, Lowik CW, Van Bockel JH, Quax PH. Prolonged In Vivo Gene Silencing by Electroporation-Mediated Plasmid Delivery of Small Interfering RNA. *Hum Gene Ther*. Sep 12 2007.
 61. Elfinger M, Geiger J, Hasenpusch G, Uzgun S, Sieverling N, Aneja MK, Maucksch C, Rudolph C. Targeting of the beta(2)-adrenoceptor increases nonviral gene delivery to pulmonary epithelial cells in vitro and lungs in vivo. *J Control Release*. May 5 2009;135(3):234-241.
 62. Elia L, Aurisicchio L, Facciabene A, Giannetti P, Ciliberto G, La Monica N, Palombo F. CD4+CD25+ regulatory T-cell-inactivation in combination with adenovirus vaccines enhances T-cell responses and protects mice from tumor challenge. *Cancer Gene Ther*. Feb 2007;14(2):201-210.
 63. Evgin L, Vaha-Koskela M, Rintoul J, Falls T, Le Boeuf F, Barrett JW, Bell JC, Stanford MM. Potent Oncolytic Activity of Raccoonpox Virus in the Absence of Natural Pathogenicity. *Mol Ther*. Feb 16 2010.
 64. Excoffon KJ, Koerber JT, Dickey DD, Murtha M, Keshavjee S, Kaspar BK, Zabner J, Schaffer DV. Directed evolution of adeno-associated virus to an infectious respiratory virus. *Proc Natl Acad Sci U S A*. Mar 10 2009;106(10):3865-3870.
 65. Fan S, Maguire CA, Ramirez SH, Bradel-Tretheway B, Sapinoro R, Sui Z, Chakraborty-Sett S, Dewhurst S. Valproic acid enhances gene expression from viral gene transfer vectors. *J Virol Methods*. April 2005;125(1):23-33.
 66. Fein DE, Limberis MP, Maloney SF, Heath JM, Wilson JM, Diamond SL. Cationic lipid formulations alter the in vivo tropism of AAV2/9 vector in lung. *Mol Ther*. Dec 2009;17(12):2078-2087.
 67. Fishbein I, Alferiev I, Bakay M, Stachelek SJ, Sobolewski P, Lai M, Choi H, Chen IW, Levy RJ. Local delivery of gene vectors from bare-metal stents by use of a biodegradable synthetic complex inhibits in-stent restenosis in rat carotid arteries. *Circulation*. Apr 22 2008;117(16):2096-2103.
 68. Fogar P, Navaglia F, Basso D, Zambon CF, Moserle L, Indraccolo S, Stranges A, Greco E, Fadi E, Padoan A, Pantano G, Sanzari MC, Pedrazzoli S, Montecucco C, Plebani M. Heat-induced transcription of diphtheria toxin A or its variants, CRM176 and CRM197: implications for pancreatic cancer gene therapy. *Cancer Gene Ther*. Jul 17 2009.

69. Fomicheva EV, Turner, II, Edwards TG, Hoff J, Arden E, D'Alecy LG, Metzger JM. Double Oxygen-sensing Vector System for Robust Hypoxia/Ischemia-regulated Gene Induction in Cardiac Muscle In Vitro and In Vivo. *Mol Ther.* Jun 24 2008.
70. Foster AE, Dotti G, Lu A, Khalil M, Brenner MK, Heslop HE, Rooney CM, Bollard CM. Antitumor activity of EBV-specific T lymphocytes transduced with a dominant negative TGF-beta receptor. *J Immunother.* Jun 2008;31(5):500-505.
71. Franc BL, Mandl SJ, Sivrashvili Z, Wender P, Contag CH. Breaching biological barriers: protein translocation domains as tools for molecular imaging and therapy. *Molecular Imaging.* October 2003;2(4):313-323.
72. Frank SJ, Wang X, He K, Yang N, Fang P, Rosenfeld RG, Hwa V, Chaudhuri TR, Deng L, Zinn KR. In Vivo Imaging of Hepatic Growth Hormone Signaling. *Mol Endocrinol.* November 1, 2006 2006;20(11):2819-2830.
73. Fujishiro J, Kawana H, Inoue S, Shimizu H, Yoshino H, Hakamata Y, Kaneko T, Murakami T, Hashizume K, Kobayashi E. Efficiency of adenovirus-mediated gene transduction in heart grafts in rats. *Transplantation proceedings.* Jan-Feb 2005;37(1):67-69.
74. Gallaher SD, Gil JS, Dorigo O, Berk AJ. Robust in vivo transduction of a genetically stable Epstein-Barr virus episome to hepatocytes in mice by a hybrid viral vector. *Journal of virology.* Apr 2009;83(7):3249-3257.
75. Gao Y, Ng SS, Chau DH, Yao H, Yang C, Man K, Huang PT, Huang C, Huang JJ, Kung HF, Lin MC. Development of recombinant adeno-associated virus and adenovirus cocktail system for efficient hTERTC27 polypeptide-mediated cancer gene therapy. *Cancer Gene Ther.* Nov 2008;15(11):723-732.
76. Garbuzenko OB, Saad M, Betigeri S, Zhang M, Vetcher AA, Soldatenkov VA, Reimer DC, Pozharov VP, Minko T. Intratracheal versus intravenous liposomal delivery of siRNA, antisense oligonucleotides and anticancer drug. *Pharmaceutical research.* Feb 2009;26(2):382-394.
77. Giering JC, Grimm D, Storm TA, Kay MA. Expression of shRNA from a tissue-specific pol II promoter is an effective and safe RNAi therapeutic. *Mol Ther.* Sep 2008;16(9):1630-1636.
78. Gould D, Yousaf N, Fatah R, Subang MC, Chernajovsky Y. Gene therapy with an improved doxycycline-regulated plasmid encoding a tumour necrosis factor-alpha inhibitor in experimental arthritis. *Arthritis research & therapy.* 2007;9(1):R7.
79. Greig JA, Buckley SM, Waddington SN, Parker AL, Bhella D, Pink R, Rahim AA, Morita T, Nicklin SA, McVey JH, Baker AH. Influence of Coagulation Factor X on In Vitro and In Vivo Gene Delivery by Adenovirus (Ad) 5, Ad35, and Chimeric Ad5/Ad35 Vectors. *Mol Ther.* Jul 14 2009.
80. Greish K, Frandsen J, Scharff S, Gustafson J, Cappello J, Li D, O'Malley BW, Jr., Ghandehari H. Silk-elasticlike protein polymers improve the efficacy of adenovirus thymidine kinase enzyme prodrug therapy of head and neck tumors. *The journal of gene medicine.* Jul 2010;12(7):572-579.
81. Griesenbach U, Meng C, Farley R, Cheng SH, Scheule RK, Davies MH, Wolstenholme-Hogg PC, ten Hove W, van der Hoeven P, Sinn PL, McCray PB, Jr., Inoue M, Geddes DM, Hasegawa M, Frankel G, Wiles S, Alton EW. In vivo imaging of gene transfer to the respiratory tract. *Biomaterials.* Apr 2008;29(10):1533-1540.
82. Griesenbach U, Meng C, Farley R, Cheng SH, Scheule RK, Davies MH, Wolstenholme-Hogg PC, Ten Hove W, van der Hoeven P, Sinn PL, McCray PB, Jr., Inoue M, Geddes DM, Hasegawa M, Frankel G, Wiles S, Alton EW. In vivo imaging of gene transfer to the respiratory tract. *Biomaterials.* Dec 19 2007.
83. Griesenbach U, Meng C, Farley R, Gardner A, Brake MA, Frankel GM, Gruenert DC, Cheng SH, Scheule RK, Alton EW. The role of doxorubicin in non-viral gene transfer in the lung. *Biomaterials.* Apr 2009;30(10):1971-1977.
84. Guse K, Dias JD, Bauerschmitz GJ, Hakkarainen T, Aavik E, Ranki T, Pisto T, Sarkioja M, Desmond RA, Kanerva A, Hemminki A. Luciferase imaging for evaluation of oncolytic adenovirus replication in vivo. *Gene therapy.* Jun 2007;14(11):902-911.
85. Guse K, Ranki T, Ala-Opas M, Bono P, Sarkioja M, Rajewski M, Kanerva A, Hakkarainen T, Hemminki A. Treatment of metastatic renal cancer with capsid-modified oncolytic adenoviruses. *Molecular cancer therapeutics.* Oct

- 2007;6(10):2728-2736.
86. Guse K, Sloniecka M, Diaconu I, Ottolino-Perry K, Tang N, Ng C, Le Boeuf F, Bell JC, McCart JA, Ristimaki A, Pesonen S, Cerullo V, Hemminki A. Antiangiogenic arming of an oncolytic vaccinia virus enhances antitumor efficacy in renal cell cancer models. *Journal of virology*. Jan 2010;84(2):856-866.
 87. Hajitou A, Lev DC, Hannay JA, Korchin B, Staquicini FI, Soghomonyan S, Alauddin MM, Benjamin RS, Pollock RE, Gelovani JG, Pasqualini R, Arap W. A preclinical model for predicting drug response in soft-tissue sarcoma with targeted AAVP molecular imaging. *Proc Natl Acad Sci U S A*. Mar 18 2008;105(11):4471-4476.
 88. Hallden G, Thorne SH, Yang J, Kirn DH. Replication-selective oncolytic adenoviruses. *Methods in molecular medicine*. 2004;90:71-90.
 89. Hamed HA, Yacoub A, Park MA, Eulitt PJ, Dash R, Sarkar D, Dmitriev IP, Lesniak MS, Shah K, Grant S, Curiel DT, Fisher PB, Dent P. Inhibition of Multiple Protective Signaling Pathways and Ad.5/3 Delivery Enhances mda-7/IL-24 Therapy of Malignant Glioma. *Mol Ther*. Feb 23 2010.
 90. Hauck ES, Zou S, Scarfo K, Nantz MH, Hecker JG. Whole animal in vivo imaging after transient, nonviral gene delivery to the rat central nervous system. *Mol Ther*. Nov 2008;16(11):1857-1864.
 91. Heidel JD, Liu JY, Yen Y, Zhou B, Heale BS, Rossi JJ, Bartlett DW, Davis ME. Potent siRNA inhibitors of ribonucleotide reductase subunit RRM2 reduce cell proliferation in vitro and in vivo. *Clin Cancer Res*. Apr 1 2007;13(7):2207-2215.
 92. Henshaw J, Mossop B, Yuan F. Relaxin treatment of solid tumors: effects on electric field-mediated gene delivery. *Molecular cancer therapeutics*. Aug 2008;7(8):2566-2573.
 93. Hickerson RP, Smith FJ, Reeves RE, Contag CH, Leake D, Leachman SA, Milstone LM, McLean WH, Kaspar RL. Single-nucleotide-specific siRNA targeting in a dominant-negative skin model. *The Journal of investigative dermatology*. Mar 2008;128(3):594-605.
 94. Hildebrandt IJ, Gambhir SS. Molecular imaging applications for immunology. *Clin Immunol*. May 2004;111(2):210-224.
 95. Hingorani M, White CL, Merron A, Peerlinck I, Gore ME, Slade A, Scott SD, Nutting CM, Pandha HS, Melcher AA, Vile RG, Vassaux G, Harrington KJ. Inhibition of repair of radiation-induced DNA damage enhances gene expression from replication-defective adenoviral vectors. *Cancer Res*. Dec 1 2008;68(23):9771-9778.
 96. Hingorani M, White CL, Zaidi S, Merron A, Peerlinck I, Gore ME, Nutting CM, Pandha HS, Melcher AA, Vile RG, Vassaux G, Harrington KJ. Radiation-mediated up-regulation of gene expression from replication-defective adenoviral vectors: implications for sodium iodide symporter gene therapy. *Clin Cancer Res*. Aug 1 2008;14(15):4915-4924.
 97. Hino A, Fukuyama S, Kataoka K, Kweon MN, Fujihashi K, Kiyono H. Nasal IL-12p70 DNA prevents and treats intestinal allergic diarrhea. *J Immunol*. Jun 1 2005;174(11):7423-7432.
 98. Hiraoka K, Kimura T, Logg CR, Kasahara N. Tumor-selective gene expression in a hepatic metastasis model after locoregional delivery of a replication-competent retrovirus vector. *Clin Cancer Res*. Dec 1 2006;12(23):7108-7116.
 99. Hiraoka K, Kimura T, Logg CR, Tai CK, Haga K, Lawson GW, Kasahara N. Therapeutic efficacy of replication-competent retrovirus vector-mediated suicide gene therapy in a multifocal colorectal cancer metastasis model. *Cancer Res*. Jun 1 2007;67(11):5345-5353.
 100. Holopainen T, Huang H, Chen C, Kim KE, Zhang L, Zhou F, Han W, Li C, Yu J, Wu J, Koh GY, Alitalo K, He Y. Angiopoietin-1 overexpression modulates vascular endothelium to facilitate tumor cell dissemination and metastasis establishment. *Cancer Res*. Jun 1 2009;69(11):4656-4664.
 101. Hoyos V, Savoldo B, Quintarelli C, Mahendravada A, Zhang M, Vera J, Heslop HE, Rooney CM, Brenner MK, Dotti G. Engineering CD19-specific T lymphocytes with interleukin-15 and a suicide gene to enhance their anti-lymphoma/leukemia effects and safety. *Leukemia*. Jun 2010;24(6):1160-1170.
 102. Huang M, Chan DA, Jia F, Xie X, Li Z, Hoyt G, Robbins RC, Chen X, Giaccia AJ, Wu

- JC. Short hairpin RNA interference therapy for ischemic heart disease. *Circulation*. Sep 30 2008;118(14 Suppl):S226-233.
103. Huang P, Kaku H, Chen J, Kashiwakura Y, Saika T, Nasu Y, Urata Y, Fujiwara T, Watanabe M, Kumon H. Potent antitumor effects of combined therapy with a telomerase-specific, replication-competent adenovirus (OBP-301) and IL-2 in a mouse model of renal cell carcinoma. *Cancer Gene Ther*. Feb 19 2010.
 104. Huang X, Guo H, Kang J, Choi S, Zhou TC, Tammana S, Lees CJ, Li ZZ, Milone M, Levine BL, Tolar J, June CH, Scott Mclvor R, Wagner JE, Blazar BR, Zhou X. Sleeping Beauty transposon-mediated engineering of human primary T cells for therapy of CD19+ lymphoid malignancies. *Mol Ther*. Mar 2008;16(3):580-589.
 105. Huang YH, Zugates GT, Peng W, Holtz D, Dunton C, Green JJ, Hossain N, Chernick MR, Padera RF, Jr., Langer R, Anderson DG, Sawicki JA. Nanoparticle-delivered suicide gene therapy effectively reduces ovarian tumor burden in mice. *Cancer Res*. Aug 1 2009;69(15):6184-6191.
 106. Hung CF, Tsai YC, He L, Coukos G, Fodor I, Qin L, Levitsky H, Wu TC. Vaccinia virus preferentially infects and controls human and murine ovarian tumors in mice. *Gene therapy*. Jan 2007;14(1):20-29.
 107. Hung CF, Tsai YC, He L, Wu TC. Control of mesothelin-expressing ovarian cancer using adoptive transfer of mesothelin peptide-specific CD8(+) T cells. *Gene therapy*. Mar 22 2007.
 108. Ibrahim A, Vande Velde G, Reumers V, Toelen J, Thiry I, Vandeputte C, Vets S, Deroose C, Bormans G, Baekelandt V, Debyser Z, Gijsbers R. Highly efficient multicistronic lentiviral vectors with peptide 2A sequences. *Hum Gene Ther*. Aug 2009;20(8):845-860.
 109. Itaka K, Osada K, Morii K, Kim P, Yun SH, Kataoka K. Polyplex nanomicelle promotes hydrodynamic gene introduction to skeletal muscle. *J Control Release*. Jan 3 2010.
 110. Iyer M, Berenji M, Templeton NS, Gambhir SS. Noninvasive imaging of cationic lipid-mediated delivery of optical and PET reporter genes in living mice. *Molecular Therapy*. October 2002;6(4):555-562.
 111. Iyer M, Salazar FB, Lewis X, Zhang L, Carey M, Wu L, Gambhir SS. Noninvasive imaging of enhanced prostate-specific gene expression using a two-step transcriptional amplification-based lentivirus vector. *Molecular Therapy*. September 2004;10(3):545-552.
 112. Iyer M, Salazar FB, Wu L, Carey M, Gambhir SS. Bioluminescence Imaging of Systemic Tumor Targeting Using a Prostate-Specific Lentiviral Vector. *Human Gene Therapy*. 2006;17(1):125-132.
 113. Iyer M, Wu L, Carey M, Wang Y, Smallwood A, Gambhir SS. Two-step transcriptional amplification as a method for imaging reporter gene expression using weak promoters. *Proc Natl Acad Sci U S A*. December 2001;98(25):14595-14600.
 114. Jearawiriyapaisarn N, Moulton HM, Buckley B, Roberts J, Sazani P, Fucharoen S, Iversen PL, Kole R. Sustained dystrophin expression induced by peptide-conjugated morpholino oligomers in the muscles of mdx mice. *Mol Ther*. Sep 2008;16(9):1624-1629.
 115. Jeon YH, Choi Y, Yoon SO, Kim CW, Chung JK. Synergistic tumoricidal effect of combined hMUC1 vaccination and hNIS radioiodine gene therapy. *Molecular cancer therapeutics*. Jul 2008;7(7):2252-2260.
 116. Jiang SN, Phan TX, Nam TK, Nguyen VH, Kim HS, Bom HS, Choy HE, Hong Y, Min JJ. Inhibition of Tumor Growth and Metastasis by a Combination of Escherichia coli-mediated Cytolytic Therapy and Radiotherapy. *Mol Ther*. Mar 2010;18(3):635-642.
 117. Jin Q, Anusaksathien O, Webb SA, Printz MA, Giannobile WV. Engineering of tooth-supporting structures by delivery of PDGF gene therapy vectors. *Mol Ther*. Apr 2004;9(4):519-526.
 118. Johnson M, Huyn S, Burton J, Sato M, Wu L. Differential Biodistribution of Adenoviral Vector In Vivo as Monitored by Bioluminescence Imaging and Quantitative Polymerase Chain Reaction. *Hum Gene Ther*. Nov 15 2006.
 119. Kanerva A, Zinn KR, Chaudhuri TR, Lam JT, Suzuki K, Uil TG, Hakkarainen T, Bauerschmitz GJ, Wang M, Liu B, Cao Z, Alvarez RD, Curiel DT, Hemminki A.

- Enhanced therapeutic efficacy for ovarian cancer with a serotype 3 receptor-targeted oncolytic adenovirus. *Molecular Therapy*. September 2003;8(3):449-458.
120. Kanerva A, Zinn KR, Peng KW, Ranki T, Kangasniemi L, Chaudhuri TR, Desmond RA, Wang M, Takayama K, Hakkarainen T, Alfthan H, Stenman UH, Curiel DT, Hemminki A. Noninvasive dual modality in vivo monitoring of the persistence and potency of a tumor targeted conditionally replicating adenovirus. *Gene therapy*. January 2005;12(1):87-94.
 121. Kesarwala AH, Prior JL, Sun J, Harpstrite SE, Sharma V, Piwnicka-Worms D. Second-generation triple reporter for bioluminescence, micro-positron emission tomography, and fluorescence imaging. *Mol Imaging*. Oct-Dec 2006;5(4):465-474.
 122. Kim D, Hung CF, Wu TC. Monitoring the trafficking of adoptively transferred antigen-specific CD8-positive T cells in vivo, using noninvasive luminescence imaging. *Hum Gene Ther*. Jul 2007;18(7):575-588.
 123. Kim D, Monie A, He L, Tsai YC, Hung CF, Wu TC. Role of IL-2 secreted by PADRE-specific CD4+ T cells in enhancing E7-specific CD8+ T-cell immune responses. *Gene therapy*. May 2008;15(9):677-687.
 124. Kim SI, Shin D, Choi TH, Lee JC, Cheon GJ, Kim KY, Park M, Kim M. Systemic and specific delivery of small interfering RNAs to the liver mediated by apolipoprotein A-I. *Mol Ther*. Jun 2007;15(6):1145-1152.
 125. Klein RL, Dayton RD, Leidenheimer NJ, Jansen K, Golde TE, Zweig RM. Efficient neuronal gene transfer with AAV8 leads to neurotoxic levels of tau or green fluorescent proteins. *Mol Ther*. Mar 2006;13(3):517-527.
 126. Klein RL, Dayton RD, Tatom JB, Henderson KM, Henning PP. AAV8, 9, Rh10, Rh43 vector gene transfer in the rat brain: effects of serotype, promoter and purification method. *Mol Ther*. Jan 2008;16(1):89-96.
 127. Ko IK, Ziady A, Lu S, Kwon YJ. Acid-degradable cationic methacrylamide polymerized in the presence of plasmid DNA as tunable non-viral gene carrier. *Biomaterials*. Oct 2008;29(28):3872-3881.
 128. Krom Y, Fallaux F, Que I, Lowik C, van Dijk K. Efficient in vivo knock-down of estrogen receptor alpha: application of recombinant adenovirus vectors for delivery of short hairpin RNA. *BMC biotechnology*. February 28 2006;6(11):1-10.
 129. Kumar S, Wan C, Ramaswamy G, Clemens TL, Ponnazhagan S. Mesenchymal Stem Cells Expressing Osteogenic and Angiogenic Factors Synergistically Enhance Bone Formation in a Mouse Model of Segmental Bone Defect. *Mol Ther*. Jan 12 2010.
 130. Kutschka I, Kofidis T, Chen IY, von Degenfeld G, Zwierzchoniewska M, Hoyt G, Arai T, Lebl DR, Hendry SL, Sheikh AY, Cooke DT, Connolly A, Blau HM, Gambhir SS, Robbins RC. Adenoviral Human BCL-2 Transgene Expression Attenuates Early Donor Cell Death After Cardiomyoblast Transplantation Into Ischemic Rat Hearts. *Circulation*. July 4, 2006 2006;114(1_suppl):I-174-180.
 131. Kwon SM, Nam HY, Nam T, Park K, Lee S, Kim K, Kwon IC, Kim J, Kang D, Park JH, Jeong SY. In vivo time-dependent gene expression of cationic lipid-based emulsion as a stable and biocompatible non-viral gene carrier. *J Control Release*. May 22 2008;128(1):89-97.
 132. Lalani AS, Chang B, Lin J, Case SS, Luan B, Wu-Prior WW, VanRoey M, Jooss K. Anti-tumor efficacy of human angiostatin using liver-mediated adeno-associated virus gene therapy. *Mol Ther*. Jan 2004;9(1):56-66.
 133. Lee JY, Kim S, Hwang DW, Jeong JM, Chung JK, Lee MC, Lee DS. Development of a Dual-Luciferase Reporter System for In Vivo Visualization of MicroRNA Biogenesis and Posttranscriptional Regulation. *J Nucl Med*. Feb 2008;49(2):285-294.
 134. Leo BM, Li X, Balian G, Anderson DG. In vivo bioluminescent imaging of virus-mediated gene transfer and transduced cell transplantation in the intervertebral disc. *Spine*. April 2004;29(8):838-844.
 135. Leyton J, Lockley M, Aerts JL, Baird SK, Aboagye EO, Lemoine NR, McNeish IA. Quantifying the Activity of Adenoviral E1A CR2 Deletion Mutants Using Renilla Luciferase Bioluminescence and 3'-Deoxy-3'-[18F]Fluorothymidine Positron Emission Tomography Imaging. *Cancer Res*. September 15, 2006 2006;66(18):9178-9185.
 136. Li H-w, Li J, Helm GA, Pan D. Highly specific expression of luciferase gene in lungs of naive nude mice directed by prostate-specific antigen promoter. *Biochemical and*

- Biophysical Research Communications*. 2005;334(4):1287-1291.
137. Li HJ, Everts M, Pereboeva L, Komarova S, Idan A, Curiel DT, Herschman HR. Adenovirus tumor targeting and hepatic untargeting by a coxsackie/adenovirus receptor ectodomain anti-carcinoembryonic antigen bispecific adapter. *Cancer Res*. Jun 1 2007;67(11):5354-5361.
 138. Li HJ, Everts M, Yamamoto M, Curiel DT, Herschman HR. Combined transductional untargeting/retargeting and transcriptional restriction enhances adenovirus gene targeting and therapy for hepatic colorectal cancer tumors. *Cancer Res*. Jan 15 2009;69(2):554-564.
 139. Li JZ, Holman D, Li H, Liu AH, Beres B, Hankins GR, Helm GA. Long-term tracing of adenoviral expression in rat and rabbit using luciferase imaging. *The journal of gene medicine*. Jun 2005;7(6):792-802.
 140. Li SD, Chen YC, Hackett MJ, Huang L. Tumor-targeted delivery of siRNA by self-assembled nanoparticles. *Mol Ther*. Jan 2008;16(1):163-169.
 141. Li YS, Davidson E, Reid CN, McHale AP. Optimising ultrasound-mediated gene transfer (sonoporation) in vitro and prolonged expression of a transgene in vivo: potential applications for gene therapy of cancer. *Cancer letters*. Jan 8 2009;273(1):62-69.
 142. Li Z, Suzuki Y, Huang M, Cao F, Xie X, Connolly AJ, Yang PC, Wu JC. Comparison of Reporter Gene and Iron Particle Labeling for Tracking Fate of Human Embryonic Stem Cells and Differentiated Endothelial Cells in Living Subjects. *Stem cells (Dayton, Ohio)*. Jan 24 2008.
 143. Liang Q, Yamamoto M, Curiel DT, Herschman HR. Noninvasive imaging of transcriptionally restricted transgene expression following intratumoral injection of an adenovirus in which the COX-2 promoter drives a reporter gene. *Mol Imaging Biol*. Nov-Dec 2004;6(6):395-404.
 144. Lillehammer T, Tveito S, Engesaeter BO, Fodstad O, Maelandsmo GM, Engebraaten O. Melanoma-specific expression in first-generation adenoviral vectors in vitro and in vivo [mdash] use of the human tyrosinase promoter with human enhancers. *Cancer Gene Ther*. 2005;12(11):864-872.
 145. Limberis MP, Bell CL, Heath J, Wilson JM. Activation of transgene-specific T cells following lentivirus-mediated gene delivery to mouse lung. *Mol Ther*. Jan 2010;18(1):143-150.
 146. Lin J, Lalani AS, Harding TC, Gonzalez M, Wu W-W, Luan B, Tu GH, Koprivnikar K, VanRoey MJ, He Y, Alitalo K, Jooss K. Inhibition of Lymphogenous Metastasis Using Adeno-Associated Virus-Mediated Gene Transfer of a Soluble VEGFR-3 Decoy Receptor. *Cancer Res*. August 1, 2005 2005;65(15):6901-6909.
 147. Lin Y, Molter J, Lee Z, Gerson SL. Bioluminescence imaging of hematopoietic stem cell repopulation in murine models. *Methods in molecular biology (Clifton, N.J.)*. 2008;430:295-306.
 148. Lipshutz GS, Gruber CA, Cao Y, Hardy J, Contag CH, Gaensler KM. In utero delivery of adeno-associated viral vectors: intraperitoneal gene transfer produces long-term expression. *Molecular Therapy*. March 2001;3(3):284-292.
 149. Lipshutz GS, Titre D, Brindle M, Bisconte AR, Contag CH, Gaensler KM. Comparison of gene expression after intraperitoneal delivery of AAV2 or AAV5 in utero. *Mol Ther*. July 2003;8(1):90-98.
 150. Liu F, Sag D, Wang J, Shollenberger LM, Niu F, Yuan X, Li SD, Thompson M, Monahan P. Sine-wave current for efficient and safe in vivo gene transfer. *Mol Ther*. Oct 2007;15(10):1842-1847.
 151. Liu S, Wang H, Yang Z, Kon T, Zhu J, Cao Y, Li F, Kirkpatrick J, Nicchitta CV, Li C-Y. Enhancement of Cancer Radiation Therapy by Use of Adenovirus-Mediated Secretable Glucose-Regulated Protein 94/gp96 Expression. *Cancer Res*. October 15, 2005 2005;65(20):9126-9131.
 152. Liu TC, Thorne SH, Kim DH. Oncolytic adenoviruses for cancer gene therapy. *Methods in molecular biology (Clifton, N.J.)*. 2008;433:243-258.
 153. Liu Y, Kon T, Li C, Zhong P. High intensity focused ultrasound-induced gene activation in solid tumors. *The Journal of the Acoustical Society of America*. Jul 2006;120(1):492-501.

154. Lobbestael E, Reumers V, Ibrahim A, Paesen K, Thiry I, Gijsbers R, Van den Haute C, Debyser Z, Baekelandt V, Taymans JM. Immunohistochemical detection of transgene expression in the brain using small epitope tags. *BMC biotechnology*. 2010;10:16.
155. Lu Y, Dang H, Middleton B, Zhang Z, Washburn L, Campbell-Thompson M, Atkinson MA, Gambhir SS, Tian J, Kaufman DL. Bioluminescent monitoring of islet graft survival after transplantation. *Mol Ther*. March 2004;9(3):428-435.
156. Marti GP, Mohebi P, Liu L, Wang J, Miyashita T, Harmon JW. KGF-1 for wound healing in animal models. *Methods in molecular biology (Clifton, N.J.)*. 2008;423:383-391.
157. Matis S, Mariani MR, Cutrona G, Cilli M, Piccardi F, Daga A, Damonte G, Millo E, Moroni M, Roncella S, Fedeli F, Boffa LC, Ferrarini M. PNAEmu can significantly reduce Burkitt's lymphoma tumor burden in a SCID mice model: cells dissemination similar to the human disease. *Cancer Gene Ther*. Apr 10 2009.
158. Matthews QL, Sibley DA, Wu H, Li J, Stoff-Khalili MA, Waehler R, Mathis JM, Curiel DT. Genetic incorporation of a herpes simplex virus type 1 thymidine kinase and firefly luciferase fusion into the adenovirus protein IX for functional display on the virion. *Mol Imaging*. Oct-Dec 2006;5(4):510-519.
159. Mayer-Kuckuk P, Gade TP, Buchanan IM, Doubrovin M, Ageyeva L, Bertino JR, Boskey AL, Blasberg RG, Koutcher JA, Banerjee D. High-resolution imaging of bone precursor cells within the intact bone marrow cavity of living mice. *Mol Ther*. Jul 2005;12(1):33-41.
160. Mays LE, Vandenberghe LH, Xiao R, Bell P, Nam HJ, Agbandje-McKenna M, Wilson JM. Adeno-associated virus capsid structure drives CD4-dependent CD8+ T cell response to vector encoded proteins. *J Immunol*. May 15 2009;182(10):6051-6060.
161. McCaffrey A, Kay MA, Contag CH. Advancing molecular therapies through in vivo bioluminescent imaging. *Mol Imaging*. April 2003;2(2):75-86.
162. McCaffrey AP, Meuse L, Karimi M, Contag CH, Kay MA. A potent and specific morpholino antisense inhibitor of hepatitis C translation in mice. *Hepatology*. August 2003;38(2):503-508.
163. McCaffrey AP, Ohashi K, Meuse L, Shen S, Lancaster AM, Lukavsky PJ, Sarnow P, Kay MA. Determinants of hepatitis C translational initiation in vitro, in cultured cells and mice. *Mol Ther*. June 2002;5(6):676-684.
164. Mi Z, Guo H, Russell MB, Liu Y, Sullenger BA, Kuo PC. RNA aptamer blockade of osteopontin inhibits growth and metastasis of MDA-MB231 breast cancer cells. *Mol Ther*. Jan 2009;17(1):153-161.
165. Mizuno T, Mohri K, Nasu S, Danjo K, Okamoto H. Dual imaging of pulmonary delivery and gene expression of dry powder inhalant by fluorescence and bioluminescence. *J Control Release*. Nov 30 2008.
166. Mocanu J, Yip K, Alajez N, Shi W, Li JH, Lunt S, Moriyama E, Wilson B, Milosevic M, Lo K, van Rooijen N, Busson P, Bastianutto C, Liu FF. Imaging the Modulation of Adenoviral Kinetics and Biodistribution for Cancer Gene Therapy. *Mol Ther*. Mar 13 2007.
167. Mocanu JD, Moriyama EH, Chia MC, Li JH, Yip KW, Huang DP, Bastianutto C, Wilson BC, Liu FF. Combined in vivo bioluminescence and fluorescence imaging for cancer gene therapy. *Mol Imaging*. October 2004;3(4):352-355.
168. Mohri K, Okuda T, Mori A, Danjo K, Okamoto H. Optimized pulmonary gene transfection in mice by spray-freeze dried powder inhalation. *J Control Release*. Jun 1 2010;144(2):221-226.
169. Monfared P, Winkeler A, Klein M, Li H, Klose A, Hoesel M, Waerzeggers Y, Korsching S, Jacobs AH. Noninvasive assessment of E2F-1-mediated transcriptional regulation in vivo. *Cancer Res*. Jul 15 2008;68(14):5932-5940.
170. Morrison J, Briggs SS, Green NK, Thoma C, Fisher KD, Kehoe S, Seymour LW. Cetuximab retargeting of adenovirus via the epidermal growth factor receptor for treatment of intraperitoneal ovarian cancer. *Hum Gene Ther*. Mar 2009;20(3):239-251.
171. Munguia A, Ota T, Miest T, Russell SJ. Cell carriers to deliver oncolytic viruses to sites of myeloma tumor growth. *Gene therapy*. May 2008;15(10):797-806.
172. Nakajima A, Seroogy CM, Sandora MR, Tarner IH, Costa GL, Taylor-Edwards C,

- Bachmann MH, Contag CH, Fathman CG. Antigen-specific T cell-mediated gene therapy in collagen-induced arthritis. *J Clin Invest*. May 2001;107(10):1293-1301.
173. Nguyen AT, Dow AC, Kupiec-Weglinski J, Busuttil RW, Lipshutz GS. Evaluation of gene promoters for liver expression by hydrodynamic gene transfer. *The Journal of surgical research*. Jul 2008;148(1):60-66.
174. Nguyen VH, Kim HS, Ha JM, Hong Y, Choy HE, Min JJ. Genetically engineered *Salmonella typhimurium* as an imageable therapeutic probe for cancer. *Cancer Res*. Jan 1 2010;70(1):18-23.
175. Niu G, Xiong Z, Cheng Z, Cai W, Gambhir SS, Xing L, Chen X. In vivo bioluminescence tumor imaging of RGD peptide-modified adenoviral vector encoding firefly luciferase reporter gene. *Mol Imaging Biol*. May-Jun 2007;9(3):126-134.
176. Nyati MK, Symon Z, Kievit E, Dornfeld KJ, Rynkiewicz SD, Ross BD, Rehemtulla A, Lawrence TS. The potential of 5-fluorocytosine/cytosine deaminase enzyme prodrug gene therapy in an intrahepatic colon cancer model *Gene therapy*. July 2002;9(13):844-849.
177. Ogura Tsuyoshi HMJMSMTOTMMUAKHHHYYSKO. Utility of intraperitoneal administration as a route of AAV serotype 5 vector-mediated neonatal gene transfer. *The journal of gene medicine*. 2006;8(8):990-997.
178. Ohlfest JR, Demorest ZL, Motooka Y, Vengco I, Oh S, Chen E, Scappaticci FA, Saplis RJ, Ekker SC, Low WC. Combinatorial Antiangiogenic Gene Therapy by Nonviral Gene Transfer Using the Sleeping Beauty Transposon Causes Tumor Regression and Improves Survival in Mice Bearing Intracranial Human Glioblastoma. *Molecular Therapy*. 2005;12(5):778-788.
179. Ohlfest JR, Lobitz PD, Perkinson SG, Largaespada DA. Integration and long-term expression in xenografted human glioblastoma cells using a plasmid-based transposon system. *Mol Ther*. Aug 2004;10(2):260-268.
180. Ong HT, Trejo TR, Pham LD, Oberg AL, Russell SJ, Peng KW. Intravascularly administered RGD-displaying measles viruses bind to and infect neovessel endothelial cells in vivo. *Mol Ther*. Jun 2009;17(6):1012-1021.
181. Ortiz-Abalia J, Sahun I, Altafaj X, Andreu N, Estivill X, Dierssen M, Fillat C. Targeting Dyrk1A with AAVshRNA attenuates motor alterations in TgDyrk1A, a mouse model of Down syndrome. *American journal of human genetics*. Oct 2008;83(4):479-488.
182. Osborn MJ, McElmurry RT, Peacock B, Tolar J, Blazar BR. Targeting of the CNS in MPS-IH using a nonviral transferrin-alpha-L-iduronidase fusion gene product. *Mol Ther*. Aug 2008;16(8):1459-1466.
183. Paneda A, Vanrell L, Mauleon I, Crettaz JS, Berraondo P, Timmermans EJ, Beattie SG, Twisk J, van Deventer S, Prieto J, Fontanellas A, Rodriguez-Pena MS, Gonzalez-Asequinolaza G. Effect of adeno-associated virus serotype and genomic structure on liver transduction and biodistribution in mice of both genders. *Hum Gene Ther*. Aug 2009;20(8):908-917.
184. Park CY, Chuck RS, Cano M, Yew M, Nguyen V, Parker J, Mori K, Gehlbach PL. Periocular triamcinolone enhances intraocular gene expression after delivery by adenovirus. *Investigative ophthalmology & visual science*. Jan 2008;49(1):399-406.
185. Park JS, Na K, Woo DG, Yang HN, Kim JM, Kim JH, Chung HM, Park KH. Non-viral gene delivery of DNA polyplexed with nanoparticles transfected into human mesenchymal stem cells. *Biomaterials*. Jan 2010;31(1):124-132.
186. Parker AL, Waddington SN, Buckley SM, Custers J, Havenga MJ, van Rooijen N, Goudsmit J, McVey JH, Nicklin SA, Baker AH. Effect of neutralizing sera on factor x-mediated adenovirus serotype 5 gene transfer. *Journal of virology*. Jan 2009;83(1):479-483.
187. Passineau MJ, Zourelis L, Machen L, Edwards PC, Benza RL. Ultrasound-assisted non-viral gene transfer to the salivary glands. *Gene therapy*. May 27 2010.
188. Patil CS, Liu M, Zhao W, Coatney DD, Li F, VanTubergen EA, D'Silva NJ, Kirkwood KL. Targeting mRNA stability arrests inflammatory bone loss. *Mol Ther*. Oct 2008;16(10):1657-1664.
189. Peng X, Tao K, Cheng T, Zhu J, Zhang X. Efficient Inhibition of wear debris-induced inflammation by locally delivered siRNA. *Biochemical and biophysical research communications*. Dec 12 2008;377(2):532-537.

190. Peterson JR, Infanger DW, Braga VA, Zhang Y, Sharma RV, Engelhardt JF, Davisson RL. Longitudinal Non-invasive Monitoring of Transcription Factor Activation in Cardiovascular Regulatory Nuclei Using Bioluminescence Imaging. *Physiol Genomics*. Jan 29 2008.
191. Pichler A, Zelcer N, Prior JL, Kuil AJ, Piwnica-Worms D. In vivo RNA Interference-Mediated Ablation of MDR1 P-Glycoprotein. *Clin Cancer Res*. June 2005;11(12):4487-4494.
192. Podetz-Pedersen KM, Bell JB, Steele TW, Wilber A, Shier WT, Belur LR, McIvor RS, Hackett PB. Gene expression in lung and liver after intravenous infusion of polyethylenimine complexes of Sleeping Beauty transposons. *Hum Gene Ther*. Feb 2010;21(2):210-220.
193. Polyak B, Fishbein I, Chorny M, Alferiev I, Williams D, Yellen B, Friedman G, Levy RJ. High field gradient targeting of magnetic nanoparticle-loaded endothelial cells to the surfaces of steel stents. *Proc Natl Acad Sci U S A*. Jan 15 2008;105(2):698-703.
194. Ponomarev V, Doubrovin M, Serganova I, Vider J, Shavrin A, Beresten T, Ivanova A, Ageyeva L, Tourkova V, Balatoni J, Bornmann W, Blasberg R, Gelovani Tjuvajev J. A novel triple-modality reporter gene for whole-body fluorescent, bioluminescent, and nuclear noninvasive imaging. *European Journal of Nuclear Medicine and Molecular Imaging*. May 2004;31(5):740-751.
195. Power AT, Wang J, Falls TJ, Paterson JM, Parato KA, Lichty BD, Stojdl DF, Forsyth PA, Atkins H, Bell JC. Carrier Cell-based Delivery of an Oncolytic Virus Circumvents Antiviral Immunity. *Mol Ther*. Jan 2007;15(1):123-130.
196. Prasad KM, Xu Y, Yang Z, Toufektsian MC, Berr SS, French BA. Topoisomerase Inhibition Accelerates Gene Expression after Adeno-associated Virus-mediated Gene Transfer to the Mammalian Heart. *Mol Ther*. Apr 2007;15(4):764-771.
197. Qiu Z, Harms JS, Zhu J, Splitter GA. Bovine herpesvirus tegument protein VP22 enhances thymidine kinase/ganciclovir suicide gene therapy for neuroblastomas compared to herpes simplex virus VP22. *Journal of virology*. April 2004;78(8):4224-4233.
198. Rabinovich PM, Komarovskaya ME, Wrzesinski SH, Alderman JL, Budak-Alpdogan T, Karpikov A, Guo H, Flavell RA, Cheung NK, Weissman SM, Bahceci E. Chimeric Receptor mRNA Transfection as a Tool to Generate Antineoplastic Lymphocytes. *Hum Gene Ther*. Oct 15 2008.
199. Raki M, Hakkarainen T, Bauerschmitz GJ, Sarkioja M, Desmond RA, Kanerva A, Hemminki A. Utility of TK/GCV in the context of highly effective oncolysis mediated by a serotype 3 receptor targeted oncolytic adenovirus. *Gene therapy*. Oct 2007;14(19):1380-1388.
200. Ranki T, Kanerva A, Ristimaki A, Hakkarainen T, Sarkioja M, Kangasniemi L, Raki M, Laakkonen P, Goodison S, Hemminki A. A heparan sulfate-targeted conditionally replicative adenovirus, Ad5.pk7-Delta24, for the treatment of advanced breast cancer. *Gene therapy*. Jan 2007;14(1):58-67.
201. Ray P, De A, Min JJ, Tsien RY, Gambhir SS. Imaging tri-fusion multimodality reporter gene expression in living subjects. *Cancer Res*. February 2004;64(4):1323-1330.
202. Ray P, Wu AM, Gambhir SS. Optical bioluminescence and positron emission tomography imaging of a novel fusion reporter gene in tumor xenografts of living mice. *Cancer Res*. March 2003;63(6):1160-1165.
203. Ray S, Paulmurugan R, Hildebrandt I, Iyer M, Wu L, Carey M, Gambhir SS. Novel bidirectional vector strategy for amplification of therapeutic and reporter gene expression. *Hum Gene Ther*. July 2004;15(7):681-690.
204. Ray S, Paulmurugan R, Patel MR, Ahn BC, Wu L, Carey M, Gambhir SS. Noninvasive imaging of therapeutic gene expression using a bidirectional transcriptional amplification strategy. *Mol Ther*. Nov 2008;16(11):1848-1856.
205. Rehemtulla A, Hall DE, Stegman LD, Prasad U, Chen G, Bhojani MS, Chenevert TL, Ross BD. Molecular imaging of gene expression and efficacy following adenoviral-mediated brain tumor gene therapy. *Mol Imaging*. Jan-March 2002;1(1):43-55.
206. Rettig GR, McAnuff M, Liu D, Kim J-S, Rice KG. Quantitative bioluminescence imaging of transgene expression in vivo. *Analytical Biochemistry*. 2006;355(1):90-94.
207. Reumers V, Deroose CM, Krylyshkina O, Nuyts J, Geraerts M, Mortelmans L,

- Gijsbers R, Van den Haute C, Debyser Z, Baekelandt V. Noninvasive and quantitative monitoring of adult neuronal stem cell migration in mouse brain using bioluminescence imaging. *Stem cells (Dayton, Ohio)*. Sep 2008;26(9):2382-2390.
208. Rodrigo-Garzon M, Berraondo P, Ochoa L, Zulueta JJ, Gonzalez-Aseguinolaza G. Antitumoral efficacy of DNA nanoparticles in murine models of lung cancer and pulmonary metastasis. *Cancer Gene Ther*. Jul 3 2009.
209. Roos AK, Eriksson F, Walters DC, Pisa P, King AD. Optimization of Skin Electroporation in Mice to Increase Tolerability of DNA Vaccine Delivery to Patients. *Mol Ther*. Jun 16 2009.
210. Sato M, Johnson M, Zhang L, Gambhir SS, Carey M, Wu L. Functionality of Androgen Receptor-Based Gene Expression Imaging in Hormone Refractory Prostate Cancer. *Clin Cancer Res*. May 15, 2005 2005;11(10):3743-3749.
211. Sato M, Johnson M, Zhang L, Zhang B, Le K, Gambhir SS, Carey M, Wu L. Optimization of adenoviral vectors to direct highly amplified prostate-specific expression for imaging and gene therapy. *Mol Ther*. November 2003;8(5):726-737.
212. Sato Y, Ajiki T, Inoue S, Fujishiro J, Yoshino H, Igarashi Y, Hakamata Y, Kaneko T, Murakami T, Kobayashi E. Gene silencing in rat-liver and limb grafts by rapid injection of small interference RNA. *Transplantation*. Jan 27 2005;79(2):240-243.
213. Sato Y, Endo H, Ajiki T, Hakamata Y, Okada T, Murakami T, Kobayashi E. Establishment of Cre/LoxP recombination system in transgenic rats. *Biochemical and biophysical research communications*. Jul 9 2004;319(4):1197-1202.
214. Scheiman J, Tseng JC, Zheng Y, Meruelo D. Multiple Functions of the 37/67-kd Laminin Receptor Make It a Suitable Target for Novel Cancer Gene Therapy. *Mol Ther*. Sep 1 2009.
215. Schuster DP, Kovacs A, Garbow J, Piwnicka-Worms D. Recent advances in imaging the lungs of intact small animals. *Am J Respir Cell Mol Biol*. February 2004;30(2):129-138.
216. Schweichel D, Steitz J, Tormo D, Gaffal E, Ferrer A, Buchs S, Speuser P, Limmer A, Tuting T. Evaluation of DNA vaccination with recombinant adenoviruses using bioluminescence imaging of antigen expression: impact of application routes and delivery with dendritic cells. *The journal of gene medicine*. Oct 2006;8(10):1243-1250.
217. Schwerdt A, Zintchenko A, Concia M, Roesen N, Fisher KD, Lindner LH, Issels RD, Wagner E, Ogris M. Hyperthermia induced targeting of thermosensitive gene carriers to tumors. *Hum Gene Ther*. Aug 26 2008.
218. Sellins K, Fradkin L, Liggitt D, Dow S. Type I interferons potently suppress gene expression following gene delivery using liposome(-)DNA complexes. *Mol Ther*. Sep 2005;12(3):451-459.
219. Sher YP, Tzeng TF, Kan SF, Hsu J, Xie X, Han Z, Lin WC, Li LY, Hung MC. Cancer targeted gene therapy of BikDD inhibits orthotopic lung cancer growth and improves long-term survival. *Oncogene*. Jul 13 2009.
220. Shih CS, Laurie N, Holzmacher J, Spence Y, Nathwani AC, Davidoff AM, Dyer MA. AAV-mediated local delivery of interferon-beta for the treatment of retinoblastoma in preclinical models. *Neuromolecular medicine*. 2009;11(1):43-52.
221. Shin S, Shea LD. Lentivirus immobilization to nanoparticles for enhanced and localized delivery from hydrogels. *Mol Ther*. Apr 2010;18(4):700-706.
222. Sims TL, Jr., Hamner JB, Bush RA, Williams RF, Zhou J, Kim SU, Aboody KS, Danks MK, Davidoff AM. Neural progenitor cell-mediated delivery of interferon beta improves neuroblastoma response to cyclophosphamide. *Annals of surgical oncology*. Nov 2008;15(11):3259-3267.
223. Sinn PL, Arias AC, Brogden KA, McCray PB, Jr. Lentivirus vector can be readministered to nasal epithelia without blocking immune responses. *Journal of virology*. Nov 2008;82(21):10684-10692.
224. Slavin AJ, Turner IH, Nakajima A, Urbanek-Ruiz I, McBride J, Contag CH, Fathman CG. Adoptive cellular gene therapy of autoimmune disease. *Autoimmun Rev*. August 2002;1(4):213-219.
225. Streck CJ, Dickson PV, Ng CYC, Zhou J, Gray JT, Nathwani AC, Davidoff AM. Adeno-Associated Virus Vector-Mediated Systemic Delivery of IFN- β Combined with Low-Dose Cyclophosphamide Affects Tumor Regression in Murine

- Neuroblastoma Models. *Clin Cancer Res.* August 15, 2005 2005;11(16):6020-6029.
226. Su YC, Chuang KH, Wang YM, Cheng CM, Lin SR, Wang JY, Hwang JJ, Chen BM, Chen KC, Roffler S, Cheng TL. Gene expression imaging by enzymatic catalysis of a fluorescent probe via membrane-anchored beta-glucuronidase. *Gene therapy.* Jan 18 2007.
227. Sun J, Hakobyan N, Valentino LA, Feldman BL, Samulski RJ, Monahan PE. Intraarticular factor IX protein or gene replacement protects against development of hemophilic synovitis in the absence of circulating factor IX. *Blood.* Dec 1 2008;112(12):4532-4541.
228. Sun W, Ziady AG. Real-time imaging of gene delivery and expression with DNA nanoparticle technologies. *Methods in molecular biology (Clifton, N.J.)* 2009;544:525-546.
229. Suzuki R, Takizawa T, Negishi Y, Utoguchi N, Sawamura K, Tanaka K, Namai E, Oda Y, Matsumura Y, Maruyama K. Tumor specific ultrasound enhanced gene transfer in vivo with novel liposomal bubbles. *J Control Release.* Jan 22 2008;125(2):137-144.
230. Szentirmai O, Baker CH, Bullain SS, Lin N, Takahashi M, Folkman J, Mulligan RC, Carter BS. Successful inhibition of intracranial human glioblastoma multiforme xenograft growth via systemic adenoviral delivery of soluble endostatin and soluble vascular endothelial growth factor receptor-2: laboratory investigation. *Journal of neurosurgery.* May 2008;108(5):979-988.
231. Takeshita F, Minakuchi Y, Nagahara S, Honma K, Sasaki H, Hirai K, Teratani T, Namatame N, Yamamoto Y, Hanai K, Kato T, Sano A, Ochiya T. Efficient delivery of small interfering RNA to bone-metastatic tumors by using atelocollagen in vivo. *Proc Natl Acad Sci U S A.* Aug 23 2005;102(34):12177-12182.
232. Tarantal AF, Lee C. Long-term Luciferase Expression Monitored by Bioluminescence Imaging after AAV-Mediated Fetal Gene Delivery in Rhesus Monkeys (*Macaca mulatta*). *Hum Gene Ther.* Sep 14 2009.
233. Tarantal AF, Lee CC, Jimenez DF, Cherry SR. Fetal Gene Transfer Using Lentiviral Vectors: In Vivo Detection of Gene Expression by microPET and Optical Imaging in Fetal and Infant Monkeys. *Hum Gene Ther.* Nov 30 2006.
234. Tarner IH, Slavin AJ, McBride J, Levicnik A, Smith R, Nolan GP, Contag CH, Fathman CG. Treatment of autoimmune disease by adoptive cellular gene therapy. *Ann N Y Acad Sci.* September 2003;998:512-519.
235. Taubes G. Firefly gene lights up lab animals from inside out. *Science.* June 1997;276(5321):1993.
236. Taymans JM, Vandenberghe LH, Haute CV, Thiry I, Deroose CM, Mortelmans L, Wilson JM, Debyser Z, Baekelandt V. Comparative analysis of adeno-associated viral vector serotypes 1, 2, 5, 7, and 8 in mouse brain. *Hum Gene Ther.* Mar 2007;18(3):195-206.
237. Thorne SH, Contag CH. Integrating the biological characteristics of oncolytic viruses and immune cells can optimize therapeutic benefits of cell-based delivery. *Gene therapy.* Mar 20 2008.
238. Toelen J, Deroose CM, Gijssbers R, Reumers V, Sbragia LN, Vets S, Chitneni SK, Bormans G, Mortelmans L, Deprest JA, Debyser Z. Fetal gene transfer with lentiviral vectors: long-term in vivo follow-up evaluation in a rat model. *American journal of obstetrics and gynecology.* Apr 2007;196(4):352 e351-356.
239. Tolar J, Osborn M, Bell S, McElmurry R, Xia L, Riddle M, Panoskaltsis-Mortari A, Jiang Y, Mclvor RS, Contag CH. Real-Time in Vivo Imaging of Stem Cells Following Transgenesis by Transposition. *Molecular Therapy.* July 2005;12(1):42-48.
240. Tolar J, Osborn M, Bell S, McElmurry R, Xia L, Riddle M, Panoskaltsis-Mortari A, Jiang Y, Mclvor RS, Contag CH, Yant SR, Kay MA, Verfaillie CM, Blazar BR. Real-time in vivo imaging of stem cells following transgenesis by transposition. *Mol Ther.* Jul 2005;12(1):42-48.
241. Tseng JC, Hurtado A, Yee H, Levin B, Boivin C, Benet M, Blank SV, Pellicer A, Meruelo D. Using sindbis viral vectors for specific detection and suppression of advanced ovarian cancer in animal models. *Cancer Res.* Sep 15 2004;64(18):6684-6692.

242. Tseng JC, Levin B, Hurtado A, Yee H, Perez de Castro I, Jimenez M, Shamamian P, Jin R, Novick RP, Pellicer A, Meruelo D. Systemic tumor targeting and killing by Sindbis viral vectors. *Nat Biotechnol*. January 2004;22(1):70-77.
243. Tseng JC, Zheng Y, Yee H, Levy DE, Meruelo D. Restricted tissue tropism and acquired resistance to Sindbis viral vector expression in the absence of innate and adaptive immunity. *Gene therapy*. Aug 2007;14(15):1166-1174.
244. Tysome JR, Briat A, Alusi G, Cao F, Gao D, Yu J, Wang P, Yang S, Dong Z, Wang S, Deng L, Francis J, Timiryasova T, Fodor I, Lemoine NR, Wang Y. Lister strain of vaccinia virus armed with endostatin-angiostatin fusion gene as a novel therapeutic agent for human pancreatic cancer. *Gene therapy*. Jul 9 2009.
245. Uchibori R, Okada T, Ito T, Urabe M, Mizukami H, Kume A, Ozawa K. Retroviral vector-producing mesenchymal stem cells for targeted suicide cancer gene therapy. *The journal of gene medicine*. May 2009;11(5):373-381.
246. van Pijkeren JP, Morrissey D, Monk IR, Cronin M, Rajendran S, O'Sullivan GC, Gahan CG, Tangney M. A novel *Listeria monocytogenes*-based DNA delivery system for cancer gene therapy. *Hum Gene Ther*. Apr 2010;21(4):405-416.
247. Vera J, Savoldo B, Vigouroux S, Biagi E, Pule M, Rossig C, Wu J, Heslop HE, Rooney CM, Brenner MK, Dotti G. T lymphocytes redirected against the kappa light chain of human immunoglobulin efficiently kill mature B lymphocyte-derived malignant cells. *Blood*. Dec 1 2006;108(12):3890-3897.
248. Vera M, Sobrevals L, Zaratiegui M, Martinez L, Palencia B, Rodriguez CM, Prieto J, Fortes P. Liver transduction with a simian virus 40 vector encoding insulin-like growth factor I reduces hepatic damage and the development of liver cirrhosis. *Gene therapy*. Feb 2007;14(3):203-210.
249. Wang Y, Iyer M, Annala AJ, Chappell S, Mauro V, Gambhir SS. Noninvasive monitoring of target gene expression by imaging reporter gene expression in living animals using improved bicistronic vectors. *J Nucl Med*. April 2005;46(4):667-674.
250. Wang Y, Liu S, Li C-Y, Yuan F. A Novel Method for Viral Gene Delivery in Solid Tumors. *Cancer Res*. September 1, 2005 2005;65(17):7541-7545.
251. Wang Y, Yang Z, Liu S, Kon T, Krol A, Li CY, Yuan F. Characterisation of systemic dissemination of nonreplicating adenoviral vectors from tumours in local gene delivery. *Br J Cancer*. Apr 25 2005;92(8):1414-1420.
252. Wangensteen KJ, Wilber A, Keng VW, He Z, Matise I, Wangensteen L, Carson CM, Chen Y, Steer CJ, Mclvor RS, Largaespada DA, Wang X, Ekker SC. A facile method for somatic, lifelong manipulation of multiple genes in the mouse liver. *Hepatology*. May 2008;47(5):1714-1724.
253. Wilber A, Frandsen JL, Wangensteen KJ, Ekker SC, Wang X, Mclvor RS. Dynamic Gene Expression After Systemic Delivery of Plasmid DNA as Determined by In Vivo Bioluminescence Imaging. *Human Gene Therapy*. 2005;16(11):1325-1332.
254. Wilber A, Linehan JL, Tian X, Woll PS, Morris JK, Belur LR, Mclvor RS, Kaufman DS. Efficient and stable transgene expression in human embryonic stem cells using transposon-mediated gene transfer. *Stem cells (Dayton, Ohio)*. Nov 2007;25(11):2919-2927.
255. Wilber A, Wangensteen KJ, Chen Y, Zhuo L, Frandsen JL, Bell JB, Chen ZJ, Ekker SC, Mclvor RS, Wang X. Messenger RNA as a Source of Transposase for Sleeping Beauty Transposon-mediated Correction of Hereditary Tyrosinemia Type I. *Mol Ther*. Apr 17 2007.
256. Wilson AA, Murphy GJ, Hamakawa H, Kwok LW, Srinivasan S, Hovav AH, Mulligan RC, Amar S, Suki B, Kotton DN. Amelioration of emphysema in mice through lentiviral transduction of long-lived pulmonary alveolar macrophages. *J Clin Invest*. Jan 2010;120(1):379-389.
257. Winkeler A, Sena-Esteves M, Paulis LE, Li H, Waerzeggers Y, Ruckriem B, Himmelreich U, Klein M, Monfared P, Rueger MA, Heneka M, Vollmar S, Hoehn M, Fraefel C, Graf R, Wienhard K, Heiss WD, Jacobs AH. Switching on the lights for gene therapy. *PLoS ONE*. 2007;2(6):e528.
258. Wu JC, Inubushi M, Sundaresan G, Schelbert HR, Gambhir SS. Optical imaging of cardiac reporter gene expression in living rats. *Circulation*. April 2002;105(14):1631-1634.

259. Wu JC, Sundaresan G, Iyer M, Gambhir SS. Noninvasive optical imaging of firefly luciferase reporter gene expression in skeletal muscles of living mice. *Mol Ther*. October 2001;4(4):297-306.
260. Xie X, Cao F, Sheikh AY, Li Z, Connolly AJ, Pei X, Li RK, Robbins RC, Wu JC. Genetic modification of embryonic stem cells with VEGF enhances cell survival and improves cardiac function. *Cloning and stem cells*. Winter 2007;9(4):549-563.
261. Xie X, Hsu JL, Choi MG, Xia W, Yamaguchi H, Chen CT, Trinh BQ, Lu Z, Ueno NT, Wolf JK, Bast RC, Jr., Hung MC. A novel hTERT promoter-driven E1A therapeutic for ovarian cancer. *Molecular cancer therapeutics*. Aug 2009;8(8):2375-2382.
262. Xing Y, Pua EC, Lu X, Zhong P. Low-amplitude ultrasound enhances hydrodynamic-based gene delivery to rat kidney. *Biochemical and biophysical research communications*. Aug 14 2009;386(1):217-222.
263. Yao H, Ng SS, Tucker WO, Tsang YK, Man K, Wang XM, Chow BK, Kung HF, Tang GP, Lin MC. The gene transfection efficiency of a folate-PEI600-cyclodextrin nanopolymer. *Biomaterials*. Oct 2009;30(29):5793-5803.
264. Yoshimitsu M, Sato T, Tao K, Walia JS, Rasaiah VI, Sleep GT, Murray GJ, Poepl AG, Underwood J, West L, Brady RO, Medin JA. Bioluminescent imaging of a marking transgene and correction of Fabry mice by neonatal injection of recombinant lentiviral vectors. *Proc Natl Acad Sci U S A*. November 2004;101(48):16909-16914.
265. Yurek DM, Fletcher AM, Smith GM, Seroogy KB, Ziady AG, Molter J, Kowalczyk TH, Padegimas L, Cooper MJ. Long-term transgene expression in the central nervous system using DNA nanoparticles. *Mol Ther*. Apr 2009;17(4):641-650.
266. Zabala M, Alzuguren P, Benavides C, Crettaz J, Gonzalez-Aseguinolaza G, Ortiz de Solorzano C, Gonzalez-Aparicio M, Kramer MG, Prieto J, Hernandez-Alcoceba R. Evaluation of bioluminescent imaging for noninvasive monitoring of colorectal cancer progression in the liver and its response to immunogene therapy. *Molecular cancer*. 2009;8:2.
267. Zhang HY, Man JH, Liang B, Zhou T, Wang CH, Li T, Li HY, Li WH, Jin BF, Zhang PJ, Zhao J, Pan X, He K, Gong WL, Zhang XM, Li AL. Tumor-targeted delivery of biologically active TRAIL protein. *Cancer Gene Ther*. Jan 15 2010.
268. Zhang KX, Moussavi M, Kim C, Chow E, Chen IS, Fazli L, Jia W, Rennie PS. Lentiviruses with trastuzumab bound to their envelopes can target and kill prostate cancer cells. *Cancer Gene Ther*. Apr 17 2009.
269. Zhang YQ, Tsai YC, Monie A, Wu TC, Hung CF. Enhancing the Therapeutic Effect Against Ovarian Cancer Through a Combination of Viral Oncolysis and Antigen-specific Immunotherapy. *Mol Ther*. Jan 19 2010.
270. Zincarelli C, Soltys S, Rengo G, Rabinowitz JE. Analysis of AAV serotypes 1-9 mediated gene expression and tropism in mice after systemic injection. *Mol Ther*. Jun 2008;16(6):1073-1080.