


## 小动物活体光学成像技术在免疫学研究中的应用

Revvity小动物活体光学成像技术已在生命科学基础研究、临床前医学研究及药物研发等领域得到广泛应用。在众多应用领域中，免疫研究是活体光学成像技术的应用热点之一。在应用活体光学成像技术进行免疫学研究中，常用的标记方法及应用领域包括：1、利用功能性探针监测免疫疾病的发生发展及相关治疗；2、利用荧光素酶基因或荧光染料标记免疫细胞，监测免疫细胞的免疫应答作用；3、利用荧光素酶作为报告基因标记疾病相关基因构建转基因动物，进行免疫疾病机理研究。下面结合一些具体实例进行阐述：

### 一.监测免疫疾病的发生发展及治疗效果


伴随免疫疾病研究的深入，目前已开发出一系列针对免疫疾病监测的功能性探针，这些探针的设计大多是基于在免疫疾病中表达的特征性分子或酶，通过对特征性分子或酶的监测而反映疾病的发生发展。利用这些探针并结合活体光学成像技术，研究者可以方便快捷地在活体动物水平监测免疫疾病的发生发展及治疗效果。

如在各种炎性疾病的发生发展中，通常会伴随大量活性氧自由基（ROS）及氧化胁迫的产生。研究者根据上述特征，设计出能够特异性探测 ROS 或氧化胁迫相关酶的功能性探针用于炎性疾病的检测。Lee 等应用其自行研发的表面结合透明质酸的金纳米颗粒（HA-AuNPs）成功检测了小鼠关节炎的发生。该纳米颗粒表面结合了荧光染料标记的透明质酸，当未被激活时，由于荧光染料及纳米颗粒本身的相互作用而处于荧光湮灭状态，而表面结合的透明质酸一旦被 ROS 或透明质酸酶剪切，便会被外界光源激发而发光。应用该探针并结合活体光学成像技术，便可在活体动物水平灵敏监测到炎性疾病的发生。




上图：应用 HA-AuNPs 及 IVIS 系统观测小鼠四肢关节炎的发生。尾静脉注射探针，不同时间点观测。

Revvity依靠强大的生物学研发团队，成功开发出多种应用于炎性疾病监测的功能性探针，如用于探测炎性细胞中髓过氧化物酶（MPO）的化学发光（Chemiluminescent）探针“XenoLight® RediJect™ Inflammation Probe”、用于探测环氧酶2（COX-2）的近红外荧光探针“XenoLight® RediJect™ COX-2 Probe”、用于探测炎性细胞中组织蛋白酶（Cathepsin）的近红外探针“ProSense 680/750”、用于探测炎性细胞中基质金属蛋白酶（Matrix Metalloproteinase）的近红外探针“MMPsense 680/750”、用于探测嗜中性粒细胞胰肽酶（Neutrophil Elastase）的近红外探针“Neutrophil Elastase 680 FAST”。这些探针已被广泛应用于各种炎性疾病的研究中。



上图：应用 IVIS 成像系统及功能性探针监测小鼠关节炎的发生。A、利用 XenoLight® RediJect™ COX-2 探针观测关节炎发生过程中 COX-2 的表达；B、利用 XenoLight® RediJect™ Inflammation Probe 观测关节炎发生过程中 MPO 的表达。




上图：应用 FMT 成像系统及功能性探针监测肺炎的发生及治疗效果。A、利用 ProSense 750 探针观测哮喘中组织蛋白酶的表达及地塞米松（dexamethasone）对哮喘的治疗效果；B、利用 Neutrophil Elastase 680 FAST 探针观测急性肺炎（COPD）的发生及 rolipram 的治疗效果。

## 二.监测免疫细胞的免疫应答

免疫应答是机体免疫系统对抗原刺激所产生的以排除抗原为目的的生理过程。这个过程是免疫系统各部分生理功能的综合体现，包括了抗原递呈、淋巴细胞活化、免疫分子形成及免疫效应发生等一系列的生理反应。通过有效的免疫应答，机体得以维护内环境的稳定。免疫细胞在机体的免疫应答中发挥着重要作用，了解免疫细胞的作用机理是免疫学研究的重要环节。活体光学成像技术已广泛应用于免疫细胞的相关研究，通过该技术可以在活体动物水平监测免疫细胞在相关疾病中的迁移、分布及功能。目前用于标记免疫细胞的主要方法包括：1、通过带有荧光素酶基因或荧光蛋白基因的病毒载体稳定转染人源或鼠源免疫细胞，使免疫细胞具有发光性质；2、直接从转基因发光小鼠中提取免疫细胞，所获得的免疫细胞即有发光性质；3、通过特定的荧光染料直接标记免疫细胞使其具有发光性质。研究者可根据具体研究，选择合适的标记方法对免疫细胞进行标记。


利用带有荧光素酶基因的病毒载体稳定转染免疫细胞的好处在于：1、生物发光灵敏度高，可以满足少量被标记免疫细胞的活体观测；2、由于是稳定转染，因此发光特性可以传代，不会因为细胞的分裂或分化而改变，可以进行数周或更长时间的观测。

如 Rabinovich 等人利用增强型萤火虫荧光素酶基因（*effluc*）标记 T 细胞，发现其发光强度比传统萤火虫荧光素酶（*ffluc*）高出至少 100 倍，将标记的 T 细胞皮下移植后，利用 IVIS 系统能够非常灵敏的观测到 3 个细胞发出的光信号。



上图：利用 IVIS 系统活体观测经荧光素酶基因标记的 T 细胞。A、分别利用 *ffluc* 或 *effluc* 标记 T 细胞的对比成像结果；B、利用 *effluc* 标记 T 细胞皮下移植成像结果。


Santos 等利用桡足类动物膜结合型荧光素酶基因 (ext GLuc) 标记 T 细胞, 发现其发光强度较萤火虫荧光素酶 (ffluc) 及海肾荧光素酶 (Rluc) 均有大幅提升。另外, 研究者将 ext GLuc 标记的 T 细胞尾静脉注射入携带淋巴瘤 A20 (OVA) 的 SCID 小鼠体内, 应用 IVIS 系统对 T 细胞在体内的分布进行长期观测, 发现 T 细胞能够有效靶向肿瘤, 并对肿瘤细胞具有杀伤作用。



上图: 利用 IVIS 系统观测 ext GLuc 标记的 T 细胞。上行: 体外及体内成像结果显示经 ext GLuc 标记的 T 细胞发光强度更高; 下行: 长期观测经 ext GLuc 标记的 T 细胞在体内的分布及对肿瘤的靶向治疗。


通过上述带有荧光素酶基因的病毒载体在体外转染免疫细胞, 需要一定的实验基础及操作成熟度, 并且与肿瘤细胞相比, 免疫细胞的转染效率较低, 没有成熟的实验操作经验较难获取荧光素酶高表达的转染细胞株。而通过从组成型 (ubiquitous) 表达或特异性 (inducible) 表达荧光素酶基因的转基因发光小鼠中提取免疫细胞的方法相对方便, 无需自行进行标记, 用提取的具有发光性质的免疫细胞即可开展实验。目前市场上已有多种商业化的发光转基因小鼠供研究者选择购买, 如用组成型表达启动子控制荧光素酶基因表达而构建的转基因小鼠 Tg( $\beta$ -actin-luc)和 Tg(CMV-luc), 其全身均有荧光素酶的表达, 因此从脾、淋巴结、胸腺、骨髓等提取的免疫细胞即具备发光性质。

Cheeran 等利用从转基因小鼠 Tg( $\beta$ -actin-luc)中提取的脾细胞及淋巴结细胞, 研究了免疫细胞对病毒感染的响应。研究者将提取的发光脾细胞及淋巴结细胞通过尾静脉注入脑室内感染巨细胞病毒的小鼠, 利用 IVIS 系统观测了上述免疫细胞在活体动物体内对感染病灶点的浸润。结果显示, 在未经病毒感染的正常小鼠体内, 移植的淋巴细胞主要聚集于脾内 (如下图-1dpi 所示), 而当小鼠脑部感染病毒后, 这些淋巴细胞会迁移至感染区域而发挥免疫作用。



上图：利用 IVIS 系统观测免疫细胞对病毒感染的免疫应答。A、小鼠背部朝上拍摄；B、小鼠腹部朝上拍摄。-1dpi 为病毒感染前 24h 尾静脉注射淋巴细胞成像结果，每张图中从左至右第一只小鼠为注射生理盐水的对照小鼠，第二、三只小鼠为脑室内感染巨细胞病毒的疾病小鼠。


Chewning 等构建了利用 T 细胞特异性表达的人 CD2 基因启动子控制荧光素酶基因表达的表达载体，并应用该载体构建出转基因小鼠，从该转基因小鼠中提取的发光细胞即为 T 细胞。研究者将提取的发光 T 细胞经尾静脉注入由 OVA 诱导而引发炎症的小鼠体内，应用 IVIS 系统观测 T 细胞的免疫应答。



上图：利用 IVIS 系统观测 T 细胞对抗原的免疫应答。

除了上述两种通过荧光素酶基因对免疫细胞进行标记的方法之外，通过膜嵌合型荧光染料标记免疫细胞的方法则更为简便，研究者只需依据染料标记的实验流程，即可自行完成对免疫细胞的标记。目前市场上主要的膜嵌合型荧光染料包括 DiD、DiL、DiR 等，它们均是通过自身亲脂性的长碳链插入目标细胞的细胞膜中而进行标记。使用此种方式进行标记的好处在于操作相对简单，但由于是膜嵌合性标记，因此被标记细胞的发光性质不会传代，只能用于观测免疫细胞在体内的短时间动态变化。

之前的研究显示，乳腺癌的发生发展通常伴随有肿瘤周边的炎症发生，而炎性细胞在肿瘤周边炎症区域的聚集是炎症发生的必要条件。Sista 等利用 DiD 染料体外标记了单核细胞，并利用 IVIS 系统观测了尾静脉注射的单核细胞对原发性乳腺癌的靶向聚集情况。




上图：利用 IVIS 系统观测经 DiD 染料标记的单核细胞对肿瘤的靶向聚集。单核细胞经尾静脉注射后不同时间点的成像结果，所用小鼠为乳腺癌原发性模型小鼠 MMTV-PyMT。

### 三.免疫疾病机理研究

对于炎症相关基因的研究，可以揭示免疫疾病的分子机理，更好的了解免疫疾病的发生发展及相关治疗。小动物活体光学成像技术已越来越多的应用于此类研究。研究者通过构建各种生物发光转基因动物，结合活体光学成像技术，在活体动物水平观测免疫疾病发展过程中相关基因的表达。目前市场上已有不少成熟的可用于炎症研究的生物发光转基因小鼠，这些小鼠的构建通常是用炎症相关基因的启动子特异性控制荧光素酶基因的表达，从而通过发光情况反映炎症的发生及治疗效果。下图列举了几种可应用于相关免疫疾病研究的生物发光转基因小鼠模型，研究者可以从美国 Taconic 公司购买。


Disease	LPTA
<input type="checkbox"/> Rheumatoid Arthritis	NFkB RE, iNOS, SAA-1
<input type="checkbox"/> Neuro-inflammatory Disease	GFAP, SAA-1, COX-2
<input type="checkbox"/> Sepsis	NFkB-RE, IκB, SAA-1
<input type="checkbox"/> IBD	pS2/TFF-1, TNF-alpha
<input type="checkbox"/> Diabetes mellitus	RIP, Ins-2
<input type="checkbox"/> Pain	GFAP, Cox-2
<input type="checkbox"/> Asthma	Epx, NFkB RE

NFκB 是一类重要的转录因子，参与免疫反应的早期和炎症反应各阶段的许多分子都受 NF-κB 的调控，因此，对于 NFκB 信号通路的研究是炎症疾病研究的一个热点。Paur 等利用 Tg(NFκB-RE-luc)转基因小鼠，观测了咖啡对 LPS 诱导的炎症的抑制作用。结果显示，咖啡能够抑制 NFκB 的表达，并激活机体的抗氧化防御，而且随着咖啡烘焙程度的提高，这种抑制效应越明显。因此，利用 Tg(NFκB-RE-luc)转基因小鼠可以观测炎症的发生发展及相关药物的治疗效果。



上图：利用 IVIS 系统观测咖啡对 NFκB 表达的抑制。A、活体成像结果：上行为未给咖啡的对照组，下行为给咖啡组；B、不同烘焙程度咖啡对 NFκB 表达抑制的定量结果。


一氧化氮 (NO) 是细胞内的一种具有生物活性的重要调节分子，它既是是组织损伤的诱发因子和各种病变的增强因子，也是免疫应答的调节性分子。在炎性条件下，体内 NO 主要由炎性细胞中的诱导性一氧化氮合成酶 (iNOS) 产生，在败血症、关节炎、心血管疾病等多种炎性疾病中起调节作用。LPS 及多种细胞因子如干扰素 IFN-γ 等可激活 iNOS 的表达。下图所示为应用 Tg(iNos-luc)转基因小鼠观测不同抗炎药物对 iNOS 表达的抑制效应。研究者利用 LPS 诱导炎症的发生，对 iNOS 特异性表达的转基因小鼠施加不同抗炎药物，应用 IVIS 系统观测了不同抗炎药物对 iNOS 表达的抑制效果。



上图：利用 IVIS 系统及 Tg(iNos-luc)转基因小鼠观测不同抗炎药物对 iNOS 表达的抑制效果。


Cordeau 等利用 Tg(GFAP-luc)转基因小鼠，监测 GFAP 在肺炎链球菌感染而引发的脑膜炎中的表达。研究者用细菌荧光素酶 (bacterial luciferase) 基因标记肺炎链球菌，以监测细菌在活体动物体内的感染情况，同时以萤火虫荧光素酶基因标记的 Tg(GFAP-luc)转基因小鼠为实验动物，观测 GFAP 在肺炎链球菌感染而引发的脑膜炎中的表达，以及经抗生素治疗后细菌的感染情况和 GFAP 的表达情况。结果显示，随着细

菌感染程度及范围的升高和扩大，GFAP 的表达量也相应升高；而经抗生素治疗后，细菌的感染情况明显被抑制，GFAP 的表达量也随之降低。对小鼠感染脑部的体外成像结果显示，细菌对于脑部局部区域的感染，能引发整个脑部 GFAP 的大量表达。



上图：应用 IVIS 系统观测肺炎链球菌在 Tg(GFAP-luc)转基因小鼠中的感染情况及 GFAP 的表达情况。A. 活体成像结果，untreated：未经抗生素治疗，treated：经抗生素治疗；B. 脑部体外成像结果。

Luo 等利用 Tg(GFAP-luc)及 Tg(SBE-luc) 转基因小鼠，结合活体光学成像技术，观测了在自身免疫性脑脊髓炎（EAE）中与 GFAP 表达相关的星形胶质细胞的聚集及与 TGF- $\beta$  信号通路相关的炎症的发生。研究者用髓磷脂少突细胞糖蛋白（MOG）免疫小鼠，引发实验性自身免疫性脑脊髓炎，随后，利用 IVIS 系统观测患病 Tg(GFAP-luc)及 Tg(SBE-luc)转基因小鼠中 GFAP 及 TGF- $\beta$  的表达情况。结果显示，小鼠免疫后的第 7 天，即可观测到 GFAP 及 TGF- $\beta$  表达量显著升高，说明在脑脊髓炎的发病初期，即伴随有星形胶质细胞的聚集以及炎症的发生。值得注意的是，脑脊髓炎的明显临床症状出现于免疫后 11 天，因此，与观察临床症状而诊断疾病发生的方法相比，通过应用活体光学成像技术观测疾病相关基因的表达，能够更早的观测到疾病的发生。



上图：应用 IVIS 系统观测自身免疫性脑脊髓炎小鼠中 GFAP 及 TGF- $\beta$  的表达情况。（上）应用 Tg(GFAP-luc)转基因小鼠观测 GFAP 的表达；（下）应用 Tg(SBE-luc)转基因小鼠观测 TGF- $\beta$  的表达。



## IVIS APPLICATIONS IN IMMUNOLOGY RESEARCH

1. Abdel-Azim H, Zhu Y, Hollis R, Wang X, Ge S, Hao QL, Smbatyan G, Kohn DB, Rosol M, Crooks GM. Expansion of multipotent and lymphoid-committed human progenitors through intracellular dimerization of Mpl. *Blood*. Apr 15 2008;111(8):4064-4074.
2. Ahlen G, Weiland M, Derk E, Jiao J, Rahbin N, Aleman S, Peterson DL, Pokrovskaja K, Grander D, Frelin L, Sallberg M. Cleavage of the IPS-1/Cardif/MAVS/VISA does not inhibit T cell-mediated elimination of hepatitis C virus non-structural 3/4A-expressing hepatocytes. *Gut*. Aug 19 2008.
3. Ahmed N, Ratnayake M, Savoldo B, Perlaky L, Dotti G, Wels WS, Bhattacharjee MB, Gilbertson RJ, Shine HD, Weiss HL, Rooney CM, Heslop HE, Gottschalk S. Regression of experimental medulloblastoma following transfer of HER2-specific T cells. *Cancer Res*. Jun 15 2007;67(12):5957-5964.
4. Ahmed N, Salsman VS, Kew Y, Shaffer D, Powell S, Zhang YJ, Grossman RG, Heslop HE, Gottschalk S. HER2-specific T cells target primary glioblastoma stem cells and induce regression of autologous experimental tumors. *Clin Cancer Res*. Jan 15 2010;16(2):474-485.
5. Ahmed N, Salsman VS, Yvon E, Louis CU, Perlaky L, Wels WS, Dishop MK, Kleinerman EE, Pule M, Rooney CM, Heslop HE, Gottschalk S. Immunotherapy for osteosarcoma: genetic modification of T cells overcomes low levels of tumor antigen expression. *Mol Ther*. Oct 2009;17(10):1779-1787.
6. Alizadeh D, Zhang L, Brown CE, Farrukh O, Jensen MC, Badie B. Induction of anti-glioma natural killer cell response following multiple low-dose intracerebral CpG therapy. *Clin Cancer Res*. Jul 1 2010;16(13):3399-3408.
7. Amante FH, Stanley AC, Randall LM, Zhou Y, Haque A, McSweeney K, Waters AP, Janse CJ, Good MF, Hill GR, Engwerda CR. A role for natural regulatory T cells in the pathogenesis of experimental cerebral malaria. *Am J Pathol*. Aug 2007;171(2):548-559.
8. Asgharpour A, Gilchrist C, Baba D, Hamano S, Houghton E. Resistance to Intestinal *Entamoeba histolytica* Infection Is Conferred by Innate Immunity and Gr-1+ Cells. *Infect Immun*. August 1, 2005 2005;73(8):4522-4529.
9. Azadniv M, Dugger K, Bowers WJ, Weaver C, Crispe IN. Imaging CD8+ T cell dynamics in vivo using a transgenic luciferase reporter. *International immunology*. Oct 2007;19(10):1165-1173.
10. Backer R, Schwandt T, Greuter M, Oosting M, Jungerkes F, Tuting T, Boon L, O'Toole T, Kraal G, Limmer A, den Haan JM. Effective collaboration between marginal metallophilic macrophages and CD8+ dendritic cells in the generation of cytotoxic T cells. *Proc Natl Acad Sci*

- U S A. Jan 5 2010;107(1):216-221.
11. Bagnara D, Ibatici A, Corselli M, Sessarego N, Tenca C, De Santanna A, Mazzarello A, Daga A, Corvo R, De Rossi G, Frassoni F, Ciccone E, Fais F. Adoptive immunotherapy mediated by ex vivo expanded natural killer T cells against CD1d-expressing lymphoid neoplasms. *Haematologica*. Jul 2009;94(7):967-974.
  12. Bahey-EI-Din M, Casey PG, Griffin BT, Gahan CG. Lactococcus lactis-expressing listeriolysin O (LLO) provides protection and specific CD8(+) T cells against *Listeria monocytogenes* in the murine infection model. *Vaccine*. Sep 26 2008;26(41):5304-5314.
  13. Barberi T, Bradbury M, Dincer Z, Panagiotakos G, Socci ND, Studer L. Derivation of engraftable skeletal myoblasts from human embryonic stem cells. *Nature medicine*. May 2007;13(5):642-648.
  14. Beilhack A, Schulz S, Baker J, Beilhack GF, Wieland CB, Herman EI, Baker EM, Cao Y-A, Contag CH, Negrin RS. In vivo analyses of early events in acute graft-versus-host disease reveal sequential infiltration of T-cell subsets. *Blood*. August 2005;106(3):1113-1122.
  15. Bins AD, Jorritsma A, Wolkers MC, Hung CF, Wu TC, Schumacher TN, Haanen JB. A rapid and potent DNA vaccination strategy defined by in vivo monitoring of antigen expression. *Nature medicine*. Aug 2005;11(8):899-904.
  16. Biswas S, Guix M, Rinehart C, Dugger TC, Chytil A, Moses HL, Freeman ML, Arteaga CL. Inhibition of TGF-beta with neutralizing antibodies prevents radiation-induced acceleration of metastatic cancer progression. *J Clin Invest*. May 2007;117(5):1305-1313.
  17. BitMansour A, Burns SM, Traver D, Akashi K, Contag CH, Weissman IL, Brown JM. Myeloid progenitors protect against invasive aspergillosis and *Pseudomonas aeruginosa* infection following hematopoietic stem cell transplantation. *Blood*. December 2002;100(13):4660-4667.
  18. Bivas-Benita M, Bar L, Gillard GO, Kaufman DR, Simmons NL, Hovav AH, Letvin NL. Efficient generation of mucosal and systemic antigen-specific CD8+ T-cell responses following pulmonary DNA immunization. *Journal of virology*. Jun 2010;84(11):5764-5774.
  19. Bogdanov AA, Jr., Lin CP, Kang HW. Optical imaging of the adoptive transfer of human endothelial cells in mice using anti-human CD31 monoclonal antibody. *Pharmaceutical research*. Jun 2007;24(6):1186-1192.
  20. Brentjens RJ, Santos E, Nikhamin Y, Yeh R, Matsushita M, La Perle K, Quintas-Cardama A, Larson SM, Sadelain M. Genetically targeted T cells eradicate systemic acute lymphoblastic leukemia xenografts. *Clin Cancer Res*. Sep 15 2007;13(18 Pt 1):5426-5435.
  21. Brody JD, Goldstein MJ, Czerwinski DK, Levy R. Immunotransplantation preferentially expands T-effector cells over T-regulatory cells and cures large lymphoma tumors. *Blood*. Jan 1

- 2009;113(1):85-94.
22. Cao YA, Bachmann MH, Beilhack A, Yang Y, Tanaka M, Swijnenburg RJ, Reeves R, Taylor-Edwards C, Schulz S, Doyle TC, Fathman CG, Robbins RC, Herzenberg LA, Negrin RS, Contag CH. Molecular imaging using labeled donor tissues reveals patterns of engraftment, rejection, and survival in transplantation. *Transplantation*. Jul 15 2005;80(1):134-139.
  23. Cao YA, Wagers AJ, Beilhack A, Dusich J, Bachmann MH, Negrin RS, Weissman IL, Contag CH. Shifting foci of hematopoiesis during reconstitution from single stem cells. *Proceedings of the National Academy of Sciences of the United States of America*. January 2004;101(1):221-226.
  24. Chang CL, Wu TC, Hung CF. Control of human mesothelin-expressing tumors by DNA vaccines. *Gene therapy*. Aug 2007;14(16):1189-1198.
  25. Chatterjea D, Burns-Guydish SM, Sciuto TE, Dvorak A, Contag CH, Galli SJ. Adoptive transfer of mast cells does not enhance the impaired survival of KitW/KitW-v mice in a model of low dose intraperitoneal infection with bioluminescent *Salmonella typhimurium*. *Immunology Letters*. 2005;99(1):122-129.
  26. Chen CH, Yeh YC, Wu GJ, Huang YH, Lai WF, Liu JY, Tzeng CR. Tracking the rejection and survival of mouse ovarian iso- and allografts in vivo with bioluminescent imaging. *Reproduction (Cambridge, England)*. Jul 2010;140(1):105-112.
  27. Chen X, Zhang X, Larson C, Chen F, Kissler H, Kaufman DB. The epididymal fat pad as a transplant site for minimal islet mass. *Transplantation*. Jul 15 2007;84(1):122-125.
  28. Chen X, Zhang X, Larson C, Xia G, Kaufman DB. Prolonging islet allograft survival using in vivo bioluminescence imaging to guide timing of antilymphocyte serum treatment of rejection. *Transplantation*. May 15 2008;85(9):1246-1252.
  29. Chen X, Zhang X, Larson CS, Baker MS, Kaufman DB. In Vivo Bioluminescence Imaging of Transplanted Islets and Early Detection of Graft Rejection. [Article]. 2006.
  30. Choy G, Choyke P, Libutti SK. Current advances in molecular imaging: noninvasive in vivo bioluminescent and fluorescent optical imaging in cancer research. *Molecular Imaging*. October 2003;2(4):303-312.
  31. Chuang CM, Monie A, Wu A, Mao CP, Hung CF. Treatment with LL-37 peptide enhances antitumor effects induced by CpG oligodeoxynucleotides against ovarian cancer. *Human gene therapy*. Apr 2009;20(4):303-313.
  32. Chuang CM, Monie A, Wu A, Pai SI, Hung CF. Combination of viral oncolysis and tumor-specific immunity to control established tumors. *Clin Cancer Res*. Jul 15 2009;15(14):4581-4588.
  33. Contag CH, Bachmann MH. Advances in in vivo bioluminescence

- imaging of gene expression. *Annual Review of Biomedical Engineering*. 2002;4:235-260.
34. Contag PR. Bioluminescence imaging to evaluate infections and host response in vivo. *Methods in molecular biology (Clifton, N.J.)*. 2008;415:101-118.
  35. Cooper LJ, Al-Kadhimi Z, Serrano LM, Pfeiffer T, Olivares S, Castro A, Chang WC, Gonzalez S, Smith D, Forman SJ, Jensen MC. Enhanced antilymphoma efficacy of CD19-redirected influenza MP1-specific CTLs by cotransfer of T cells modified to present influenza MP1. *Blood*. February 2005;105(4):1622-1631.
  36. Costa GL, Sandora MR, Nakajima A, Nguyen EV, Taylor-Edwards C, Slavin AJ, Contag CH, Fathman CG, Benson JM. Adoptive immunotherapy of experimental autoimmune encephalomyelitis via T cell delivery of the IL-12 p40 subunit. *Journal of Immunology*. August 2001;167(4):2379-2387.
  37. Creusot RJ, Yaghoubi SS, Kodama K, Dang DN, Dang VH, Breckpot K, Thielemans K, Gambhir SS, Fathman CG. Tissue-targeted therapy of autoimmune diabetes using dendritic cells transduced to express IL-4 in NOD mice. *Clinical immunology (Orlando, Fla.)*. Mar 10 2008.
  38. Cubillos-Ruiz JR, Engle X, Scarlett UK, Martinez D, Barber A, Elgueta R, Wang L, Nesbeth Y, Durant Y, Gewirtz AT, Sentman CL, Kedl R, Conejo-Garcia JR. Polyethylenimine-based siRNA nanocomplexes reprogram tumor-associated dendritic cells via TLR5 to elicit therapeutic antitumor immunity. *J Clin Invest*. Aug 2009;119(8):2231-2244.
  39. de Mare A, Lambeck AJ, Regts J, van Dam GM, Nijman HW, Snippe H, Wilschut J, Daemen T. Viral vector-based prime-boost immunization regimens: a possible involvement of T-cell competition. *Gene therapy*. Mar 2008;15(6):393-403.
  40. Dobrenkov K, Olszewska M, Likar Y, Shenker L, Gunset G, Cai S, Pillarsetty N, Hricak H, Sadelain M, Ponomarev V. Monitoring the efficacy of adoptively transferred prostate cancer-targeted human T lymphocytes with PET and bioluminescence imaging. *J Nucl Med*. Jul 2008;49(7):1162-1170.
  41. Durai M, Krueger C, Ye Z, Cheng L, Mackensen A, Oelke M, Schneck JP. In vivo functional efficacy of tumor-specific T cells expanded using HLA-Ig based artificial antigen presenting cells (aAPC). *Cancer Immunol Immunother*. Feb 2009;58(2):209-220.
  42. Echchannaoui H, Frei K, Schnell C, Leib SL, Zimmerli W, Landmann R. Toll-like receptor 2-deficient mice are highly susceptible to *Streptococcus pneumoniae* meningitis because of reduced bacterial clearing and enhanced inflammation. *Journal of Infectious Disease*. September 2002;186(6):798-806.
  43. Edinger M, Cao YA, Verneris MR, Bachmann MH, Contag CH, Negrin RS. Revealing lymphoma growth and the efficacy of immune cell

- therapies using in vivo bioluminescence imaging. *Blood*. January 2003;101(2):640-648.
44. Elia L, Aurisicchio L, Facciabene A, Giannetti P, Ciliberto G, La Monica N, Palombo F. CD4+CD25+ regulatory T-cell-inactivation in combination with adenovirus vaccines enhances T-cell responses and protects mice from tumor challenge. *Cancer gene therapy*. Feb 2007;14(2):201-210.
  45. Elinav E, Waks T, Eshhar Z. Redirection of regulatory T cells with predetermined specificity for the treatment of experimental colitis in mice. *Gastroenterology*. Jun 2008;134(7):2014-2024.
  46. Fernandez I, Zeiser R, Karsunky H, Kambham N, Beilhack A, Soderstrom K, Negrin RS, Engleman E. CD101 Surface Expression Discriminates Potency Among Murine FoxP3+ Regulatory T Cells. *J Immunol*. Sep 1 2007;179(5):2808-2814.
  47. Fowler M, Virostko J, Chen Z, Poffenberger G, Radhika A, Brissova M, Shiota M, Nicholson WE, Shi Y, Hirshberg B, Harlan DM, Jansen ED, Powers AC. Assessment of pancreatic islet mass after islet transplantation using in vivo bioluminescence imaging. *Transplantation*. April 2005;79(7):768-776.
  48. Gade TPF, Hassen W, Santos E, Gunset G, Saudemont A, Gong MC, Brentjens R, Zhong X-S, Stephan M, Stefanski J, Lyddane C, Osborne JR, Buchanan IM, Hall SJ, Heston WD, Riviere I, Larson SM, Koutcher JA, Sadelain M. Targeted Elimination of Prostate Cancer by Genetically Directed Human T Lymphocytes. *Cancer Res*. October 1, 2005 2005;65(19):9080-9088.
  49. Gardner CL, Burke CW, Tesfay MZ, Glass PJ, Klimstra WB, Ryman KD. Eastern and Venezuelan equine encephalitis viruses differ in their ability to infect dendritic cells and macrophages: impact of altered cell tropism on pathogenesis. *Journal of virology*. Nov 2008;82(21):10634-10646.
  50. Gaymalov ZZ, Yang Z, Pisarev VM, Alakhov VY, Kabanov AV. The effect of the nonionic block copolymer pluronic P85 on gene expression in mouse muscle and antigen-presenting cells. *Biomaterials*. Feb 2009;30(6):1232-1245.
  51. Geiben-Lynn R, Greenland JR, Frimpong-Boateng K, van Rooijen N, Hovav AH, Letvin NL. CD4+ T lymphocytes mediate in vivo clearance of plasmid DNA vaccine antigen expression and potentiate CD8+ T-cell immune responses. *Blood*. Dec 1 2008;112(12):4585-4590.
  52. Gerdoni E, Gallo B, Casazza S, Musio S, Bonanni I, Pedemonte E, Mantegazza R, Frassoni F, Mancardi G, Pedotti R, Uccelli A. Mesenchymal stem cells effectively modulate pathogenic immune response in experimental autoimmune encephalomyelitis. *Annals of neurology*. Mar 2007;61(3):219-227.
  53. Giannoukakis N, Phillips B, Trucco M. Toward a cure for type 1 diabetes mellitus: diabetes-suppressive dendritic cells and beyond. *Pediatric diabetes*. Jun 2008;9(3 Pt 2):4-13.

54. Hakamata Y, Murakami T, Kobayashi E. "Firefly Rats" as an Organ/Cellular Source for Long-Term In Vivo Bioluminescent Imaging. [Article]. *Transplantation*. 2006.
55. Hardy J, Edinger M, Bachmann MH, Negrin RS, Fathman CG, Contag CH. Bioluminescence imaging of lymphocyte trafficking in vivo. *Exp Hematol*. December 2001;29(12):1353-1360.
56. Hemrajani C, Marches O, Wiles S, Girard F, Dennis A, Dziva F, Best A, Phillips AD, Berger CN, Mousnier A, Crepin VF, Kruidenier L, Woodward MJ, Stevens MP, La Ragione RM, MacDonald TT, Frankel G. Role of NleH, a type III secreted effector from attaching and effacing pathogens, in colonization of the bovine, ovine, and murine gut. *Infect Immun*. Nov 2008;76(11):4804-4813.
57. Highfill SL, Kelly RM, O'Shaughnessy MJ, Zhou Q, Xia L, Panoskaltis-Mortari A, Taylor PA, Tolar J, Blazar BR. Multipotent adult progenitor cells can suppress graft-versus-host disease via prostaglandin E2 synthesis and only if localized to sites of allopriming. *Blood*. Jul 16 2009;114(3):693-701.
58. Hildebrandt IJ, Gambhir SS. Molecular imaging applications for immunology. *Clinical immunology (Orlando, Fla)*. May 2004;111(2):210-224.
59. Hino A, Fukuyama S, Kataoka K, Kweon MN, Fujihashi K, Kiyono H. Nasal IL-12p70 DNA prevents and treats intestinal allergic diarrhea. *J Immunol*. Jun 1 2005;174(11):7423-7432.
60. Ho TY, Chen YS, Hsiang CY. Noninvasive nuclear factor-kappaB bioluminescence imaging for the assessment of host-biomaterial interaction in transgenic mice. *Biomaterials*. Oct 2007;28(30):4370-4377.
61. Huang P, Kaku H, Chen J, Kashiwakura Y, Saika T, Nasu Y, Urata Y, Fujiwara T, Watanabe M, Kumon H. Potent antitumor effects of combined therapy with a telomerase-specific, replication-competent adenovirus (OBP-301) and IL-2 in a mouse model of renal cell carcinoma. *Cancer gene therapy*. Feb 19 2010.
62. Huang X, Guo H, Kang J, Choi S, Zhou TC, Tammana S, Lees CJ, Li ZZ, Milone M, Levine BL, Tolar J, June CH, Scott Mclvor R, Wagner JE, Blazar BR, Zhou X. Sleeping Beauty transposon-mediated engineering of human primary T cells for therapy of CD19+ lymphoid malignancies. *Mol Ther*. Mar 2008;16(3):580-589.
63. Huarte E, Cubillos-Ruiz JR, Nesbeth YC, Scarlett UK, Martinez DG, Buckanovich RJ, Benencia F, Stan RV, Keler T, Sarobe P, Sentman CL, Conejo-Garcia JR. Depletion of dendritic cells delays ovarian cancer progression by boosting antitumor immunity. *Cancer Res*. Sep 15 2008;68(18):7684-7691.
64. Hung CF, Tsai YC, He L, Wu TC. Control of mesothelin-expressing ovarian cancer using adoptive transfer of mesothelin peptide-specific

- CD8(+) T cells. *Gene therapy*. Mar 22 2007.
65. Hutchens M, Luker KE, Sottile P, Sonstein J, Lukacs NW, Nunez G, Curtis JL, Luker GD. TLR3 Increases Disease Morbidity and Mortality from Vaccinia Infection. *J Immunol*. Jan 1 2008;180(1):483-491.
  66. Ilkovitch D, Lopez DM. The liver is a site for tumor-induced myeloid-derived suppressor cell accumulation and immunosuppression. *Cancer Res*. Jul 1 2009;69(13):5514-5521.
  67. Inglefield JR, Dumitru CD, Alkan SS, Gibson SJ, Lipson KE, Tomai MA, Larson CJ, Vasilakos JP. TLR7 agonist 852A inhibition of tumor cell proliferation is dependent on plasmacytoid dendritic cells and type I IFN. *J Interferon Cytokine Res*. Apr 2008;28(4):253-263.
  68. Izmailova ES, Paz N, Alencar H, Chun M, Schopf L, Hepperle M, Lane JH, Harriman G, Xu Y, Ocain T, Weissleder R, Mahmood U, Healy AM, Jaffee B. Use of molecular imaging to quantify response to IKK-2 inhibitor treatment in murine arthritis. *Arthritis and rheumatism*. Jan 2007;56(1):117-128.
  69. Jenq RR, King CG, Volk C, Suh D, Smith OM, Rao UK, Yim NL, Holland AM, Lu SX, Zakrzewski JL, Goldberg GL, Diab A, Alpdogan O, Penack O, Na IK, Kappel LW, Wolchok JD, Houghton AN, Perales MA, van den Brink MR. Keratinocyte growth factor enhances DNA plasmid tumor vaccine responses after murine allogeneic bone marrow transplantation. *Blood*. Feb 12 2009;113(7):1574-1580.
  70. Jeon YH, Choi Y, Kang JH, Kim CW, Jeong JM, Lee DS, Chung JK. Immune Response to Firefly Luciferase as a Naked DNA. *Cancer Biol Ther*. May 13 2007;6(5).
  71. Jing W, Gershan JA, Johnson BD. Depletion of CD4 T cells enhances immunotherapy for neuroblastoma after syngeneic HSCT but compromises development of antitumor immune memory. *Blood*. Apr 30 2009;113(18):4449-4457.
  72. Kalchenko V, Shvitiel S, Malina V, Lapid K, Haramati S, Lapidot T, Brill A, Harmelin A. Use of lipophilic near-infrared dye in whole-body optical imaging of hematopoietic cell homing. *Journal of biomedical optics*. Sep-Oct 2006;11(5):050507.
  73. Keller SA, Hernandez-Hopkins D, Vider J, Ponomarev V, Hyjek E, Schattner EJ, Cesarman E. NF-kappaB is essential for the progression of KSHV- and EBV-infected lymphomas in vivo. *Blood*. Apr 15 2006;107(8):3295-3302.
  74. Kim D, Hung CF, Wu TC. Monitoring the trafficking of adoptively transferred antigen- specific CD8-positive T cells in vivo, using noninvasive luminescence imaging. *Human gene therapy*. Jul 2007;18(7):575-588.
  75. Kim D, Monie A, He L, Tsai YC, Hung CF, Wu TC. Role of IL-2 secreted by PADRE-specific CD4+ T cells in enhancing E7-specific CD8+ T-cell immune responses. *Gene therapy*. May 2008;15(9):677-687.

76. Kim MH, Liu W, Borjesson DL, Curry FR, Miller LS, Cheung AL, Liu FT, Isseroff RR, Simon SI. Dynamics of Neutrophil Infiltration during Cutaneous Wound Healing and Infection Using Fluorescence Imaging. *J Invest Dermatol*. Jan 10 2008.
77. Kim WI, Wiesner SM, Largaespada DA. Vav promoter-tTA conditional transgene expression system for hematopoietic cells drives high level expression in developing B and T cells. *Exp Hematol*. Aug 2007;35(8):1231-1239.
78. Kim YJ, Dubey P, Ray P, Gambhir SS, Witte ON. Multimodality imaging of lymphocytic migration using lentiviral-based transduction of a tri-fusion reporter gene. *Molecular Imaging and Biology*. Sep-October 2004;6(5):331-340.
79. Kimura T, Koya RC, Anselmi L, Sternini C, Wang HJ, Comin-Anduix B, Prins RM, Faure-Kumar E, Rozengurt N, Cui Y, Kasahara N, Stripecke R. Lentiviral Vectors with CMV or MHCII Promoters Administered In Vivo: Immune Reactivity Versus Persistence of Expression. *Mol Ther*. May 1 2007.
80. Kolka JA, Vreede AP, Roessler BJ. Lipopolysaccharide recognition protein, MD-2, facilitates cellular uptake of E. coli-derived plasmid DNA in synovium. *J Gene Med*. March 2005;7(7):956-964.
81. Kruschinski A, Moosmann A, Poschke I, Norell H, Chmielewski M, Seliger B, Kiessling R, Blankenstein T, Abken H, Charo J. Engineering antigen-specific primary human NK cells against HER-2 positive carcinomas. *Proc Natl Acad Sci U S A*. Nov 11 2008;105(45):17481-17486.
82. Kung AL, Zabludoff SD, France DS, Freedman SJ, Tanner EA, Vieira A, Cornell-Kennon S, Lee J, Wang B, Wang J, Memmert K, Naegeli HU, Petersen F, Eck MJ, Bair KW, Wood AW, Livingston DM. Small molecule blockade of transcriptional coactivation of the hypoxia-inducible factor pathway. *Cancer Cell*. July 2004;6(1):33-43.
83. Kuns RD, Morris ES, Macdonald KP, Markey KA, Morris HM, Raffelt NC, Banovic T, Don AL, Rowe V, Burman AC, Clouston AD, Farah C, Besra GS, Illarionov PA, Smyth MJ, Porcelli SA, Hill GR. Invariant natural killer T cell-natural killer cell interactions dictate transplantation outcome after alpha-galactosylceramide administration. *Blood*. Jun 4 2009;113(23):5999-6010.
84. Kutschka I, Chen IY, Kofidis T, Arai T, von Degenfeld G, Sheikh AY, Hendry SL, Pearl J, Hoyt G, Sista R, Yang PC, Blau HM, Gambhir SS, Robbins RC. Collagen Matrices Enhance Survival of Transplanted Cardiomyoblasts and Contribute to Functional Improvement of Ischemic Rat Hearts. *Circulation*. July 4, 2006 2006;114(1\_suppl):I-167-173.
85. Kutschka I, Kofidis T, Chen IY, von Degenfeld G, Zwierzchoniewska M, Hoyt G, Arai T, Lebl DR, Hendry SL, Sheikh AY, Cooke DT, Connolly A, Blau HM, Gambhir SS, Robbins RC. Adenoviral Human BCL-2



- Transgene Expression Attenuates Early Donor Cell Death After Cardiomyoblast Transplantation Into Ischemic Rat Hearts. *Circulation*. July 4, 2006 2006;114(1\_suppl):I-174-180.
86. Lapteva N, Seethammagari MR, Hanks BA, Jiang J, Levitt JM, Slawin KM, Spencer DM. Enhanced activation of human dendritic cells by inducible CD40 and Toll-like receptor-4 ligation. *Cancer Res*. Nov 1 2007;67(21):10528-10537.
  87. Lauritzen KH, Moldestad O, Eide L, Carlsen H, Nesse G, Storm JF, Mansuy IM, Bergersen LH, Klungland A. Mitochondrial DNA toxicity in forebrain neurons causes apoptosis, neurodegeneration, and impaired behavior. *Molecular and cellular biology*. Mar 2010;30(6):1357-1367.
  88. Le Devedec SE, van Roosmalen W, Maria N, Grimbergen M, Pont C, Lalai R, van de Water B. An improved model to study tumor cell autonomous metastasis programs using MTLn3 cells and the Rag2(-/-) gammac (-/-) mouse. *Clin Exp Metastasis*. 2009;26(7):673-684.
  89. Lee SH, Hao E, Savinov AY, Geron I, Strongin AY, Itkin-Ansari P. Human beta-cell precursors mature into functional insulin-producing cells in an immunoisolation device: implications for diabetes cell therapies. *Transplantation*. Apr 15 2009;87(7):983-991.
  90. Leo BM, Li X, Balian G, Anderson DG. In vivo bioluminescent imaging of virus-mediated gene transfer and transduced cell transplantation in the intervertebral disc. *Spine*. April 2004;29(8):838-844.
  91. Li D, Gromov K, Soballe K, Puzas JE, O'Keefe RJ, Awad H, Drissi H, Schwarz EM. Quantitative mouse model of implant-associated osteomyelitis and the kinetics of microbial growth, osteolysis, and humoral immunity. *J Orthop Res*. Jan 2008;26(1):96-105.
  92. Li Z, Suzuki Y, Huang M, Cao F, Xie X, Connolly AJ, Yang PC, Wu JC. Comparison of Reporter Gene and Iron Particle Labeling for Tracking Fate of Human Embryonic Stem Cells and Differentiated Endothelial Cells in Living Subjects. *Stem Cells*. Jan 24 2008.
  93. Liang Q, Yamamoto M, Curiel DT, Herschman HR. Noninvasive imaging of transcriptionally restricted transgene expression following intratumoral injection of an adenovirus in which the COX-2 promoter drives a reporter gene. *Mol Imaging Biol*. Nov-Dec 2004;6(6):395-404.
  94. Liang X, Romo-Vivar A, Schapiro NE, Loughran P, Thorne SH, Amoscato AA, Zeh HJ, Beer-Stolz D, Lotze MT, de Vera ME. Ethyl pyruvate administration inhibits hepatic tumor growth. *Journal of leukocyte biology*. Jul 7 2009.
  95. Liang Y, Liu C, Djeu JY, Zhong B, Peters T, Scharffetter-Kochanek K, Anasetti C, Yu XZ. Beta2 integrins separate graft-versus-host disease and graft-versus-leukemia effects. *Blood*. Jan 15 2008;111(2):954-962.
  96. Limberis MP, Bell CL, Heath J, Wilson JM. Activation of transgene-specific T cells following lentivirus-mediated gene delivery to mouse lung. *Mol Ther*. Jan 2010;18(1):143-150.

97. Ling X, Konopleva M, Zeng Z, Ruvolo V, Stephens LC, Schober W, McQueen T, Dietrich M, Madden TL, Andreeff M. The novel triterpenoid C-28 methyl ester of 2-cyano-3, 12-dioxoolen-1, 9-dien-28-oic acid inhibits metastatic murine breast tumor growth through inactivation of STAT3 signaling. *Cancer Res.* May 1 2007;67(9):4210-4218.
98. Liu C, Yu S, Kappes J, Wang J, Grizzle WE, Zinn KR, Zhang HG. Expansion of spleen myeloid suppressor cells represses NK cell cytotoxicity in tumor-bearing host. *Blood.* May 15 2007;109(10):4336-4342.
99. Lu Y, Dang H, Middleton B, Zhang Z, Washburn L, Campbell-Thompson M, Atkinson MA, Gambhir SS, Tian J, Kaufman DL. Bioluminescent monitoring of islet graft survival after transplantation. *Mol Ther.* March 2004;9(3):428-435.
100. Lundqvist A, Su S, Rao S, Childs R. Cutting edge: Bortezomib-treated tumors sensitized to NK cell apoptosis paradoxically acquire resistance to antigen-specific T cells. *J Immunol.* Feb 1 2010;184(3):1139-1142.
101. Luo J, Ho P, Steinman L, Wyss-Coray T. Bioluminescence in vivo imaging of autoimmune encephalomyelitis predicts disease. *Journal of neuroinflammation.* 2008;5:6.
102. Luo J, Ho PP, Buckwalter MS, Hsu T, Lee LY, Zhang H, Kim DK, Kim SJ, Gambhir SS, Steinman L, Wyss-Coray T. Glia-dependent TGF-beta signaling, acting independently of the TH17 pathway, is critical for initiation of murine autoimmune encephalomyelitis. *J Clin Invest.* Nov 1 2007;117(11):3306-3315.
103. Ma L, Xiang Z, Sherrill TP, Wang L, Blackwell TS, Williams P, Chong A, Chari R, Yin DP. Bioluminescence Imaging Visualizes Activation of Nuclear Factor-kappaB in Mouse Cardiac Transplantation. *Transplantation.* Mar 27 2008;85(6):903-910.
104. Maes W, Deroose C, Reumers V, Krylyshkina O, Gijsbers R, Baekelandt V, Ceuppens J, Debyser Z, Van Gool SW. In vivo bioluminescence imaging in an experimental mouse model for dendritic cell based immunotherapy against malignant glioma. *Journal of neuro-oncology.* Jan 2009;91(2):127-139.
105. Malley R, Henneke P, Morse SC, Cieslewicz MJ, Lipsitch M, Thompson CM, Kurt-Jones E, Paton JC, Wessels MR, Golenbock DT. Recognition of pneumolysin by Toll-like receptor 4 confers resistance to pneumococcal infection. *Proc Natl Acad Sci U S A.* February 2003;100(4):1966-1971.
106. Mandl S, Schimmelpfennig C, Edinger M, Negrin RS, Contag CH. Understanding immune cell trafficking patterns via in vivo bioluminescence imaging. *J Cell Biochem Suppl.* 2002;39:239-248.
107. Mandl SJ, Mari C, Edinger M, Negrin RS, Tait JF, Contag CH, Blankenberg FG. Multi-modality imaging identifies key times for annexin V imaging as an early predictor of therapeutic outcome. *Mol Imaging.*

- January 2004;3(1):1-8.
108. Markey KA, Burman AC, Banovic T, Kuns RD, Raffelt NC, Rowe V, Olver SD, Don AL, Morris ES, Pettit AR, Wilson YA, Robb RJ, Randall LM, Korner H, Engwerda CR, Clouston AD, Macdonald KP, Hill GR. Soluble lymphotoxin is an important effector molecule in GVHD and GVL. *Blood*. Jan 7 2010;115(1):122-132.
  109. Marttila-Ichihara F, Turja R, Miiluniemi M, Karikoski M, Maksimow M, Niemela J, Martinez-Pomares L, Salmi M, Jalkanen S. Macrophage mannose receptor on lymphatics controls cell trafficking. *Blood*. Jul 1 2008;112(1):64-72.
  110. Mason KM, Munson RS, Jr., Bakaletz LO. A mutation in the sap operon attenuates survival of nontypeable *Haemophilus influenzae* in a chinchilla model of otitis media. *Infect Immun*. January 2005;73(1):599-608.
  111. Mays LE, Vandenberghe LH, Xiao R, Bell P, Nam HJ, Agbandje-McKenna M, Wilson JM. Adeno-associated virus capsid structure drives CD4-dependent CD8+ T cell response to vector encoded proteins. *J Immunol*. May 15 2009;182(10):6051-6060.
  112. Min JJ AY, Moon S, Kim YS, Park JE, Kim SM, Le UN, Wu JC, Joo SY, Hong MH, Yang DH, Jeong MH, Song CH, Jeong YH, Yoo KY, Kang KS, Bom HS. In vivo bioluminescence imaging of cord blood derived mesenchymal stem cell transplantation into rat myocardium. *Annals of Nuclear Medicine*. 2006;20(3):165-170.
  113. Mirmonsef P, Tan G, Zhou G, Morino T, Noonan K, Borrello I, Levitsky HI. Escape from suppression: tumor-specific effector cells outcompete regulatory T cells following stem-cell transplantation. *Blood*. Feb 15 2008;111(4):2112-2121.
  114. Mitsiades CS, Mitsiades NS, McMullan CJ, Poulaki V, Shringarpure R, Akiyama M, Hideshima T, Chauhan D, Joseph M, Libermann TA, Garcia-Echeverria C, Pearson MA, Hofmann F, Anderson KC, Kung AL. Inhibition of the insulin-like growth factor receptor-1 tyrosine kinase activity as a therapeutic strategy for multiple myeloma, other hematologic malignancies, and solid tumors. *Cancer Cell*. March 2004;5(3):221-230.
  115. Moriyama Y, Moriyama EH, Blackmore K, Akens MK, Lilge L. In vivo Study of the Inflammatory Modulating Effects of Low Level Laser Therapy on iNOS Expression Using Bioluminescence Imaging. *Photochem Photobiol*. Feb 1 2005.
  116. Murai K, Sakai D, Nakamura Y, Nakai T, Igarashi T, Seo N, Murakami T, Kobayashi E, Mochida J. Primary immune system responders to nucleus pulposus cells: evidence for immune response in disc herniation. *European cells & materials*. 2010;19:13-21.
  117. Murakami T, Sato A, Chun NA, Hara M, Naito Y, Kobayashi Y, Kano Y, Ohtsuki M, Furukawa Y, Kobayashi E. Transcriptional modulation using

- HDACi depsipeptide promotes immune cell-mediated tumor destruction of murine B16 melanoma. *J Invest Dermatol.* Jun 2008;128(6):1506-1516.
118. Nakajima A, Kamijo S, Yoshioka T. Visualization of antigen-specific T cell in living arthritic mice. *Journal of Nippon Medical School = Nihon Ika Daigaku zasshi.* Dec 2006;73(6):298-299.
  119. Nakajima A, Seroogy CM, Sandora MR, Tarner IH, Costa GL, Taylor-Edwards C, Bachmann MH, Contag CH, Fathman CG. Antigen-specific T cell-mediated gene therapy in collagen-induced arthritis. *J Clin Invest.* May 2001;107(10):1293-1301.
  120. Nervi B, Rettig MP, Ritchey JK, Wang HL, Bauer G, Walker J, Bonyhadi ML, Berenson RJ, Prior JL, Piwnica-Worms D, Nolte JA, Dipersio JF. Factors affecting human T cell engraftment, trafficking, and associated xenogeneic graft-vs-host disease in NOD/SCID beta2m(null) mice. *Exp Hematol.* Aug 29 2007.
  121. Ni J, Nolte B, Arnold A, Fournier P, Schirmacher V. Targeting anti-tumor DNA vaccines to dendritic cells via a short CD11c promoter sequence. *Vaccine.* Sep 4 2009;27(40):5480-5487.
  122. Nishimura R, Baker J, Beilhack A, Zeiser R, Olson JA, Segal EI, Karimi M, Negrin RS. In vivo trafficking and survival of cytokine-induced killer cells resulting in minimal GVHD with retention of antitumor activity. *Blood.* Sep 15 2008;112(6):2563-2574.
  123. O'Mahony C, Scully P, O'Mahony D, Murphy S, O'Brien F, Lyons A, Sherlock G, MacSharry J, Kiely B, Shanahan F, O'Mahony L. Commensal-induced regulatory T cells mediate protection against pathogen-stimulated NF-kappaB activation. *PLoS pathogens.* 2008;4(8):e1000112.
  124. Ogunniyi AD, Paton JC, Kirby AC, McCullers JA, Cook J, Hyodo M, Hayakawa Y, Karaolis DK. c-di-GMP is an effective immunomodulator and vaccine adjuvant against pneumococcal infection. *Vaccine.* Aug 26 2008;26(36):4676-4685.
  125. Ojogho O, Li A, Zuppan CW, Escher A. Bioluminescent mammalian cells grown in sponge matrices to monitor immune rejection. *Mol Imaging.* Sep-Oct 2007;6(5):323-330.
  126. Okada S, Ishii K, Yamane J, Iwanami A, Ikegami T, Katoh H, Iwamoto Y, Nakamura M, Miyoshi H, Okano HJ, Contag CH, Toyama Y, Okano H. In vivo imaging of engrafted neural stem cells: its application in evaluating the optimal timing of transplantation for spinal cord injury. *Faseb J.* Sep 1 2005.
  127. Olson JA, Zeiser R, Beilhack A, Goldman JJ, Negrin RS. Tissue-specific homing and expansion of donor NK cells in allogeneic bone marrow transplantation. *J Immunol.* Sep 1 2009;183(5):3219-3228.
  128. Pham W, Xie J, Gore JC. Tracking the migration of dendritic cells by in vivo optical imaging. *Neoplasia (New York, N.Y.)* Dec

- 2007;9(12):1130-1137.
129. Phillips B, Nylander K, Harnaha J, Machen J, Lakomy R, Styche A, Gillis K, Brown L, Lafreniere D, Gallo M, Knox J, Hogeland K, Trucco M, Giannoukakis N. A microsphere-based vaccine prevents and reverses new-onset autoimmune diabetes. *Diabetes*. Jun 2008;57(6):1544-1555.
  130. Prins R, Bruhn KW, Craft N, Lin JW, Kim CH, Odesa SK, Miller JF, Liau LM. Central nervous system tumor immunity generated by a recombinant listeria monocytogenes vaccine targeting tyrosinase related protein-2 and real-time imaging of intracranial tumor burden. *Neurosurgery*. 2006;58(1):169-178.
  131. Prins RM, Craft N, Bruhn KW, Khan-Farooqi H, Koya RC, Stripecke R, Miller JF, Liau LM. The TLR-7 agonist, imiquimod, enhances dendritic cell survival and promotes tumor antigen-specific T cell priming: relation to central nervous system antitumor immunity. *J Immunol*. Jan 1 2006;176(1):157-164.
  132. Prins RM, Shu CJ, Radu CG, Vo DD, Khan-Farooqi H, Soto H, Yang MY, Lin MS, Shelly S, Witte ON, Ribas A, Liau LM. Anti-tumor activity and trafficking of self, tumor-specific T cells against tumors located in the brain. *Cancer Immunol Immunother*. Sep 2008;57(9):1279-1289.
  133. Prins RM, Vo DD, Khan-Farooqi H, Yang MY, Soto H, Economou JS, Liau LM, Ribas A. NK and CD4 cells collaborate to protect against melanoma tumor formation in the brain. *J Immunol*. Dec 15 2006;177(12):8448-8455.
  134. Rabinovich BA, Ye Y, Etto T, Chen JQ, Levitsky HI, Overwijk WW, Cooper LJ, Gelovani J, Hwu P. Visualizing fewer than 10 mouse T cells with an enhanced firefly luciferase in immunocompetent mouse models of cancer. *Proc Natl Acad Sci U S A*. Sep 23 2008;105(38):14342-14346.
  135. Rabinovich PM, Komarovskaya ME, Wrzesinski SH, Alderman JL, Budak-Alpdogan T, Karpikov A, Guo H, Flavell RA, Cheung NK, Weissman SM, Bahceci E. Chimeric Receptor mRNA Transfection as a Tool to Generate Antineoplastic Lymphocytes. *Human gene therapy*. Oct 15 2008.
  136. Ramakrishna E, Woller N, Mundt B, Knocke S, Gurlevik E, Saborowski M, Malek N, Manns MP, Wirth T, Kuhnel F, Kubicka S. Antitumoral immune response by recruitment and expansion of dendritic cells in tumors infected with telomerase-dependent oncolytic viruses. *Cancer Res*. Feb 15 2009;69(4):1448-1458.
  137. Ramirez-Montagut T, Chow A, Kochman AA, Smith OM, Suh D, Sindhi H, Lu S, Borsotti C, Grubin J, Patel N, Terwey TH, Kim TD, Heller G, Murphy GF, Liu C, Alpdogan O, van den Brink MR. IFN-gamma and Fas ligand are required for graft-versus-tumor activity against renal cell carcinoma in the absence of lethal graft-versus-host disease. *J Immunol*. Aug 1 2007;179(3):1669-1680.

138. Ramphal R, Balloy V, Jyot J, Verma A, Si-Tahar M, Chignard M. Control of *Pseudomonas aeruginosa* in the lung requires the recognition of either lipopolysaccharide or flagellin. *J Immunol.* Jul 1 2008;181(1):586-592.
139. Rauch D, Gross S, Harding J, Niewiesk S, Lairmore M, Piwnica-Worms D, Ratner L. Imaging spontaneous tumorigenesis: Inflammation precedes development of peripheral NK tumors. *Blood.* Oct 29 2008.
140. Real PJ, Tosello V, Palomero T, Castillo M, Hernando E, de Stanchina E, Sulis ML, Barnes K, Sawai C, Homminga I, Meijerink J, Aifantis I, Basso G, Cordon-Cardo C, Ai W, Ferrando A. Gamma-secretase inhibitors reverse glucocorticoid resistance in T cell acute lymphoblastic leukemia. *Nature medicine.* Jan 2009;15(1):50-58.
141. Reichardt W, Durr C, von Elverfeldt D, Juttner E, Gerlach UV, Yamada M, Smith B, Negrin RS, Zeiser R. Impact of mammalian target of rapamycin inhibition on lymphoid homing and tolerogenic function of nanoparticle-labeled dendritic cells following allogeneic hematopoietic cell transplantation. *J Immunol.* Oct 1 2008;181(7):4770-4779.
142. Rivera R, Hutchens M, Luker KE, Sonstein J, Curtis JL, Luker GD. Murine alveolar macrophages limit replication of vaccinia virus. *Virology.* Feb 26 2007.
143. Roth DJ, Jansen ED, Powers AC, Wang TG. A Novel Method of Monitoring Response to Islet Transplantation: Bioluminescent Imaging of an NF- $\kappa$ B Transgenic Mouse Model. [Article]. *Transplantation.* 2006.
144. Ryan A, Scrable H. Visualization of the dynamics of gene expression in the living mouse. *Mol Imaging.* January 2004;3(1):33-42.
145. Sato A, Ohtsuki M, Hata M, Kobayashi E, Murakami T. Antitumor activity of IFN- $\lambda$  in murine tumor models. *J Immunol.* Jun 15 2006;176(12):7686-7694.
146. Scheffold C, Kornacker M, Scheffold YC, Contag CH, Negrin RS. Visualization of effective tumor targeting by CD8<sup>+</sup> natural killer T cells redirected with bispecific antibody F(ab')<sub>2</sub>HER2xCD3. *Cancer Res.* October 2002;62(20):5785-5791.
147. Schimmelpfennig CH, Schulz S, Arber C, Baker J, Tarner I, McBride J, Contag CH, Negrin RS. Ex Vivo Expanded Dendritic Cells Home to T-Cell Zones of Lymphoid Organs and Survive in Vivo after Allogeneic Bone Marrow Transplantation. *Am J Pathol.* November 1, 2005 2005;167(5):1321-1331.
148. Schrepfer S, Deuse T, Reichenspurner H, Fischbein MP, Robbins RC, Pelletier MP. Stem cell transplantation: the lung barrier. *Transplantation proceedings.* Mar 2007;39(2):573-576.
149. Schweichel D, Steitz J, Tormo D, Gaffal E, Ferrer A, Buchs S, Speuser P, Limmer A, Tuting T. Evaluation of DNA vaccination with recombinant adenoviruses using bioluminescence imaging of antigen expression: impact of application routes and delivery with dendritic cells. *J Gene*

- Med.* Oct 2006;8(10):1243-1250.
150. Singh H, Serrano LM, Pfeiffer T, Olivares S, McNamara G, Smith DD, Al-Kadhimi Z, Forman SJ, Gillies SD, Jensen MC, Colcher D, Raubitschek A, Cooper LJ. Combining adoptive cellular and immunocytokine therapies to improve treatment of B-lineage malignancy. *Cancer Res.* Mar 15 2007;67(6):2872-2880.
  151. Slavin AJ, Tarner IH, Nakajima A, Urbanek-Ruiz I, McBride J, Contag CH, Fathman CG. Adoptive cellular gene therapy of autoimmune disease. *Autoimmun Rev.* August 2002;1(4):213-219.
  152. Stephan MT, Ponomarev V, Brentjens RJ, Chang AH, Dobrenkov KV, Heller G, Sadelain M. T cell-encoded CD80 and 4-1BBL induce auto- and transcostimulation, resulting in potent tumor rejection. *Nature medicine.* Dec 2007;13(12):1440-1449.
  153. Suffner J, Hochweller K, Kuhnle MC, Li X, Kroczeck RA, Garbi N, Hammerling GJ. Dendritic cells support homeostatic expansion of Foxp3+ regulatory T cells in Foxp3.LuciDTR mice. *J Immunol.* Feb 15 2010;184(4):1810-1820.
  154. Tanaka M, Swijnenburg R-J, Gunawan F, Cao Y-A, Yang Y, Caffarelli AD, de Bruin JL, Contag CH, Robbins RC. In Vivo Visualization of Cardiac Allograft Rejection and Trafficking Passenger Leukocytes Using Bioluminescence Imaging. *Circulation.* August 30, 2005 2005;112(9\_suppl):I-105-110.
  155. Tarner IH, Slavin AJ, McBride J, Levicnik A, Smith R, Nolan GP, Contag CH, Fathman CG. Treatment of autoimmune disease by adoptive cellular gene therapy. *Ann N Y Acad Sci.* September 2003;998:512-519.
  156. Terwey TH, Kim TD, Kochman AA, Hubbard VM, Lu S, Zakrzewski JL, Ramirez-Montagut T, Eng JM, Muriglian SJ, Heller G, Murphy GF, Liu C, Budak-Alpdogan T, Alpdogan O, van den Brink MR. CCR2 is required for CD8-induced graft-versus-host disease. *Blood.* Nov 1 2005;106(9):3322-3330.
  157. Tian X, Woll PS, Morris JK, Linehan JL, Kaufman DS. Hematopoietic engraftment of human embryonic stem cell-derived cells is regulated by recipient innate immunity. *Stem Cells.* February 2, 2006 2006:2005-0340.
  158. Tolar J, O'Shaughnessy MJ, Panoskaltsis-Mortari A, McElmurry RT, Bell S, Riddle M, Mclvor RS, Yant SR, Kay MA, Krause D, Verfaillie CM, Blazar BR. Host factors that impact the biodistribution and persistence of multipotent adult progenitor cells. *Blood.* January 12, 2006 2006:2005-2008-3289.
  159. Tolar J, Wang X, Braunlin E, McElmurry RT, Nakamura Y, Bell S, Xia L, Zhang J, Hu Q, Panoskaltsis-Mortari A, Zhang J, Blazar BR. The host immune response is essential for the beneficial effect of adult stem cells after myocardial ischemia. *Exp Hematol.* Apr 2007;35(4):682-690.
  160. Tseng CW, Hung CF, Alvarez RD, Trimble C, Huh WK, Kim D, Chuang

- CM, Lin CT, Tsai YC, He L, Monie A, Wu TC. Pretreatment with cisplatin enhances E7-specific CD8+ T-Cell-mediated antitumor immunity induced by DNA vaccination. *Clin Cancer Res.* May 15 2008;14(10):3185-3192.
- 161.** Ushiki T, Kizaka-Kondoh S, Ashihara E, Tanaka S, Masuko M, Hirai H, Kimura S, Aizawa Y, Maekawa T, Hiraoka M. Noninvasive tracking of donor cell homing by near-infrared fluorescence imaging shortly after bone marrow transplantation. *PLoS one.* 2010;5(6):e11114.
- 162.** Varol C, Landsman L, Fogg DK, Greenshtein L, Gildor B, Margalit R, Kalchenko V, Geissmann F, Jung S. Monocytes give rise to mucosal, but not splenic, conventional dendritic cells. *The Journal of experimental medicine.* Jan 22 2007;204(1):171-180.
- 163.** Vera J, Savoldo B, Vigouroux S, Biagi E, Pule M, Rossig C, Wu J, Heslop HE, Rooney CM, Brenner MK, Dotti G. T lymphocytes redirected against the kappa light chain of human immunoglobulin efficiently kill mature B lymphocyte-derived malignant cells. *Blood.* Dec 1 2006;108(12):3890-3897.
- 164.** Verra NC, Jorritsma A, Weijer K, Ruizendaal JJ, Voordouw A, Weder P, Hooijberg E, Schumacher TN, Haanen JB, Spits H, Luiten RM. Human telomerase reverse transcriptase-transduced human cytotoxic T cells suppress the growth of human melanoma in immunodeficient mice. *Cancer Res.* March 2004;64(6):2153-2161.
- 165.** Vikis HG, Jackson EN, Krupnick AS, Franklin A, Gelman AE, Chen Q, Piwnicka-Worms D, You M. Strain-specific susceptibility for pulmonary metastasis of sarcoma 180 cells in inbred mice. *Cancer Res.* Jun 15 2010;70(12):4859-4867.
- 166.** Vishwanath RP, Brown CE, Wagner JR, Meechoovet HB, Naranjo A, Wright CL, Olivares S, Qian D, Cooper LJ, Jensen MC. A quantitative high-throughput chemotaxis assay using bioluminescent reporter cells. *Journal of immunological methods.* Jul 2005;302(1-2):78-89.
- 167.** von Kockritz-Blickwede M, Rohde M, Oehmcke S, Miller LS, Cheung AL, Herwald H, Foster S, Medina E. Immunological mechanisms underlying the genetic predisposition to severe *Staphylococcus aureus* infection in the mouse model. *Am J Pathol.* Dec 2008;173(6):1657-1668.
- 168.** Vuong C, Kocianova S, Yu J, Kadurugamuwa JL, Otto M. Development of real-time in vivo imaging of device-related *Staphylococcus epidermidis* infection in mice and influence of animal immune status on susceptibility to infection. *The Journal of infectious diseases.* Jul 15 2008;198(2):258-261.
- 169.** Waldman E, Lu SX, Hubbard VM, Kochman AA, Eng JM, Terwey TH, Muriglan SJ, Kim TD, Heller G, Murphy GF, Liu C, Alpdogan O, van den Brink MR. Absence of beta7 integrin results in less graft-versus-host disease because of decreased homing of alloreactive T cells to intestine. *Blood.* Feb 15 2006;107(4):1703-1711.



170. Walensky LD, Kung AL, Escher I, Malia TJ, Barbuto S, Wright RD, Wagner G, Verdine GL, Korsmeyer SJ. Activation of apoptosis in vivo by a hydrocarbon-stapled BH3 helix. *Science*. September 2004;305(5689):1466-1470.
171. Wan YY, Flavell RA. Regulatory T-cell functions are subverted and converted owing to attenuated Foxp3 expression. *Nature*. Jan 14 2007.
172. Wang X, Rosol M, Ge S, Peterson D, McNamara G, Pollack H, Kohn DB, Nelson MD, Crooks GM. Dynamic tracking of human hematopoietic stem cell engraftment using in vivo bioluminescence imaging. *Blood*. November 2003;102(10):3478-3482.
173. Woll PS, Grzywacz B, Tian X, Marcus RK, Knorr DA, Verneris MR, Kaufman DS. Human embryonic stem cells differentiate into a homogeneous population of natural killer cells with potent in vivo antitumor activity. *Blood*. Jun 11 2009;113(24):6094-6101.
174. Xie D, Gore C, Liu J, Pong RC, Mason R, Hao G, Long M, Kabbani W, Yu L, Zhang H, Chen H, Sun X, Boothman DA, Min W, Hsieh JT. Role of DAB2IP in modulating epithelial-to-mesenchymal transition and prostate cancer metastasis. *Proc Natl Acad Sci U S A*. Feb 9 2010;107(6):2485-2490.
175. Zaharoff DA, Rogers CJ, Hance KW, Schlom J, Greiner JW. Chitosan solution enhances the immunoadjuvant properties of GM-CSF. *Vaccine*. Dec 17 2007;25(52):8673-8686.
176. Zakrzewski JL, Kochman AA, Lu SX, Terwey TH, Kim TD, Hubbard VM, Muriglan SJ, Suh D, Smith OM, Grubin J, Patel N, Chow A, Cabrera-Perez J, Radhakrishnan R, Diab A, Perales MA, Rizzuto G, Menet E, Pamer EG, Heller G, Zuniga-Pflucker JC, Alpdogan O, van den Brink MR. Adoptive transfer of T-cell precursors enhances T-cell reconstitution after allogeneic hematopoietic stem cell transplantation. *Nature medicine*. Sep 2006;12(9):1039-1047.
177. Zakrzewski JL, Suh D, Markley JC, Smith OM, King C, Goldberg GL, Jenq R, Holland AM, Grubin J, Cabrera-Perez J, Brentjens RJ, Lu SX, Rizzuto G, Sant'Angelo DB, Riviere I, Sadelain M, Heller G, Zuniga-Pflucker JC, Lu C, van den Brink MR. Tumor immunotherapy across MHC barriers using allogeneic T-cell precursors. *Nature biotechnology*. Apr 2008;26(4):453-461.
178. Zangani M, Carlsen H, Kielland A, Os A, Hauglin H, Blomhoff R, Munthe LA, Bogen B. Tracking early autoimmune disease by bioluminescent imaging of NF-kappaB activation reveals pathology in multiple organ systems. *Am J Pathol*. Apr 2009;174(4):1358-1367.
179. Zeiser R, Leveson-Gower DB, Zambricki EA, Kambham N, Beilhack A, Loh J, Hou JZ, Negrin RS. Differential impact of mammalian target of rapamycin inhibition on CD4+CD25+Foxp3+ regulatory T cells compared with conventional CD4+ T cells. *Blood*. Jan 1 2008;111(1):453-462.

180. Zeiser R, Nguyen VH, Beilhack A, Buess M, Schulz S, Baker J, Contag CH, Negrin RS. Inhibition of CD4+CD25+ regulatory T-cell function by calcineurin-dependent interleukin-2 production. *Blood*. July 1, 2006 2006;108(1):390-399.
181. Zeiser R, Nguyen VH, Hou J-Z, Beilhack A, Zambricki EA, Buess M, Contag CH, Negrin RS. Early CD30 signaling is critical for adoptively transferred CD4+CD25+ regulatory T cells in prevention of acute graft versus host disease. *Blood*. October 26, 2006 2006:blood-2006-2007-038455.
182. Zeiser R, Nguyen VH, Hou JZ, Beilhack A, Zambricki E, Buess M, Contag CH, Negrin RS. Early CD30 signaling is critical for adoptively transferred CD4+CD25+ regulatory T cells in prevention of acute graft-versus-host disease. *Blood*. Mar 1 2007;109(5):2225-2233.
183. Zeiser R, Zambricki EA, Leveson-Gower D, Kambham N, Beilhack A, Negrin RS. Host-derived interleukin-18 differentially impacts regulatory and conventional T cell expansion during acute graft-versus-host disease. *Biol Blood Marrow Transplant*. Dec 2007;13(12):1427-1438.
184. Zhang GJ, Safran M, Wei W, Sorensen E, Lassota P, Zhelev N, Neuberg DS, Shapiro G, Kaelin WG, Jr. Bioluminescent imaging of Cdk2 inhibition in vivo. *Nature medicine*. June 2004;10(6):643-648.
185. Zhang N, Weber A, Li B, Lyons R, Contag PR, Purchio AF, West DB. An inducible nitric oxide synthase-luciferase reporter system for in vivo testing of anti-inflammatory compounds in transgenic mice. *J Immunol*. June 2003;170(12):6307-6319.
186. Zhang YQ, Tsai YC, Monie A, Wu TC, Hung CF. Enhancing the Therapeutic Effect Against Ovarian Cancer Through a Combination of Viral Oncolysis and Antigen-specific Immunotherapy. *Mol Ther*. Jan 19 2010.